




# **CGK Surgical**

## **Product Catalog**





## Electrosurgery


|                                      |             |
|--------------------------------------|-------------|
| BipoJET® Bipolar Scissors .....      | 1168        |
| Non-Irrigating Bipolar Forceps ..... | 1169 - 1190 |
| Irrigating Bipolar Forceps .....     | 1191 - 1194 |
| Disposable Bipolar Forceps .....     | 1195 - 1197 |
| Bipolar Connecting Cords .....       | 1198 - 1200 |

# Electrosurgery


## BipoJET® Bipolar Scissors

The versatility of BipoJET Dissecting Bipolar Scissors enables:


■ Parallel Cutting and Coagulating



■ Precoagulation



■ Pinpoint Coagulation



■ Surface Coagulation



For more information about Bipolar Scissors see our brochure number DOC456.  
See page 1199-1200 for appropriate connecting cords.



BC691R-BC693R



BC691R-BC693R

### METZENBAUM Bipolar Dissecting Scissors

| Catalog Number | Tip    | Length (in.) |
|----------------|--------|--------------|
| BC691R         | X-FINE | 4 3/4        |
| BC692R         | X-FINE | 5 1/2        |
| BC693R         | X-FINE | 6 1/4        |



BC694R-BC696R



BC694R-BC696R

### STEVENS Bipolar Dissecting Scissors

| Catalog Number | Tip      | Length (in.) |
|----------------|----------|--------------|
| BC694R         | TENOTOMY | 4 3/4        |
| BC695R         | TENOTOMY | 5 1/2        |
| BC696R         | TENOTOMY | 6 1/4        |



BC685R-BC690R



BC685R




BC686R

### METZENBAUM Bipolar Dissecting Scissors

| Catalog Number | Tip      | Length (in.) |
|----------------|----------|--------------|
| BC685R         | STANDARD | 7            |
| BC687R         | STANDARD | 9            |
| BC689R         | STANDARD | 11           |
| BC686R         | FINE     | 7            |
| BC688R         | FINE     | 9            |
| BC690R         | FINE     | 11           |

# Electrosurgery

## Non-Irrigating Bipolar Forceps




### YASARGIL Bipolar Forceps Features:

- High precision, tissue preserving preparation instrument
- Stronger tension of the forceps allows for blunt spreading of the tissue.

### Tip Size Reference Chart

|         |  |        |  |
|---------|--|--------|--|
| 0.2 mm  |  | 0.8 mm |  |
| 0.25 mm |  | 0.9 mm |  |
| 0.4 mm  |  | 1.0 mm |  |
| 0.5 mm  |  | 1.3 mm |  |
| 0.6 mm  |  | 1.5 mm |  |
| 0.7 mm  |  | 2.0 mm |  |



### YASARGIL Style Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| GK746R         | 0.4           | BAYONET | INS*       | NON-IRRIG* | 3/4                  | 5 1/2        |
| GK747R         | 0.7           | BAYONET | INS        | NON-IRRIG  | 3/4                  | 5 1/2        |
| GK748R         | 1.0           | BAYONET | INS        | NON-IRRIG  | 3/4                  | 5 1/2        |
| GK749R         | 1.3           | BAYONET | INS        | NON-IRRIG  | 3/4                  | 5 1/2        |
| GK760R         | 0.4           | BAYONET | INS        | NON-IRRIG  | 1 3/8                | 6 1/8        |
| GK761R         | 0.7           | BAYONET | INS        | NON-IRRIG  | 1 3/8                | 6 1/8        |
| GK762R         | 1.0           | BAYONET | INS        | NON-IRRIG  | 1 3/8                | 6 1/8        |
| GK763R         | 1.3           | BAYONET | INS        | NON-IRRIG  | 1 3/8                | 6 1/8        |
| GK764R         | 0.4           | BAYONET | INS        | NON-IRRIG  | 2 1/4                | 6 7/8        |
| GK765R         | 0.7           | BAYONET | INS        | NON-IRRIG  | 2 1/4                | 6 7/8        |
| GK766R         | 1.0           | BAYONET | INS        | NON-IRRIG  | 2 1/4                | 6 7/8        |
| GK767R         | 1.3           | BAYONET | INS        | NON-IRRIG  | 2 1/4                | 6 7/8        |

\* Key:

INS = Insulating

NON-IRRIG = Non-Irrigating






IRRIG = Irrigating

Electrosurgery

Non-Irrigating Bipolar Forceps

YASARGIL Style Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| GK768R         | 0.4           | BAYONET | INS        | NON-IRRIG  | 3                    | 7 3/4        |
| GK769R         | 0.7           | BAYONET | INS        | NON-IRRIG  | 3                    | 7 3/4        |
| GK770R         | 1.0           | BAYONET | INS        | NON-IRRIG  | 3                    | 7 3/4        |
| GK771R         | 1.3           | BAYONET | INS        | NON-IRRIG  | 3                    | 7 3/4        |
| GK772R         | 0.4           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |
| GK773R         | 0.7           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |
| GK774R         | 1.0           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |
| GK775R         | 1.3           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |









ELECTRO  
SURGERY

Electrosurgery  
Non-Irrigating Bipolar Forceps

YASARGIL Style Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| GK788R         | 0.4           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 9 1/4        |
| GK789R         | 0.7           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 9 1/4        |
| GK790R         | 1.0           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 9 1/4        |
| GK791R         | 1.3           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 9 1/4        |
| GK793R         | 0.7           | BAYONET | INS        | NON-IRRIG  | 5 1/4                | 10           |



ELECTRO  
SURGERY






INCHES

Electrosurgery

Non-Irrigating Bipolar Forceps



YASARGIL Style Bipolar Forceps, Stainless Steel

| Catalog Number | Tip      | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|----------|---------------|---------|------------|------------|----------------------|--------------|
| GK777R         | 45° UP   | 0.7           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |
| GK780R         | 90° UP   | 0.4           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |
| GK781R         | 90° UP   | 0.7           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |
| GK785R         | 90° DOWN | 0.7           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |



YASARGIL Style Bipolar Forceps





| Catalog Number | Tip  | Tip Size (mm) | Working Length (in.) | Length (in.) |
|----------------|------|---------------|----------------------|--------------|
| GK831R         | RING | 2/3           | 4 1/2                | 9 1/4        |
| GK832R         | RING | 3/4           | 4 1/2                | 9 1/4        |
| GK833R         | RING | 4/5           | 4 1/2                | 9 1/4        |





Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| GK660R         | 0.6           | BAYONET | INS        | NON-IRRIG  | 1 1/2                | 6 1/4        |
| GK670R         | 0.6           | BAYONET | INS        | NON-IRRIG  | 3 1/8                | 7 7/8        |
| GK671R         | 0.9           | BAYONET | INS        | NON-IRRIG  | 3 1/8                | 7 7/8        |
| GK675R         | 0.6           | BAYONET | INS        | NON-IRRIG  | 4                    | 8 3/4        |
| GK677R         | 0.6           | BAYONET | INS        | NON-IRRIG  | 4 3/4                | 9 1/2        |
| GK679R         | 0.9           | BAYONET | INS        | NON-IRRIG  | 4 3/4                | 9 1/2        |





# Electrosurgery

## Non-Irrigating Bipolar Forceps

### Bayonet Forceps, Rose Gold, Insulated

■ Rose Gold plated tips reduce tissue adhesion



#### Non-Stick Rose Gold forceps

The Rose Gold forceps enable heat to be quickly discharged away from the tip.



#### Standard forceps

The heat cannot be discharged fast enough. This leads to overheating and sticking of tissue.

### Rose Gold Bipolar Forceps, Bayonet-Shaped

| Catalog Number | Tip Size (mm) | Style    | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|----------|------------|------------|----------------------|--------------|
| US403          | 0.5           | HARDY    | INS        | NON-IRRG   | 2                    | 6            |
| US404          | 0.7           | HARDY    | INS        | NON-IRRG   | 2                    | 6            |
| US405          | 1.0           | HARDY    | INS        | NON-IRRG   | 2                    | 6            |
| US406          | 1.5           | HARDY    | INS        | NON-IRRG   | 2                    | 6            |
| US407          | 2.0           | HARDY    | INS        | NON-IRRG   | 2                    | 6            |
| US605          | 0.5           | SCOVILLE | INS        | NON-IRRG   | 2 3/4                | 7 3/4        |
| US606          | 0.7           | SCOVILLE | INS        | NON-IRRG   | 2 3/4                | 7 3/4        |
| US607          | 1.0           | SCOVILLE | INS        | NON-IRRG   | 2 3/4                | 7 3/4        |
| US608          | 1.5           | SCOVILLE | INS        | NON-IRRG   | 2 3/4                | 7 3/4        |
| US609          | 2.0           | SCOVILLE | INS        | NON-IRRG   | 2 3/4                | 7 3/4        |
| US600          | 0.5           | CUSHING  | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US601          | 0.7           | CUSHING  | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US602          | 1.0           | CUSHING  | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US603          | 1.5           | CUSHING  | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US604          | 2.0           | CUSHING  | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US610          | 0.5           | HARDY    | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |
| US611          | 0.7           | HARDY    | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |
| US612          | 1.0           | HARDY    | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |
| US613          | 1.5           | HARDY    | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |
| US614          | 2.0           | HARDY    | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |
| US615          | 0.5           | HARDY    | INS        | NON-IRRG   | 5                    | 9            |
| US616          | 0.7           | HARDY    | INS        | NON-IRRG   | 5                    | 9            |
| US617          | 1.0           | HARDY    | INS        | NON-IRRG   | 5                    | 9            |
| US618          | 1.5           | HARDY    | INS        | NON-IRRG   | 5                    | 9            |
| US619          | 2.0           | HARDY    | INS        | NON-IRRG   | 5                    | 9            |
| US530          | 0.5           | HARDY    | INS        | NON-IRRG   | 6                    | 10           |
| US620          | 0.7           | HARDY    | INS        | NON-IRRG   | 6                    | 10           |
| US621          | 1.0           | HARDY    | INS        | NON-IRRG   | 6                    | 10           |
| US173          | 1.5           | HARDY    | INS        | NON-IRRG   | 6                    | 10           |
| US531          | 2.0           | HARDY    | INS        | NON-IRRG   | 6                    | 10           |

US173-US621

US403, US530, US600, US610, US615

US404, US601, US611, US616, US621

US405, US602, US612, US617

US173, US406, US603, US613, US618

US407, US531, US604, US614, US619

# Electrosurgery

## Non-Irrigating Bipolar Forceps

### Rose Gold Microform Bipolar Forceps

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US526          | 0.5           | BAYONET | INS        | NON-IRRIG  | 1 3/8                | 6            |
| US527          | 0.7           | BAYONET | INS        | NON-IRRIG  | 1 3/8                | 6            |
| US528          | 1.0           | BAYONET | INS        | NON-IRRIG  | 1 3/8                | 6            |
| US529          | 1.3           | BAYONET | INS        | NON-IRRIG  | 1 3/8                | 6            |
| US570          | 0.5           | BAYONET | INS        | NON-IRRIG  | 2 3/8                | 7            |
| US571          | 0.7           | BAYONET | INS        | NON-IRRIG  | 2 3/8                | 7            |
| US572          | 1.0           | BAYONET | INS        | NON-IRRIG  | 2 3/8                | 7            |
| US573          | 1.3           | BAYONET | INS        | NON-IRRIG  | 2 3/8                | 7            |
| US574          | 0.5           | BAYONET | INS        | NON-IRRIG  | 3                    | 7 3/4        |
| US575          | 0.7           | BAYONET | INS        | NON-IRRIG  | 3                    | 7 3/4        |
| US576          | 1.0           | BAYONET | INS        | NON-IRRIG  | 3                    | 7 3/4        |
| US577          | 1.3           | BAYONET | INS        | NON-IRRIG  | 3                    | 7 3/4        |
| US578          | 0.5           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |
| US579          | 0.7           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |
| US580          | 1.0           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |
| US581          | 1.3           | BAYONET | INS        | NON-IRRIG  | 3 3/4                | 8 1/2        |
| US582          | 0.5           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 9 1/4        |
| US583          | 0.7           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 9 1/4        |
| US584          | 1.0           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 9 1/4        |
| US585          | 1.3           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 9 1/4        |
| US586          | 0.5           | BAYONET | INS        | NON-IRRIG  | 5 1/2                | 10 1/4       |
| US587          | 0.7           | BAYONET | INS        | NON-IRRIG  | 5 1/2                | 10 1/4       |
| US588          | 1.0           | BAYONET | INS        | NON-IRRIG  | 5 1/2                | 10 1/4       |
| US589          | 1.3           | BAYONET | INS        | NON-IRRIG  | 5 1/2                | 10 1/4       |



US526-US589

Electrosurgery

Non-Irrigating Bipolar Forceps


Rose Gold Reverse Action Bipolar Forceps

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US547          | 0.7           | BAYONET | INS        | NON-IRRIG  | 6                    | 10           |
| US548          | 1.0           | BAYONET | INS        | NON-IRRIG  | 6                    | 10           |



Rose Gold Minimally Invasive (MIS) Forceps

| Catalog Number | Tip    | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|--------|---------------|---------|------------|------------|----------------------|--------------|
| US541          | STR    | 0.7           | BAYONET | INS        | NON-IRRIG  | 5 1/4                | 9 1/2        |
| US542          | STR    | 1.2           | BAYONET | INS        | NON-IRRIG  | 5 1/4                | 9 1/2        |
| US543          | STR    | 0.7           | BAYONET | INS        | NON-IRRIG  | 6 1/4                | 10 1/2       |
| US544          | STR    | 1.2           | BAYONET | INS        | NON-IRRIG  | 6 1/4                | 10 1/2       |
| US545          | 45° UP | 0.7           | BAYONET | INS        | NON-IRRIG  | 5 1/4                | 9 1/2        |
| US546          | 45° UP | 0.7           | BAYONET | INS        | NON-IRRIG  | 6 1/4                | 10 1/2       |



Rose Gold SEMKIN Bipolar Forceps

| Catalog Number | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|---------------|-------|------------|------------|--------------|
| US413          | 1.5           | STR   | INS        | NON-IRRIG  | 6            |



# Electrosurgery

## Non-Irrigating Bipolar Forceps

### SCOVILLE Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US534          | 1.0           | BAYONET | INS        | NON-IRRG   | 2 3/4                | 7 3/4        |
| US301          | 1.5           | BAYONET | INS        | NON-IRRG   | 2 3/4                | 7 3/4        |
| US300          | 1.5           | BAYONET | NON-INS    | NON-IRRG   | 2 3/4                | 7 3/4        |





### GERALD Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US342          | 1.0           | BAYONET | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US307          | 1.0           | BAYONET | NON-INS    | NON-IRRG   | 3 1/2                | 7 1/2        |



### CUSHING Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US308          | 0.7           | BAYONET | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US182          | 1.0           | BAYONET | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US183          | 1.5           | BAYONET | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US305          | 2.0           | BAYONET | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |



ELECTRO SURGERY






INCHES

Electrosurgery

Non-Irrigating Bipolar Forceps

HARDY Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US310          | 0.5           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 8 1/2        |
| US536          | 0.7           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 8 1/2        |
| US184          | 1.0           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 8 1/2        |
| US185          | 1.5           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 8 1/2        |
| US186          | 2.0           | BAYONET | INS        | NON-IRRIG  | 4 1/2                | 8 1/2        |
| US187          | 0.5           | BAYONET | INS        | NON-IRRIG  | 6                    | 10           |
| US188          | 1.0           | BAYONET | INS        | NON-IRRIG  | 6                    | 10           |
| US189          | 1.5           | BAYONET | INS        | NON-IRRIG  | 6                    | 10           |



HARDY Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US309          | 0.5           | BAYONET | NON-INS    | NON-IRRIG  | 4 1/2                | 8 1/2        |



ELECTRO SURGERY

# Electrosurgery

## Non-Irrigating Bipolar Forceps

### BARRÉ BipoJET® Nerve/Tissue Sparing Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Length (in.) |
|----------------|---------------|---------|------------|------------|--------------|
| GK685R         | 3.0           | BAYONET | INS        | NON-IRRIG  | 10 1/4       |



### Strong Tension Bipolar Forceps, Stainless Steel


| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US314          | 0.5           | BAYONET | NON-INS    | NON-IRRIG  | 2 3/8                | 6 1/4        |
| US315          | 0.7           | BAYONET | NON-INS    | NON-IRRIG  | 2 3/8                | 6 1/4        |
| US316          | 1.0           | BAYONET | NON-INS    | NON-IRRIG  | 2 3/8                | 6 1/4        |
| US317          | 1.3           | BAYONET | NON-INS    | NON-IRRIG  | 2 3/8                | 6 1/4        |
| US318          | 0.5           | BAYONET | NON-INS    | NON-IRRIG  | 3 1/4                | 7 1/8        |
| US319          | 0.7           | BAYONET | NON-INS    | NON-IRRIG  | 3 1/4                | 7 1/8        |
| US320          | 1.0           | BAYONET | NON-INS    | NON-IRRIG  | 3 1/4                | 7 1/8        |
| US321          | 1.3           | BAYONET | NON-INS    | NON-IRRIG  | 3 1/4                | 7 1/8        |
| US322          | 0.5           | BAYONET | NON-INS    | NON-IRRIG  | 4                    | 7 7/8        |
| US323          | 0.7           | BAYONET | NON-INS    | NON-IRRIG  | 4                    | 7 7/8        |
| US324          | 1.0           | BAYONET | NON-INS    | NON-IRRIG  | 4                    | 7 7/8        |
| US325          | 1.3           | BAYONET | NON-INS    | NON-IRRIG  | 4                    | 7 7/8        |



INCHES

Electrosurgery

Non-Irrigating Bipolar Forceps







CASPAR Style Bipolar Forceps:

- Developed for precision neurosurgery with a special focus on spine surgery.
- Parallel guiding handle for exact tip alignment.

CASPAR Style Bipolar Forceps, Stainless Steel







| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| GK899R         | 0.2           | BAYONET | INS        | NON-IRRIG  | 2 3/4                | 6 1/2        |
| GK900R         | 0.5           | BAYONET | INS        | NON-IRRIG  | 2 3/4                | 6 1/2        |
| GK910R         | 1.0           | BAYONET | INS        | NON-IRRIG  | 2 3/4                | 6 1/2        |
| GK920R         | 2.0           | BAYONET | INS        | NON-IRRIG  | 2 3/4                | 6 1/2        |
| GK929R         | 0.2           | BAYONET | INS        | NON-IRRIG  | 4                    | 7 3/4        |
| GK930R         | 0.5           | BAYONET | INS        | NON-IRRIG  | 4                    | 7 3/4        |
| GK940R         | 1.0           | BAYONET | INS        | NON-IRRIG  | 4                    | 7 3/4        |
| GK950R         | 2.0           | BAYONET | INS        | NON-IRRIG  | 4                    | 7 3/4        |
| GK959R         | 0.2           | BAYONET | INS        | NON-IRRIG  | 5                    | 8 3/4        |
| GK960R         | 0.5           | BAYONET | INS        | NON-IRRIG  | 5                    | 8 3/4        |
| GK970R         | 1.0           | BAYONET | INS        | NON-IRRIG  | 5                    | 8 3/4        |
| GK980R         | 2.0           | BAYONET | INS        | NON-IRRIG  | 5                    | 8 3/4        |





CASPAR Style Bipolar Forceps, Stainless Steel

| Catalog Number | Tip      | Shape   | Tip Size (mm) | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|----------|---------|---------------|------------|------------|----------------------|--------------|
| GK976R         | 90° UP   | BAYONET | 1.0           | INS        | NON-IRRIG  | 4 3/4                | 8 1/2        |
| GK978R         | 90° DOWN | BAYONET | 1.0           | INS        | NON-IRRIG  | 4 3/4                | 8 1/2        |
| GK972R         | 45° UP   | BAYONET | 1.0           | INS        | NON-IRRIG  | 5                    | 8 3/4        |
| GK974R         | 45° DOWN | BAYONET | 1.0           | INS        | NON-IRRIG  | 5                    | 8 3/4        |




Electrosurgery

Non-Irrigating Bipolar Forceps

SINTRAM Bipolar Forceps

| Catalog Number | Tip     | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------|---------------|---------|------------|------------|----------------------|--------------|
| GK644R         | STR     | 0.7 x 6       | BAYONET | INS        | NON-IRRIG  | 2 3/8                | 6 1/4        |
| GK645R         | STR     | 1.0 x 6       | BAYONET | INS        | NON-IRRIG  | 3 1/4                | 7 1/4        |
| GK646R         | STR     | 1.0 x 6       | BAYONET | INS        | NON-IRRIG  | 4                    | 7 7/8        |
| GK614R         | STR     | 0.7 x 23      | BAYONET | INS        | NON-IRRIG  | 2 3/8                | 6 1/4        |
| GK615R         | STR     | 1.0 x 23      | BAYONET | INS        | NON-IRRIG  | 3 1/4                | 7 1/4        |
| GK616R         | STR     | 1.0 x 23      | BAYONET | INS        | NON-IRRIG  | 4                    | 7 7/8        |
| GK604R         | STR     | 0.9 x 23      | STR     | INS        | NON-IRRIG  | -                    | 6 1/4        |
| GK605R         | STR     | 0.9 x 23      | STR     | INS        | NON-IRRIG  | -                    | 7 1/4        |
| GK606R         | STR     | 0.9 x 23      | STR     | INS        | NON-IRRIG  | -                    | 7 7/8        |
| GK624R         | ANG 35° | 0.9 x 23      | STR     | INS        | NON-IRRIG  | -                    | 6 1/4        |
| GK625R         | ANG 35° | 0.9 x 23      | STR     | INS        | NON-IRRIG  | -                    | 7 7/8        |



MALIS Style Titanium Bipolar Forceps

| Catalog Number | Tip      | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|----------|---------------|---------|------------|------------|----------------------|--------------|
| US331          | STR      | 0.5           | BAYONET | NON-INS    | NON-IRRG   | 3 3/4                | 8 3/4        |
| US332          | STR      | 1.0           | BAYONET | NON-INS    | NON-IRRG   | 3 3/4                | 8 3/4        |
| US330          | STR      | 1.5           | BAYONET | NON-INS    | NON-IRRG   | 3 3/4                | 8 3/4        |
| US333          | ANG-UP   | 0.5           | BAYONET | NON-INS    | NON-IRRG   | 3 3/4                | 8 3/4        |
| US175          | ANG-DOWN | 0.5           | BAYONET | NON-INS    | NON-IRRG   | 3 3/4                | 8 3/4        |



US330-US332



US330-US332



US175, US333



US175, US333

MALIS Style Rose Gold Bipolar Forceps

| Catalog Number | Tip      | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|----------|---------------|---------|------------|------------|----------------------|--------------|
| US161          | STR      | 0.5           | BAYONET | INS        | NON-IRRG   | 3 3/4                | 8 3/4        |
| US162          | STR      | 1.0           | BAYONET | INS        | NON-IRRG   | 3 3/4                | 8 3/4        |
| US163          | STR      | 1.5           | BAYONET | INS        | NON-IRRG   | 3 3/4                | 8 3/4        |
| US164          | ANG-UP   | 0.5           | BAYONET | INS        | NON-IRRG   | 3 3/4                | 8 3/4        |
| US165          | ANG-DOWN | 0.5           | BAYONET | INS        | NON-IRRG   | 3 3/4                | 8 3/4        |




US161-US165




US161



US162




US163



US164



US165



# Electrosurgery

## Non-Irrigating Bipolar Forceps

### Titanium Bipolar Forceps

| Catalog Number | Tip | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|-----|---------------|-------|------------|------------|--------------|
| US326          | STR | 0.5           | STR   | NON-INS    | NON-IRRG   | 6            |
| US328          | STR | 0.5           | STR   | NON-INS    | NON-IRRG   | 7 1/8        |
| US327          | CVD | 0.5           | STR   | NON-INS    | NON-IRRG   | 6            |
| US329          | CVD | 0.5           | STR   | NON-INS    | NON-IRRG   | 7 1/8        |



US326, US328



US327, US329

### Micro Bipolar Forceps, Stainless Steel

| Catalog Number | Tip | Tip Size (mm) | Shape | Irrigation | Insulation | Length (in.) |
|----------------|-----|---------------|-------|------------|------------|--------------|
| US345          | STR | 0.7           | STR   | NON-INS    | NON-IRRG   | 6            |
| US347          | STR | 0.7           | STR   | NON-INS    | NON-IRRG   | 7 1/8        |
| US346          | CVD | 0.7           | STR   | NON-INS    | NON-IRRG   | 6            |
| US348          | CVD | 0.7           | STR   | NON-INS    | NON-IRRG   | 7 1/8        |



US345, US347



US346, US348

### Classic Bipolar Forceps, Stainless Steel

| Catalog Number | Tip      | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|----------|---------------|-------|------------|------------|--------------|
| GK700R         | STR      | 0.6           | STR   | INS        | NON-IRRG   | 6 1/4        |
| GK701R         | STR      | 0.9           | STR   | INS        | NON-IRRG   | 6 1/4        |
| GK703R         | STR      | 0.9           | STR   | INS        | NON-IRRG   | 7 7/8        |
| GK711R         | ANG-DOWN | 0.9           | STR   | INS        | NON-IRRG   | 6 1/4        |
| GK713R         | ANG-DOWN | 0.9           | STR   | INS        | NON-IRRG   | 7 7/8        |



GK700R-GK703R



GK700R, GK701R, GK703R




GK711R-GK713R



GK711R, GK713R


# Electrosurgery

## Non-Irrigating Bipolar Forceps



Classic Bipolar Forceps, Stainless Steel

| Catalog Number | Tip      | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|----------|---------------|-------|------------|------------|--------------|
| GK690R         | STR      | 0.4           | STR   | INS        | NON-IRRIG  | 4 3/4        |
| GK720R         | STR      | 0.6           | STR   | INS        | NON-IRRIG  | 6 1/4        |
| GK722R         | STR      | 0.9           | STR   | INS        | NON-IRRIG  | 6 1/4        |
| GK730R         | STR      | 0.9           | STR   | INS        | NON-IRRIG  | 7 7/8        |
| GK723R         | ANG-DOWN | 0.9           | STR   | INS        | NON-IRRIG  | 6 1/4        |



Jewelers Bipolar Forceps, Stainless Steel

| Catalog Number | Tip | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|-----|---------------|-------|------------|------------|--------------|
| US337          | STR | 0.4           | STR   | INS        | NON-IRRIG  | 4            |
| US335          | STR | 0.4           | STR   | NON-INS    | NON-IRRIG  | 4            |
| US159          | CVD | 0.4           | STR   | INS        | NON-IRRIG  | 4            |
| US336          | CVD | 0.4           | STR   | NON-INS    | NON-IRRIG  | 4            |




TENZEL Bipolar Forceps, Stainless Steel

| Catalog Number | Tip | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|-----|---------------|-------|------------|------------|--------------|
| US178          | CVD | 0.4           | STR   | INS        | NON-IRRIG  | 4            |
| US338          | CVD | 0.4           | STR   | NON-INS    | NON-IRRIG  | 4            |



Electrosurgery

Non-Irrigating Bipolar Forceps



US339-US340

ADSON Bipolar Forceps, Stainless Steel

| Catalog Number | Tip | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|-----|---------------|-------|------------|------------|--------------|
| US340          | STR | 1.0           | STR   | INS        | NON-IRRIG  | 5            |
| US339          | STR | 1.0           | STR   | NON-INS    | NON-IRRIG  | 5            |



US179-US344

SEMKIN Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|---------------|-------|------------|------------|--------------|
| US179          | 0.7           | STR   | INS        | NON-IRRIG  | 5 1/2        |
| US311          | 0.7           | STR   | NON-INS    | NON-IRRIG  | 5 1/2        |
| US344          | 2.0           | STR   | INS        | NON-IRRIG  | 6            |
| US196          | 2.0           | STR   | NON-INS    | NON-IRRIG  | 6            |



US197-US304

CUSHING Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|---------------|-------|------------|------------|--------------|
| US197          | 0.7           | STR   | INS        | NON-IRRIG  | 6            |
| US302          | 0.7           | STR   | NON-INS    | NON-IRRIG  | 6            |
| US304          | 1.5           | STR   | INS        | NON-IRRIG  | 7            |
| US303          | 1.5           | STR   | NON-INS    | NON-IRRIG  | 7            |

# Electrosurgery

## Non-Irrigating Bipolar Forceps

### GERALD Bipolar Forceps, Stainless Steel


| Catalog Number | Tip | Shape | Insulation | Irrigation | Length (in.) |
|----------------|-----|-------|------------|------------|--------------|
| US341          | STR | STR   | INS        | NON-IRRG   | 7            |
| US306          | STR | STR   | NON-INS    | NON-IRRG   | 7            |




US306-US341

### Rose Gold Bipolar Forceps

| Catalog Number | Tip Size (mm) | Shape | Style    | Insulation | Irrigation | Length (in.) |
|----------------|---------------|-------|----------|------------|------------|--------------|
| US171          | 1.0           | STR   | ADSON    | INS        | NON-IRRG   | 5            |
| US172          | 0.4           | STR   | JEWELERS | INS        | NON-IRRG   | 4            |




US171



US171



US172



US172

### MALIS Style Titanium Bipolar Forceps

| Catalog Number | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|---------------|-------|------------|------------|--------------|
| US176          | 0.25          | STR   | NON-INS    | NON-IRRG   | 6            |
| US177          | 0.5           | STR   | NON-INS    | NON-IRRG   | 6            |



US176-US177




US176-US177

### Micro and Laboratory Forceps, Stainless Steel


| Catalog Number | Tip | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|-----|---------------|-------|------------|------------|--------------|
| GK596R         | STR | 0.3           | STR   | INS        | NON-IRRG   | 5 3/4        |
| GK595R         | STR | 0.6           | STR   | INS        | NON-IRRG   | 4 3/4        |
| GK597R         | CVD | 0.6           | STR   | INS        | NON-IRRG   | 4 1/4        |



GK595R-GK596R



GK595R



GK596R



GK597R



GK597R



# Electrosurgery

## Non-Irrigating Bipolar Forceps

### BiProtect® Ceramic Insulated Tip, Insulated


- Ceramic tip insulation helps ensure that heat concentrates on the inner and distal part of the forceps
  - Less collateral damage and reduced thermal spread
  - Partially insulated tip allows for precise and safe surface coagulation

YASARGIL BiProtect Ceramic Insulated Tip Forceps


| Catalog Number | Tip | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|-----|---------------|---------|------------|------------|----------------------|--------------|
| GL470R         | STR | 0.4           | BAYONET | INS        | NON-IRRG   | 3 1/2                | 7 3/4        |
| GL471R         | STR | 0.8           | BAYONET | INS        | NON-IRRG   | 3 1/2                | 7 3/4        |
| GL480R         | STR | 0.4           | BAYONET | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |
| GL481R         | STR | 0.8           | BAYONET | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |



GL470R-GL481R



GL470R, GL480R



GL471R, GL481R

Transnasal Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|-------|------------|------------|----------------------|--------------|
| GK719R         | 1.0           | ANG   | INS        | NON-IRRG   | 5 1/8                | 8 1/2        |



GK719R



GK719R



MIN Instruments, Pivot-Point Bipolar Forceps, Transnasal

| Catalog Number | Tip | Tip Size (mm) | Shape   | Working Length (in.) | Length (in.) |
|----------------|-----|---------------|---------|----------------------|--------------|
| GK822R         | STR | 0.7           | BAYONET | 3 3/4                | 8 1/2        |
| GK824R         | STR | 1.0           | BAYONET | 3 3/4                | 8 1/2        |
| GK826R         | STR | 0.7           | BAYONET | 5 1/4                | 10           |
| GK828R         | STR | 1.0           | BAYONET | 5 1/4                | 10           |
| GK823R         | CVD | 0.7           | BAYONET | 3 3/4                | 8 1/2        |
| GK825R         | CVD | 1.0           | BAYONET | 3 3/4                | 8 1/2        |
| GK827R         | CVD | 0.7           | BAYONET | 5 1/4                | 10           |
| GK829R         | CVD | 1.0           | BAYONET | 5 1/4                | 10           |



Special pin between the branches opens the tip of the forceps by additional compression of the handle, allowing secure coagulation in narrow and deep seated surgical fields.


95 mm, 3 3/4"

GK822R-GK824R



135 mm, 5 1/4"

GK826R-GK828R




95 mm, 3 3/4"

GK823R-GK825R



135 mm, 5 1/4"

GK827R-GK829R




1/1

GK822R, GK823R, GK826R, GK827R

1/1

GK824R, GK825R, GK828R, GK829R




Electrosurgery

Non-Irrigating Bipolar Forceps

LANDOLT Bipolar Forceps, Stainless Steel

| Catalog Number | Tip      | Insulation | Irrigation | Description                           | Working Length (in.) |
|----------------|----------|------------|------------|---------------------------------------|----------------------|
| GK560R         | ANG 90°  | INS        | NON-IRRG   | COAGULATION FORCEPS F/ HYPOPHYSECTOMY | 6                    |
| GK561R         | 90°      | -          | -          | SPARE INNER PART F/ GK560R            |                      |
| GK580R         | ANG 120° | INS        | NON-IRRG   | COAGULATION FORCEPS F/ HYPOPHYSECTOMY | 6                    |
| GK581R         | 120°     | -          | -          | SPARE INNER PART F/ GK580R            |                      |



ELECTRO  
SURGERY


# Electrosurgery

## Irrigating Bipolar Forceps

■ Irrigating coagulation forceps facilitate gentle coagulation without tissue adhesion.

### Rose Gold Irrigating Bipolar Forceps


| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US622          | 0.5           | BAYONET | INS        | IRRG       | 3 1/2                | 7 1/2        |
| US623          | 0.7           | BAYONET | INS        | IRRG       | 3 1/2                | 7 1/2        |
| US624          | 1.0           | BAYONET | INS        | IRRG       | 3 1/2                | 7 1/2        |
| US625          | 1.5           | BAYONET | INS        | IRRG       | 3 1/2                | 7 1/2        |
| US626          | 2.0           | BAYONET | INS        | IRRG       | 3 1/2                | 7 1/2        |
| US627          | 0.5           | BAYONET | INS        | IRRG       | 3                    | 7 3/4        |
| US628          | 0.7           | BAYONET | INS        | IRRG       | 3                    | 7 3/4        |
| US629          | 1.0           | BAYONET | INS        | IRRG       | 3                    | 7 3/4        |
| US630          | 1.5           | BAYONET | INS        | IRRG       | 3                    | 7 3/4        |
| US631          | 2.0           | BAYONET | INS        | IRRG       | 3                    | 7 3/4        |
| US632          | 0.5           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US633          | 0.7           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US634          | 1.0           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US635          | 1.5           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US636          | 2.0           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US637          | 0.5           | BAYONET | INS        | IRRG       | 5                    | 9            |
| US638          | 0.7           | BAYONET | INS        | IRRG       | 5                    | 9            |
| US639          | 1.0           | BAYONET | INS        | IRRG       | 5                    | 9            |
| US640          | 1.5           | BAYONET | INS        | IRRG       | 5                    | 9            |
| US641          | 2.0           | BAYONET | INS        | IRRG       | 5                    | 9            |
| US532          | 0.5           | BAYONET | INS        | IRRG       | 6                    | 10           |
| US642          | 0.7           | BAYONET | INS        | IRRG       | 6                    | 10           |
| US643          | 1.0           | BAYONET | INS        | IRRG       | 6                    | 10           |
| US174          | 1.5           | BAYONET | INS        | IRRG       | 6                    | 10           |
| US533          | 2.0           | BAYONET | INS        | IRRG       | 6                    | 10           |



US174-US533

### CUSHING Irrigating Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US369          | 0.7           | BAYONET | INS        | IRRG       | 3 1/2                | 7 1/2        |
| US190          | 1.0           | BAYONET | INS        | IRRG       | 3 1/2                | 7 1/2        |
| US191          | 1.5           | BAYONET | INS        | IRRG       | 3 1/2                | 7 1/2        |
| US370          | 2.0           | BAYONET | INS        | IRRG       | 3 1/2                | 7 1/2        |



US190-US370



US190-US370

# Electrosurgery

## Irrigating Bipolar Forceps

### HARDY Irrigating Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US366          | 0.5           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US382          | 1.0           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US385          | 1.5           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US192          | 2.0           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US193          | 0.5           | BAYONET | INS        | IRRG       | 6                    | 10           |
| US194          | 1.0           | BAYONET | INS        | IRRG       | 6                    | 10           |
| US195          | 1.5           | BAYONET | INS        | IRRG       | 6                    | 10           |
| US558          | 0.5           | BAYONET | NON-INS    | IRRG       | 4 1/2                | 8 1/2        |

US192-US195, US382, US385, US558

US192-US195, US382, US385, US558

### SCOVILLE Irrigating Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US367          | 1.5           | BAYONET | INS        | IRRG       | 2 3/4                | 7 3/4        |
| US556          | 1.5           | BAYONET | NON-INS    | IRRG       | 2 3/4                | 7 3/4        |

US367, US556

### GERALD Irrigating Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US368          | 1.0           | BAYONET | INS        | IRRG       | 3 1/2                | 7 1/2        |
| US557          | 1.0           | BAYONET | NON-INS    | IRRG       | 3 1/2                | 7 1/2        |


US368, US557

# Electrosurgery

## Irrigating Bipolar Forceps

### MALIS Style Titanium Irrigating Bipolar Forceps

| Catalog Number | Tip      | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|----------|---------------|---------|------------|------------|----------------------|--------------|
| US360          | STR      | 1.5           | BAYONET | INS        | IRRG       | 3 3/4                | 8 3/4        |
| US361          | STR      | 0.5           | BAYONET | INS        | IRRG       | 3 3/4                | 8 3/4        |
| US362          | STR      | 1.0           | BAYONET | INS        | IRRG       | 3 3/4                | 8 3/4        |
| US363          | ANG-UP   | 0.5           | BAYONET | INS        | IRRG       | 3 3/4                | 8 3/4        |
| US372          | ANG-DOWN | 0.5           | BAYONET | INS        | IRRG       | 3 3/4                | 8 3/4        |



US360-US362



US360-US362



US363, US372



US363, US372

### MALIS Style Rose Gold Irrigating Bipolar Forceps

| Catalog Number | Tip      | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|----------|---------------|---------|------------|------------|----------------------|--------------|
| US166          | STR      | 0.5           | BAYONET | INS        | IRRG       | 3 3/4                | 8 3/4        |
| US167          | STR      | 1.0           | BAYONET | INS        | IRRG       | 3 3/4                | 8 3/4        |
| US168          | STR      | 1.5           | BAYONET | INS        | IRRG       | 3 3/4                | 8 3/4        |
| US169          | ANG-UP   | 0.5           | BAYONET | INS        | IRRG       | 3 3/4                | 8 3/4        |
| US170          | ANG-DOWN | 0.5           | BAYONET | INS        | IRRG       | 3 3/4                | 8 3/4        |



US166-US168



US169-US170









INCHES

Electrosurgery

Irrigating Bipolar Forceps

Irrigating Bipolar Forceps, Stainless Steel

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| GK838R         | 0.7           | BAYONET | INS        | IRRG       | 2                    | 6 1/2        |
| GK840R         | 0.4           | BAYONET | INS        | IRRG       | 2 3/4                | 7            |
| GK841R         | 0.7           | BAYONET | INS        | IRRG       | 2 3/4                | 7            |
| GK842R         | 1.3           | BAYONET | INS        | IRRG       | 2 3/4                | 7            |
| GK843R         | 0.4           | BAYONET | INS        | IRRG       | 3 1/2                | 8            |
| GK844R         | 0.7           | BAYONET | INS        | IRRG       | 3 1/2                | 8            |
| GK845R         | 1.3           | BAYONET | INS        | IRRG       | 3 1/2                | 8            |
| GK846R         | 0.4           | BAYONET | INS        | IRRG       | 4 1/2                | 9            |
| GK847R         | 0.7           | BAYONET | INS        | IRRG       | 4 1/2                | 9            |
| GK848R         | 1.3           | BAYONET | INS        | IRRG       | 4 1/2                | 9            |



# Electrosurgery

## Disposable Bipolar Forceps

### Stainless Steel Disposable Bipolar Forceps, Box of 10


| Catalog Number | Tip Size (mm) | Shape   | Style    | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|----------|------------|------------|----------------------|--------------|
| US153SU        | 1.5           | BAYONET | SCOVILLE | INS        | NON-IRRIG  | 2 3/4                | 7 3/4        |
| US567SU        | 2.0           | BAYONET | CUSHING  | INS        | NON-IRRIG  | 2 3/4                | 7 3/4        |
| US562SU        | 0.7           | BAYONET | CUSHING  | INS        | NON-IRRIG  | 3 1/2                | 7 1/2        |
| US560SU        | 1.0           | BAYONET | GERALD   | INS        | NON-IRRIG  | 3 1/2                | 7 1/2        |
| US154SU        | 0.5           | BAYONET | HARDY    | INS        | NON-IRRIG  | 4 1/2                | 8 1/2        |
| US155SU        | 1.0           | BAYONET | HARDY    | INS        | NON-IRRIG  | 4 1/2                | 8 1/2        |
| US156SU        | 1.5           | BAYONET | HARDY    | INS        | NON-IRRIG  | 4 1/2                | 8 1/2        |
| US157SU        | 2.0           | BAYONET | HARDY    | INS        | NON-IRRIG  | 4 1/2                | 8 1/2        |



US153SU-US567SU



US153SU-US567SU



US158SU



US158SU

### ADSON Disposable Bipolar Forceps, Stainless Steel, Box of 10

| Catalog Number | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|---------------|-------|------------|------------|--------------|
| US158SU        | 1.0           | STR   | INS        | NON-IRRIG  | 5            |

### Disposable Jewelers Bipolar Forceps, Stainless Steel, Box of 10

| Catalog Number | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|---------------|-------|------------|------------|--------------|
| US568SU        | 0.4           | STR   | NON-INS    | NON-IRRIG  | 4            |
| US569SU        | 0.4           | STR   | INS        | NON-IRRIG  | 4            |



US568SU-US569SU



Electrosurgery  
Disposable Bipolar Forceps

Rose Gold Disposable Jewelers Bipolar Forceps, Box of 5

| Catalog Number | Tip Size (mm) | Shape | Insulation | Irrigation | Length (in.) |
|----------------|---------------|-------|------------|------------|--------------|
| US408SU        | 0.4           | STR   | INS        | NON-IRRG   | 4            |
| US409SU        | 1.0           | STR   | INS        | NON-IRRG   | 5            |



US408SU-US409SU

Rose Gold Disposable Bipolar Forceps, Box of 5

| Catalog Number | Tip Size (mm) | Shape   | Style    | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|----------|------------|------------|----------------------|--------------|
| US565SU        | 0.7           | STR     | SEMKin   | INS        | NON-IRRG   |                      | 5 1/2        |
| US559SU        | 1.5           | BAYONET | SCOVILLE | INS        | NON-IRRG   | 2 3/4                | 7 3/4        |
| US201SU        | 0.5           | BAYONET | CUSHING  | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US561SU        | 0.7           | BAYONET | CUSHING  | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US202SU        | 1.0           | BAYONET | CUSHING  | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US203SU        | 1.5           | BAYONET | CUSHING  | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US566SU        | 2.0           | BAYONET | CUSHING  | INS        | NON-IRRG   | 3 1/2                | 7 1/2        |
| US204SU        | 0.5           | BAYONET | HARDY    | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |
| US563SU        | 0.7           | BAYONET | HARDY    | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |
| US205SU        | 1.0           | BAYONET | HARDY    | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |
| US206SU        | 1.5           | BAYONET | HARDY    | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |
| US564SU        | 2.0           | BAYONET | HARDY    | INS        | NON-IRRG   | 4 1/2                | 8 1/2        |



US201SU-US564SU



US565SU



# Electrosurgery

## Disposable Bipolar Forceps

### Rose Gold Disposable Irrigating Bipolar Forceps, Box of 5


| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US410SU        | 0.5           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US411SU        | 1.0           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US412SU        | 1.5           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |



US410SU-US412SU

### Stainless Steel Disposable Irrigating Bipolar Forceps, Box of 5

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Working Length (in.) | Length (in.) |
|----------------|---------------|---------|------------|------------|----------------------|--------------|
| US523SU        | 0.5           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US524SU        | 1.0           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |
| US525SU        | 1.5           | BAYONET | INS        | IRRG       | 4 1/2                | 8 1/2        |



US523SU-US525SU

### Specialty Disposable Non-Stick Bipolar Forceps

| Catalog Number | Tip Size (mm) | Shape   | Insulation | Irrigation | Length (in.) |
|----------------|---------------|---------|------------|------------|--------------|
| US501SU        | 0.5           | BAYONET | INS        | NON-IRRG   | 7            |
| US502SU        | 1.0           | BAYONET | INS        | NON-IRRG   | 7            |
| US503SU        | 1.5           | BAYONET | INS        | NON-IRRG   | 7            |
| US504SU        | 0.5           | BAYONET | INS        | NON-IRRG   | 8            |
| US505SU        | 1.0           | BAYONET | INS        | NON-IRRG   | 8            |
| US506SU        | 1.5           | BAYONET | INS        | NON-IRRG   | 8            |
| US507SU        | 0.5           | BAYONET | INS        | NON-IRRG   | 9            |
| US508SU        | 1.0           | BAYONET | INS        | NON-IRRG   | 9            |
| US509SU        | 1.5           | BAYONET | INS        | NON-IRRG   | 9            |



US501SU-US509SU



US501SU-US509SU





INCHES

# Electrosurgery


## Bipolar Connecting Cords

Bipolar Cords for BipoJET®


| Catalog Number | New Cord                                                                           | Pin Distance | Compatible Generators                 |
|----------------|------------------------------------------------------------------------------------|--------------|---------------------------------------|
| GN145          |   | 22 MM        | EMC<br>DOLLEX<br>ERBE Int'l<br>MARTIN |
| GN085          |   | 29 MM        | AESCLAP<br>VALLEYLAB<br>ERBE          |
| GN146          |  | 38 MM        | CODMAN CMC II / III<br>ERBE<br>MARTIN |



GK195



GN219



GN224

LANDOLT Bipolar Cord and Adapters

| Catalog Number | Pin Distance (mm) | Description                                |
|----------------|-------------------|--------------------------------------------|
| GK195          | 28.6              | REUSABLE CORD F/ INSULATED LANDOLT FORCEPS |
| GN219          | 32                | ADAPTER F/ CONMED                          |
| GN224          | 38                | ADAPTER F/ CODMAN                          |

# Electrosurgery

## Bipolar Connecting Cords




Dual Bipolar Cord, Sterile Disposable

| Catalog Number | Packaging | Length (ft) |
|----------------|-----------|-------------|
| US349SP        | 10/BOX    | 12          |




Round Bipolar Cord, Reusable

| Catalog Number | Material | Length (ft) |
|----------------|----------|-------------|
| US351          | SILICONE | 12          |



Dual Bipolar Cord, Reusable


| Catalog Number | Length (ft) |
|----------------|-------------|
| US359          | 12          |



INCHES

Electrosurgery

Bipolar Connecting Cords



| Combined Bipolar and Irrigation Cord |          |           |                              |
|--------------------------------------|----------|-----------|------------------------------|
| Catalog Number                       | Pin Type | Packaging | Compatible w/                |
| US352                                | DUAL     | BOX OF 10 | -                            |
| US353                                | ROUND    | BOX OF 10 | -                            |
| US355                                | DUAL     | BOX OF 10 | MALIS 1000 IRRIGATION SYSTEM |
| US356                                | ROUND    | BOX OF 10 | MALIS 1000 IRRIGATION SYSTEM |