



CGK Surgical

Product Catalog





General Surgery


Scissors	2 - 52
Needle Holders and Cases	53 - 76
Ligature Needles	77 - 78
Forceps and Clamps	79 - 189
Skin Hooks	190 - 192
Retractors	193 - 254
Universal Ring Retractor System	255 - 271
Self Retaining Retractors	272 - 283
Rectal Specula Instruments	284 - 291
Probes, Dilators and Scoops	292 - 299
Scalpel Blades and Handles	300 - 307
Amputation and Autopsy Instruments	308 - 310
Suction Instruments	311 - 322
Rulers, Calipers and Accessories	323 - 324
Metal Bowls and Cups	325

General Surgery


Scissors




BC374R-BC474R



BC374R



BC464R




BC474R

DEAVER Suture Scissors


Catalog Number	Tip	Blade	Length (in.)
BC374R	S/B	STR	5 3/4
BC474R	S/B	CVD	5 3/4
BC464R	B/B	CVD	5 3/4



BC130R-BC131R



BC130R




BC131R

Suture Scissors


Catalog Number	Tip	Blade	Length (in.)
BC130R	S/S	STR	4 3/4
BC341R	S/S	STR	4 1/2
BC343R	S/S	STR	5 1/8
BC344R	S/S	STR	5 3/4
BC345R	S/S	STR	6
BC346R	S/S	STR	6 1/2
BC131R	S/S	CVD	4 3/4




BC341R-BC346R



BC341R




BC343R



BC344R




BC345R




BC346R




BC441R-BC445R



BC441R



BC444R




BC443R




BC445R

Suture Scissors


Catalog Number	Tip	Blade	Length (in.)
BC441R	S/S	CVD	4 1/2
BC443R	S/S	CVD	5 1/8
BC444R	S/S	CVD	5 3/4
BC445R	S/S	CVD	6



BC303R-BC304R



BC303R




BC304R

Suture Scissors

Catalog Number	Tip	Blade	Length (in.)
BC303R	S/B	DEL, STR	5 1/8
BC304R	S/B	DEL, STR	5 3/4



BC496R




BC496R

Scissors


Catalog Number	Tip	Blade	Length (in.)
BC496R	S/B	STR/1-SERR	5 1/8

General Surgery


Scissors




BC320R-BC330R




BC320R



BC321R




BC323R




BC324R




BC325R




BC326R




BC327R




BC328R




BC330R



BC403R-BC404R



BC403R



BC404R

Suture Scissors


Catalog Number	Tip	Blade	Length (in.)
BC320R	S/B	STR	4 1/8
BC321R	S/B	STR	4 1/2
BC323R	S/B	STR	5 1/8
BC324R	S/B	STR	5 3/4
BC325R	S/B	STR	6
BC326R	S/B	STR	6 1/2
BC327R	S/B	STR	6 7/8
BC328R	S/B	STR	7 1/4
BC330R	S/B	STR	7 7/8
BC403R	S/B	CVD	5 1/8
BC404R	S/B	CVD	5 3/4




BC421R-BC428R

Suture Scissors


Catalog Number	Tip	Blade	Length (in.)
BC421R	S/B	CVD	4 1/2
BC423R	S/B	CVD	5 1/8
BC424R	S/B	CVD	5 3/4
BC425R	S/B	CVD	6
BC426R	S/B	CVD	6 1/2
BC428R	S/B	CVD	7 1/4




BC421R



BC423R




BC424R



BC425R



BC426R




BC428R



BC135R

Suture Scissors


Catalog Number	Tip	Blade	Length (in.)
BC135R	B/B	CVD, DEL	4 3/4




BC135R

General Surgery


Scissors




BC313R-BC317R




BC313R




BC314R



BC315R



BC316R




BC317R




BC411R-BC417R




BC411R




BC413R




BC414R



BC415R



BC416R




BC417R

COOPER Suture Scissors

Catalog Number	Tip	Blade	Length (in.)
BC313R	B/B	STR	5 1/8
BC314R	B/B	STR	5 3/4
BC315R	B/B	STR	6
BC316R	B/B	STR	6 1/2
BC317R	B/B	STR	6 7/8
BC411R	B/B	CVD	4 1/2
BC413R	B/B	CVD	5 1/8
BC414R	B/B	CVD	5 3/4
BC415R	B/B	CVD	6
BC416R	B/B	CVD	6 1/2
BC417R	B/B	CVD	6 7/8


LITTAUER Ligature Scissors

Catalog Number	Tip	Blade	Length (in.)
BC803R	B/1-NOTCH	STR	5 1/4




NORTHBENT Ligature Scissors

Catalog Number	Tip	Blade	Length (in.)
BC805R	B/1-NOTCH	LAT-CVD	4 3/4




SPENCER Suture Scissors

Catalog Number	Tip	Blade	Length (in.)
BC801R-BC802R	B/1-NOTCH	STR	3 1/2
BC802R	B/1-NOTCH	STR	4 1/2



KNOWLES Bandage Scissors

Catalog Number	Tip	Blade	Length (in.)
BC846R	1-PROBE-POINT	STR	5 3/4
BC847R	1-PROBE-POINT	LAT-ANG	5 1/2



General Surgery


Scissors




BC848R-BC863R

LISTER Bandage Scissors


Catalog Number	Blade	Length (in.)
BC848R	LONG	6 1/8
BC849R	LONG	7
BC860R	STANDARD	4 1/2
BC861R	STANDARD	5 1/2
BC862R	STANDARD	7
BC863R	STANDARD	7 7/8




BC848R 1/1




BC849R 1/1




BC860R 1/1




BC861R 1/1



BC862R 1/1



BC863R 1/1



MG379R



MG379R 1/1

LISTER Bandage Scissors, Gold Handle

Catalog Number	Length (in.)
MG379R	5 1/2





LISTER Bandage Scissors, W/ One Large Ring

Catalog Number	Length (in.)
MB223	8



LISTER Bandage Scissors, Serrated

Catalog Number	Length (in.)
MB146R	3 1/2
BC876R	6 1/4



General Surgery

Scissors

BC890R

1 1/2

BC890R

1 1/1

Nurse's Bandage Scissors

Catalog Number	Tip	Blade	Handle Color	Length (in.)
BC890R	B/B	STR	RED	5 3/4

BC885R-BC979R

1 1/2

BC885R-BC979R

1 1/1

LISTER Universal Bandage Scissors, Serrated

Catalog Number	Handle Color	Length (in.)
BC885R	BLACK	5 3/4
BC970R	RED	5 3/4
BC971R	BLUE	5 3/4
BC973R	YELLOW	5 3/4
BC975R	ORANGE	5 3/4
BC977R	PINK	5 3/4
BC978R	LIGHT GREEN	5 3/4
BC979R	VIOLET	5 3/4

BC881R-BC989R

1 1/2

BC881R-BC989R


1 1/1

Universal Bandage Scissors, W/ One Large Ring, Serrated

Catalog Number	Handle Color	Length (in.)
BC881R	BLACK	7
BC980R	RED	7
BC981R	BLUE	7
BC983R	YELLOW	7
BC985R	ORANGE	7
BC987R	PINK	7
BC988R	LIGHT GREEN	7
BC989R	VIOLET	7




BC580R-BC581R



BC580R

STILLE Dissecting Scissors, Round Blades


Catalog Number	Tip	Blade	Length (in.)
BC580R	B/B	STR	6
BC581R	B/B	CVD	6




BC581R



BC584R-BC589R




BC584R




BC585R

MAYO Dissecting Scissors, Round Blade


Catalog Number	Tip	Blade	Length (in.)
BC584R	B/B	STR	6 1/2
BC585R	B/B	STR	7 7/8
BC587R	B/B	CVD	6 1/2
BC588R	B/B	CVD	7 7/8
BC589R	B/B	CVD	8 1/4



BC587R




BC588R




BC589R

General Surgery


Scissors



BC560R-BC561R




BC560R




BC561R

MAYO-HARRINGTON Dissecting Scissors


Catalog Number	Tip	Blade	Length (in.)
BC560R	B/B	STR	9
BC561R	B/B	CVD	9



BC567R-BC568R




BC568R




BC567R

MAYO-HARRINGTON Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC568R	B/B	STR	11
BC567R	B/B	CVD	11



MB148R



MB148R

MAYO Dissecting Scissors, W/ One Large Ring


Catalog Number	Tip	Blade	Length (in.)
MB148R	B/B	HEAVY	14

MAYO Dissecting Scissors, Beveled
Blade


Catalog Number	Tip	Blade	Length (in.)
BC544R-BC557R	B/B	STR	5 1/2
BC545R	B/B	STR	6 1/8
BC547R	B/B	STR	6 3/4
BC562R	B/B	STR	7 1/2
BC564R	B/B	STR	8 1/2
BC554R	B/B	CVD	5 1/2
BC555R	B/B	CVD	6 1/8
BC557R	B/B	CVD	6 3/4
BC563R	B/B	CVD	7 1/2
BC565R	B/B	CVD	8 1/2




BC544R-BC557R




BC544R




BC545R




BC547R



BC554R




BC555R




BC557R




BC562R-BC565R




BC562R



BC564R



BC563R



BC565R

General Surgery

Scissors



BC466R-BC467R

1/2

1/1

BC466R

1/1

BC467R

MAYO Operating Scissors

Catalog Number	Tip	Blade	Length (in.)
BC466R	S/S	STR	5 1/2
BC467R	S/S	CVD	5 1/2

BC594R

1/2

1/1

BC594R

MAYO Dissecting Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC594R	B/B	STR	5 1/2

BC804R

1/2

1/1

BC804R

NEW'S Suture Scissors


Catalog Number	Tip	Blade	Length (in.)
BC804R	B/B	ANG	5 1/2

Suture Scissors

Catalog Number	Tip	Blade	Length (in.)
BC806R	B/B	STR	4 1/2




BC806R




BC806R

LEXER Dissecting Scissors, Delicate


Catalog Number	Tip	Blade	Length (in.)
BC590R	B/B	STR	6 1/2
BC591R	B/B	CVD	6 1/2



BC590R-BC591R




BC590R




BC591R

MAYO-LEXER Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC592R	B/B	CVD	6 1/2




BC592R




BC592R

MAYO-SIMS Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
MB149R	B/B	STR	7 3/4




MB149R




MB149R

General Surgery


Scissors



BC566R-BC569R



BC566R




BC569R

MAYO–NOBLE Dissecting Scissors, Broad Blade


Catalog Number	Tip	Blade	Length (in.)
BC566R	B/B	STR	6 3/4
BC569R	B/B	CVD	6 3/4




BC576R-BC579R




BC576R



BC578R



BC577R




BC579R

MAYO–STILLE Dissecting Scissors, Beveled Blade


Catalog Number	Tip	Blade	Length (in.)
BC576R	B/B	STR	5 1/2
BC578R	B/B	STR	6 1/2
BC577R	B/B	CVD	5 1/2
BC579R	B/B	CVD	6 1/2




BC570R-BC575R




BC570R



BC572R



BC573R







BC575R

MAYO–STILLE Dissecting Scissors, Flat Blade

Catalog Number	Tip	Blade	Length (in.)
BC570R	B/B	STR	5 1/2
BC572R	B/B	STR	6 3/4
BC573R	B/B	CVD	5 1/2
BC575R	B/B	CVD	6 3/4



TOENNIS Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC633R	B/B	X-DEL, STR	7
BC634R	B/B	DEL, STR	7
BC635R	B/B	DEL, CVD	7






Nerve Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC628R	S/S	CVD	6



TOENNIS-ADSON Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC630R	X-DEL	CVD	6 7/8
BC631R	DEL	CVD	6 7/8




General Surgery

Scissors




BC601R-BC607R



BC601R

BC602R

BC604R



BC605R

BC606R


BC607R

METZENBAUM Dissecting Scissors


Catalog Number	Tip	Blade	Length (in.)
BC601R	B/B	STR	5 3/4
BC602R	B/B	STR	7
BC604R	B/B	STR	7 7/8
BC605R	B/B	CVD	5 3/4
BC606R	B/B	CVD	7
BC607R	B/B	CVD	7 7/8



BC600R-BC603R




BC600R




BC603R

Baby METZENBAUM Dissecting Scissors, Delicate


Catalog Number	Tip	Blade	Length (in.)
BC600R	B/B	STR	5 3/4
BC603R	B/B	CVD	5 3/4



BC022R-BC023R




BC022R



BC023R

Baby METZENBAUM Dissecting Scissors


Catalog Number	Tip	Blade	Length (in.)
BC022R	B/B	STR	5 3/4
BC023R	B/B	CVD	5 3/4



MA464R-MA465R



MA465R



MA465R



METZENBAUM Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
MA464R	B/B	CVD	5 1/2
MA465R	B/B	CVD	8


NELSON-METZENBAUM Dissecting
Scissors

Catalog Number	Tip	Blade	Length (in.)
BC611R	B/B	STR	9
BC612R	B/B	STR	9 3/4
BC613R	B/B	STR	11 1/4
BC614R	B/B	STR	12
BC615R	B/B	CVD	9
BC616R	B/B	CVD	9 3/4
BC618R	B/B	CVD	11 1/4
BC619R	B/B	CVD	12
MB926R	B/B	CVD	16


BC611R-BC619R, MB926R




BC611R




BC612R




BC613R




BC614R




BC615R




BC616R



BC618R



BC619R



MB926R



General Surgery

Scissors





NELSON-METZENBAUM Dissecting
Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC625R	B/B	CVD	13 3/4




WELLER Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC620R	B/B	HEAVY-CVD	11




Abdominal Scissors

Catalog Number	Tip	Blade	Length (in.)
BC774R	B/B	STR	7



Rectal Scissors

Catalog Number	Tip	Blade	Length (in.)
BC617R	B/B	HEAVY-CVD	12 3/4




General Surgery

Scissors


LLOYD-DAVIES Rectal Scissors

Catalog Number	Tip	Blade	Length (in.)
BC621R	B/B	CVD	10 5/8




STELZNER Rectal Scissors

Catalog Number	Tip	Blade	Length (in.)
BC622R	B/B	STR	12 1/4
BC623R	B/B	CVD	12 1/4




MUELLER Rectal Scissors

Catalog Number	Tip	Blade	Length (in.)
BC627R	B/B	CVD	12 3/4




HARVEY Wire Scissors

Catalog Number	Tip	Blade	Length (in.)
BC815R	B/B	SERR	5










Wire Cutting Scissors

Catalog Number	For Wire Size	Length (in.)
DP551R	SOFT, UP TO 1.2MM AND HARD, UP TO 0.6 MM	4 3/4



General Surgery

Delicate Scissors





Iris Scissors

Catalog Number	Tip	Blade	Length (in.)
BC100R	S/S	STR	3 1/2
BC110R	S/S	STR	4 1/4
BC106R	S/S	STR	4 1/2
BC064R	S/S	STR	4 3/4
BC101R	S/S	CVD	3 1/2
BC111R	S/S	CVD	4 1/4
BC107R	S/S	CVD	4 1/2
BC065R	S/S	CVD	4 3/4
BC112R	S/B	STR	4 1/4
BC113R	S/B	CVD	4 1/4
BC117R	S/S	LAT-ANG	4 1/4




REYNOLDS-JAMESON Tenotomy
Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
BC015R	B/B	CVD	SQUARE	5 1/2





Strabismus Scissors

Catalog Number	Tip	Blade	Length (in.)
BC164R	B/B	STR	4 1/2
BC165R	B/B	CVD	4 1/2



General Surgery

Delicate Scissors





RAGNELL Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC026R	B/B	STR	5 1/4
BC027R	B/B	CVD	5 1/4
BC029R	B/B	CVD	6



CHADWICK Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC153R	S/S	CVD	4 1/2






LA GRANGE Scissors

Catalog Number	Tip	Blade	Length (in.)
BC154R	S/S	1-SERR	4 1/4




KELLY Scissors

Catalog Number	Tip	Blade	Length (in.)
BC160R	S/S	STR	6 7/8
BC161R	S/S	CVD	6 7/8



General Surgery

Delicate Scissors



ROCHESTER-FERGUSON Scissors

Catalog Number	Tip	Blade	Length (in.)
MG220R	B/B	ANG FLAT	7

DUROTIP Tungsten Carbide Scissors


- Excellent cutting accuracy
- Long working life
- Gold-plated rings

DUROTIP TC Suture Scissors

Catalog Number	Tip	Blade	Length (in.)
BC228R	S/S	STR	5 3/4
BC224R	S/B	STR	5 3/4
BC234R	S/B	CVD	5 3/4




BC228R




BC228R



BC224R-BC234R




BC224R




BC234R

DUROTIP TC Operating Scissors

Catalog Number	Tip	Blade	Length (in.)
BC236R	S/S	STR	5 3/4
MB171R	S/S	CVD	5 3/4



BC236R-MB171R






BC236R





MB171R

General Surgery
DUROTIP™ TC Scissors



DUROTIP TC Wire Cutting Scissors

Catalog Number	Blade	For Wire Size	Length (in.)
DP552R	1-SERR	SOFT, UP TO 1.2 MM AND HARD, UP TO 0.8 MM	4 3/4



DUROTIP TC MAYO-LEXER Dissecting Scissors, Round Blades

Catalog Number	Tip	Blade	Length (in.)
BC284R	B/B	CVD	6 1/2







BC240R-BC243R

**DUROTIP TC MAYO Dissecting Scissors,
Round Blades**


Catalog Number	Tip	Blade	Length (in.)
BC240R	B/B	STR	5 1/2
BC242R	B/B	STR	6 3/4
BC241R	B/B	CVD	5 1/2
BC243R	B/B	CVD	6 3/4




BC240R




BC241R



BC242R




BC243R




BC250R-BC253R

**DUROTIP TC MAYO Dissecting Scissors,
Beveled Blades**


Catalog Number	Tip	Blade	Length (in.)
BC250R	B/B	STR	5 1/2
BC252R	B/B	STR	6 3/4
BC251R	B/B	CVD	5 1/2
BC253R	B/B	CVD	6 3/4




BC250R



BC251R




BC252R




BC253R

General Surgery
DUROTIP™ TC Scissors




DUROTIP TC MAYO Dissecting Scissors


Catalog Number	Tip	Blade	Length (in.)
BC558R	B/B	STR	9
MG215R	B/B	STR	14
BC559R	B/B	CVD	9
MB726R	B/B	CVD	14



BC558R




BC559R



MB726R



MG215R





MG215R








DUROTIP TC METZENBAUM Dissecting Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
MB205R	S/S	STR	9 1/4
MB207R	S/S	CVD	5 3/4
MD496	S/S	CVD	7
MB206R	S/S	CVD	9 1/4




DUROTIP TC METZENBAUM Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC260R	B/B	STR	5 3/4
BC262R	B/B	STR	7
BC261R	B/B	CVD	5 3/4
BC263R	B/B	CVD	7
BC265R	B/B	CVD	7 7/8




General Surgery
DUROTIP™ TC Scissors




DUROTIP TC METZENBAUM Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC287R	B/B	CVD	9



DUROTIP TC METZENBAUM Dissecting Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC270R	B/B	STR	7
BC274R	B/B	STR	7 7/8
BC276R	B/B	STR	9
BC271R	B/B	CVD	7
BC275R	B/B	CVD	7 7/8
BC277R	B/B	CVD	9







BC258R-BC259R

**DUROTIP TC Baby METZENBAUM
Dissecting Scissors, Delicate**


Catalog Number	Tip	Blade	Length (in.)
BC258R	B/B	STR	5 3/4
BC259R	B/B	CVD	5 3/4



BC258R




BC259R




BC266R-BC285R

**DUROTIP TC NELSON-METZENBAUM
Dissecting Scissors**


Catalog Number	Tip	Blade	Length (in.)
BC266R	B/B	STR	9
BC267R	B/B	CVD	9
BC279R	B/B	CVD	10 1/4
BC281R	B/B	CVD	11
BC282R	B/B	CVD	11 1/4
BC289R	B/B	CVD	11 3/4
BC285R	B/B	CVD	13 3/4




BC266R




BC267R




BC279R




BC281R



BC282R





BC289R



BC285R



General Surgery

DUROTIP™ TC Scissors




DUROTIP TC TOENNIS-ADSON
Dissecting Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC273R	B/B	CVD	6 7/8



DUROTIP TC LEXER Dissecting Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC283R	B/B	CVD	6 1/2






Wavecut

- Wavecut scissors have one serrated cutting edge to prevent tissue from slipping. The W after the item number indicates the scissors are the wavecut style, e.g. BC211W.




DUROTIP TC MAYO Wavecut Dissecting Scissors, Beveled Blades

Catalog Number	Tip	Blade	Length (in.)
BC253W	B/B	CVD	6 3/4




DUROTIP TC Baby METZENBAUM Wavecut Dissecting Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC259W	B/B	CVD	5 3/4




General Surgery


DUROTIP™ TC Scissors



BC261W-BC265W




BC261W




BC265W

DUROTIP TC METZENBAUM Wavecut Dissecting Scissors


Catalog Number	Tip	Blade	Length (in.)
BC261W	B/B	CVD	5 3/4
BC263W	B/B	CVD	7
BC265W	B/B	CVD	7 7/8




BC263W



BC271W-BC277W




BC271W




BC277W

DUROTIP TC METZENBAUM Wavecut Dissecting Scissors, Delicate


Catalog Number	Tip	Blade	Length (in.)
BC271W	B/B	CVD	7
BC275W	B/B	CVD	7 7/8
BC277W	B/B	CVD	9



BC275W



BC287W




BC287W

DUROTIP TC METZENBAUM Wavecut Dissecting Scissors, Heavy

Catalog Number	Tip	Blade	Length (in.)
BC287W	B/B	CVD	9





DUROTIP TC TOENNIS-ADSON Wavecut Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC273W	B/B	CVD	6 7/8





DUROTIP TC NELSON-METZENBAUM Wavecut Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC267W	B/B	CVD	9
BC279W	B/B	CVD	10 1/4
BC281W	B/B	CVD	11




General Surgery
DUROTIP™ TC Scissors



DUROTIP TC MAYO-LEXER Wavecut
Scissors

Catalog Number	Tip	Blade	Length (in.)
BC284W	B/B	CVD	6 1/2




DUROTIP TC Ligature Scissors

- Carbide inserts guarantee permanent sharp cutting edge (also in the distal area of the blades).
- Wavecut prevents suture materials from slipping.
- Identified by 1 gold-plated ring.




DUROTIP TC Ligature Wavecut Scissors


Catalog Number	Tip	Blade	Length (in.)
BC294W	B/B	CVD	5 3/4
BC295W	B/B	CVD	7
BC296W	B/B	CVD	9
BC297W	B/B	CVD	10 1/4
BC298W	B/B	CVD	13 3/4




BC294W




BC295W




BC296W



BC297W



BC298W



General Surgery

SuperCut Scissors

SuperCut Scissors


Your everyday black handle scissors

- Microtome Edge
 - Razor sharp cutting edge facilitates gentle cutting with minimal force
- Wavecut Blade
 - Prevents tissue from slipping out during the cutting process


For more information about SuperCut Scissors see brochure number **DOC821**.



BC900R-BC901R



BC900R




BC901R

Iris SuperCut Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC900R	S/S	STR	4 1/4
BC901R	S/S	CVD	4 1/4



BC903R




BC903R

SuperCut Plastic Standard Scissors


Catalog Number	Tip	Blade	Length (in.)
BC903R	S/S	CVD	5



BC904R-BC905R



BC904R



BC905R

STEVENS SuperCut Scissors

Catalog Number	Tip	Tip Size (mm)	Blade	Length (in.)
BC904R	B/B	12	STR	4 1/8
BC905R	B/B	12	CVD	4 1/8




General Surgery

SuperCut Scissors

GENERAL
SURGERY


STEVENS SuperCut Scissors (continued)

Catalog Number	Tip	Tip Size (mm)	Blade	Length (in.)
BC906R	B/B	13	STR	4 1/4
BC907R	B/B	13	CVD	4 1/4



GOLDMANN-FOX SuperCut Scissors



Catalog Number	Tip	Blade	Length (in.)
BC909R	S/S	CVD	5 1/8



INCHES


General Surgery

SuperCut Scissors




JAMESON SuperCut Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC911R	B/B	CVD	5 1/8
BC915R	B/B	CVD	6
BC919R	B/B	CVD	7



JOSEPH SuperCut Scissors, Delicate


Catalog Number	Tip	Blade	Length (in.)
BC916R	S/S	STR	6
BC917R	S/S	CVD	6



General Surgery

SuperCut Scissors

GENERAL
SURGERY




GORNEY SuperCut Scissors, Delicate


Catalog Number	Tip	Blade	Length (in.)
BC920R	B/B	STR	7 3/4
BC921R	B/B	CVD	7 3/4



BC920R-BC921R




BC920R




BC921R

KAYE SuperCut Facelift Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC933R	B/B	CVD	7 3/4




BC933R




BC933R

MAYO SuperCut Scissors, Delicate


Catalog Number	Tip	Blade	Length (in.)
BC926R	B/B	STR, BEVELED	5 3/4
BC927R	B/B	CVD, BEVELED	5 3/4
BC924R	B/B	STR, ROUND	6 1/2
BC925R	B/B	CVD, ROUND	6 1/2



BC926R-BC927R




BC926R




BC927R



BC924R-BC925R



BC924R



BC925R

INCHES


General Surgery

SuperCut Scissors




MAYO SuperCut Scissors, Beveled


Catalog Number	Tip	Blade	Length (in.)
MG216R	B/B	STR	9
BC922R	B/B	STR	10
BC938R	B/B	STR	12
MG087R	B/B	CVD	9
BC929R	B/B	CVD	11




BC922R-MG216R




MG216R




BC922R



BC938R



MG087R



BC929R

General Surgery

SuperCut Scissors


GENERAL
SURGERY




BC930R-BC931R

Baby METZENBAUM SuperCut Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC930R	B/B	STR	4 1/2
BC931R	B/B	CVD	4 1/2



BC930R




BC931R




BC934R-BC937R

METZENBAUM SuperCut Scissors, Delicate


Catalog Number	Tip	Blade	Length (in.)
BC934R	B/B	STR	5 3/4
BC936R	B/B	STR	7
BC935R	B/B	CVD	5 3/4
BC937R	B/B	CVD	7




BC934R



BC935R




BC936R



BC937R

INCHES


General Surgery
SuperCut Scissors




BC945R-BC953R

METZENBAUM SuperCut Scissors,
Delicate


Catalog Number	Tip	Blade	Length (in.)
BC945R	B/B	CVD	9
BC947R	B/B	CVD	10 1/4
BC949R	B/B	CVD	11 1/4
BC951R	B/B	CVD	11 3/4
BC953R	B/B	CVD	14




BC945R




BC947R




BC949R



BC951R



BC953R



BC957R-BC959R

POTTS-SMITH SuperCut Scissors

Catalog Number	Tip	Blade Angle	Blade Size (mm)	Length (in.)
BC957R	B/B	25°	22	7 1/2
BC958R	B/B	45°	20	7 1/4
BC959R	B/B	60°	19	7



BC957R




BC958R




BC959R

General Surgery

Noir™ Supreme Cut Scissors



Microtome Edge



Wavecut

Noir Supreme Cut Scissors

■ Noir

The combination of the renowned Aesculap Tungsten Carbide technology with the Aesculap Noir coating results in an extremely high surface hardness and therefore high wear resistance in daily hospital use.

■ Microtome edge

The extremely sharp cutting edge of the scissors blade facilitates gentle cutting of tissue with minimal use of force.

■ Wavecut prevents tissue from slipping out during the cutting process.

■ The WB after the item number indicates the Noir scissors are the wavecut style, e.g. BC211WB.

For more information about Noir Supreme Cut Scissors see our brochure number DOC637.



BC252B-BC253WB



BC252B



BC253B, BC253WB

Noir MAYO Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC252B	B/B	STR	6 3/4
BC253B	B/B	CVD	6 3/4
BC253WB	B/B	CVD	6 3/4



BC284B-BC284WB



BC284B-BC284WB

Noir MAYO-LEXER Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC284B	B/B	CVD	6 1/2
BC284WB	B/B	CVD	6 1/2




General Surgery

Noir™ Supreme Cut Scissors




Noir LEXER Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC283B	B/B	DEL, CVD	6 1/2
BC283WB	B/B	DEL, CVD	6 1/2



Noir Iris Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC210B	S/S	STR	4 1/4
BC211B	S/S	CVD	4 1/4
BC211WB	S/S	CVD	4 1/4



Noir TOENNIS-ADSON Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC273B	B/B	DEL, CVD	6 7/8
BC273WB	B/B	DEL, CVD	6 7/8

General Surgery

Noir™ Supreme Cut Scissors

GENERAL
SURGERY

Noir METZENBAUM Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC261B	B/B	CVD	5 3/4
BC261WB	B/B	CVD	5 3/4
BC263B	B/B	CVD	7
BC263WB	B/B	CVD	7
BC271B	B/B	DEL, CVD	7
BC271WB	B/B	DEL, CVD	7
BC265B	B/B	CVD	7 7/8
BC265WB	B/B	CVD	7 7/8
BC275B	B/B	DEL, CVD	7 7/8
BC275WB	B/B	DEL, CVD	7 7/8
BC267B	B/B	CVD	9
BC267WB	B/B	CVD	9
BC277B	B/B	DEL, CVD	9
BC277WB	B/B	DEL, CVD	9



BC261B-BC277WB



1/1

BC261B, BC261WB



1/1

BC263B, BC263WB



1/1

BC271B, BC271WB



1/1

BC265B, BC265WB



1/1

BC275B, BC275WB



1/1

BC267B, BC267WB



1/1

BC277B, BC277WB


INCHES

General Surgery

Noir™ Supreme Cut Scissors

Noir NELSON-METZENBAUM Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC279B	B/B	DEL, CVD	10 1/4
BC279WB	B/B	DEL, CVD	10 1/4
BC281B	B/B	DEL, CVD	11
BC281WB	B/B	DEL, CVD	11



BC279B-BC281WB



BC279B, BC279WB

1/1




BC281B, BC281WB

1/1




CONVERSE Needle Holder

Catalog Number	Jaw	Length (in.)
BM200R	LONG SERR	5 1/8




HALSEY Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM201R	DEL	SMOOTH	5 1/8
BM205R	DEL	X-SERR	5 1/8




WEBSTER Needle Holder

Catalog Number	Jaw	Length (in.)
BM202R	SMOOTH	5



DERF Needle Holder


Catalog Number	Jaw	Length (in.)
BM204R	SERR	5



General Surgery

Needle Holders





BM210R


$\frac{1}{2}$



BM210R $\frac{1}{1}$

DOYEN Needle Holder

Catalog Number	Jaw	Length (in.)
BM210R	X-SERR	5



MG366R

$\frac{1}{2}$



MG366R $\frac{1}{1}$

COLLIER Needle Holder, Fenestrated

Catalog Number	Jaw	Length (in.)
MG366R	SERR	5



BM215R, MG072R

$\frac{1}{2}$



MG072R $\frac{1}{1}$




BM215R $\frac{1}{1}$

BAUMGARTNER Needle Holder


Catalog Number	Jaw	Length (in.)
MG072R	SERR, GROOVED	5 $\frac{1}{4}$
BM215R	X-SERR	5 $\frac{3}{4}$

CRILE-WOOD Needle Holder


Catalog Number	Jaw	Length (in.)
BM218R	X-SERR	6
BM220R	X-SERR	6 1/4



BM218R-BM220R




BM218R



BM220R

CRILE-MURRAY Needle Holder

Catalog Number	Jaw	Length (in.)
BM219R	SERR	6




BM219R



BM219R

ADSON Needle Holder, Fenestrated

Catalog Number	Jaw	Length (in.)
BM224R	SERR	7



BM224R




BM224R



1/1

General Surgery

Needle Holders



BM235R-BM237R



BM235R




BM236R



BM237R

MAYO-HEGAR Needle Holder

Catalog Number	Jaw	Length (in.)
BM235R	X-SERR	6
BM236R	X-SERR	7 1/4
BM237R	X-SERR	7 7/8



BM245R-BM255R



BM245R



BM247R




BM249R




BM255R

HEGAR Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM245R	HEAVY	X-SERR	6
BM247R	HEAVY	X-SERR	6 7/8
BM249R	HEAVY	X-SERR	8
BM255R	HEAVY	X-SERR	8



BM259R



BM259R

HEANEY Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM259R	CVD	SERR	7 7/8

JAMESON Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM270R	DEL	X-SERR	9

BM270R



WANGENSTEEN Needle Holder

Catalog Number	Jaw	Length (in.)
BM272R	X-SERR	10 3/4

BM272R



MASSON Needle Holder


Catalog Number	Jaw	Length (in.)
BM275R	X-SERR	10 1/2

BM275R





General Surgery


Needle Holders



BM290R-BM291R



BM290R




BM291R

HEGAR-OLSEN Needle Holder/Scissors

Catalog Number	Jaw	Length (in.)
BM290R	SERR	5 1/2
BM291R	SERR	7 1/4



BM299R



BM299R

GILLIES Needle Holder/Scissors



Catalog Number	Jaw	Length (in.)
BM299R	SERR	6 1/4

General Surgery

DUROGRIP™ TC Needle Holders

Aesculap features Needle Holders with Tungsten Carbide (TC) inserts.

A smooth jaw and three different serration patterns are available, depending on the size of the suture needle.



smooth



0.2
(16,000 teeth/in.²)



0.4
(4,000 teeth/in.²)



0.5
(2,5000 teeth/in.²)

We recommend carbide inserts as follows:

Jaw Surface

Smooth
0.2 mm pitch of serration
0.4 mm pitch of serration
0.5 mm pitch of serration

Suture Material


9/0 - 11/0
6/0 - 10/0
4/0 - 6/0
up to 3/0




BM320R

DUROGRIP TC HELLBERG-KUPKA Needle Holder W/ Ergonomic Handle

Catalog Number	Jaw	Length (in.)
BM320R	SERR	8 1/4




BM320R



BM008R

DUROGRIP TC DERF Needle Holder


Catalog Number	Jaw	Length (in.)
BM008R	SERR	4 3/4



BM008R


General Surgery

DUROGRIP™ TC Needle Holders



DUROGRIP TC MASING Needle Holder

Catalog Number	Jaw	Length (in.)
BM009R	SERR	5 1/8



DUROGRIP TC HALSEY Needle Holder


Catalog Number	Jaw	Length (in.)
BM010R	SMOOTH	4 3/4
BM012R	SERR	5 1/8

General Surgery

DUROGRIP™ TC Needle Holders



DUROGRIP TC CONVERSE Needle Holder

Catalog Number	Jaw	Length (in.)
BM011R	SERR	5 1/8





DUROGRIP TC WEBSTER Needle Holder

Catalog Number	Jaw	Length (in.)
BM014R	SMOOTH	5 1/8
BM024R	SERR	5 1/8




DUROGRIP TC BAUMGARTNER Needle Holder

Catalog Number	Jaw	Length (in.)
BM015R	SERR	5 3/4




General Surgery

DUROGRIP™ TC Needle Holders




DUROGRIP TC NEW ORLEANS Needle Holder

Catalog Number	Jaw	Length (in.)
MG522R	SERR	7



DUROGRIP TC SAROT Needle Holder

Catalog Number	Jaw	Length (in.)
BM020R	SERR	7
BM023R	SERR	10 1/4



DUROGRIP TC Baby CRILE-WOOD Needle Holder

Catalog Number	Jaw	Length (in.)
BM013R	SERR	6


General Surgery

DUROGRIP™ TC Needle Holders







GENERAL
SURGERY

DUROGRIP TC CRILE-WOOD Needle Holder



Catalog Number	Jaw	Length (in.)
BM016R	SERR	5 3/4
BM017R	SERR	7 1/4
BM019R	SERR	7 7/8
BM030R	SERR	10 5/8
BM021R	SERR	12
MB604R	SERR	14
BM029R	SMOOTH	6



BM016R, BM021R, MB604R



BM030R



BM029R



INCHES

General Surgery

DUROGRIP™ TC Needle Holders



BM022R

1/16 0.4
BM022R

1/2

**DUROGRIP TC MAYO-HEGAR-SEELEY
Needle Holder**

Catalog Number	Jaw	Length (in.)
BM022R	SERR	8

BM025R-BM027R

1/16 0.2
BM025R
1/16 0.2
BM026R
1/16 0.2
BM027R

1/2

DUROGRIP TC Fine Needle Holder











Catalog Number	Tip	Jaw	Length (in.)
BM025R	X-DEL	SERR	7
BM026R	X-DEL	SERR	7 7/8
BM027R	X-DEL	SERR	9

General Surgery

DUROGRIP™ TC Needle Holders


DUROGRIP TC DEBAKEY Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM032R	DEL	SERR	6
BM033R	DEL	SERR	6 1/2
BM035R	DEL	SERR	7
BM034R	DEL	SERR	8 1/4
BM036R	DEL	SERR	9
BM037R	DEL	SERR	10 1/4
BM038R	DEL	SERR	12
BM063R	X-DEL	SERR	6 1/2
BM064R	X-DEL	SERR	8 1/4
BM040R	X-DEL	SERR	12




General Surgery

DUROGRIP™ TC Needle Holders



BM054R-BM058R




BM054R




BM055R




BM056R




BM057R




BM058R




BM044R-BM047R




BM044R




BM045R




BM046R




BM047R




BM074R-BM076R



BM074R



BM075R








BM076R

DUROGRIP TC RYDER Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM054R	X-DEL	SERR	5 1/4
BM055R	X-DEL	SERR	6 1/8
BM056R	X-DEL	SERR	6 7/8
BM057R	X-DEL	SERR	8 1/4
BM058R	X-DEL	SERR	9
BM044R	DEL	SERR	5 1/4
BM045R	DEL	SERR	6 7/8
BM046R	DEL	SERR	7 3/4
BM047R	DEL	SERR	10 1/2
BM074R	DEL	SMOOTH	6
BM075R	DEL	SMOOTH	7
BM076R	DEL	SMOOTH	8 7/8

DUROGRIP TC FRENCH Eye Needle Holder

Catalog Number	Jaw	Length (in.)
MB606R	SERR	6 1/4
MB607R	SERR	7 1/4
MB611R	SERR	8
MB612R	SERR	9 1/4



General Surgery


DUROGRIP™ TC Needle Holders




BM065R-BM070R

DUROGRIP TC MAYO-HEGAR Needle Holder

Catalog Number	Jaw	Length (in.)
BM065R	SERR	6
BM066R	SERR	7 1/4
BM067R	SERR	8
BM068R	SERR	9 1/4
BM070R	SERR	10 1/2
BM069R	SERR	12
MB617R	SERR	14



BM065R



BM066R




BM067R




BM068R




BM070R



BM069R



BM069R






MB617R



MB617R



DUROGRIP TC HEGAR Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM077R	HEAVY	SERR	8
BM078R	HEAVY	SERR	9 1/2



DUROGRIP TC WANGENSTEEN Needle Holder

Catalog Number	Jaw	Length (in.)
BM042R	SERR	10 3/4




General Surgery

DUROGRIP™ TC Needle Holders

DUROGRIP TC MASSON Needle Holder

Catalog Number	Jaw	Length (in.)
BM090R	SERR	10 1/2




BM090R

1/2

BM090R

DUROGRIP TC MAYO-MASSON Needle Holder

Catalog Number	Jaw	Length (in.)
MB918R	SERR	16



MB918R

MB918R


General Surgery

DUROGRIP™ TC Needle Holders

GENERAL
SURGERY


DUROGRIP TC HEANEY Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM081R	CVD	SERR	7 7/8
BM084R	CVD	SERR	10 1/4



DUROGRIP TC MUELLER Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM089R	CVD	SERR	14 1/4





INCHES

General Surgery

DUROGRIP™ TC Needle Holders




DUROGRIP TC FINOCHIETTO Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM092R	ANG	SERR	10 1/2




DUROGRIP TC STRATTE Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM095R	CVD	SERR	8 1/4
BM096R	CVD	SERR	9





DUROGRIP TC Wire Twister


Catalog Number	Tip	Jaw	Length (in.)
BM110R	ROUND	SERR	6 1/8
BM111R	SQUARE	SERR	6



BM110R




BM111R




DUROGRIP TC BAUMGARTNER Wire Twister

Catalog Number	Tip	Jaw	Length (in.)
BM112R	SQUARE	SERR	7 7/8




BM112R




General Surgery

DUROGRIP™ TC Needle Holders



DUROGRIP TC BERRY Sternal Wire Twister

Catalog Number	Tip	Jaw	Length (in.)
MD456R	ROUND	SERR	7 1/4



DUROGRIP TC RUBIO Wire Twister

Catalog Number	Tip	Jaw	Length (in.)
MB550R	ROUND	SERR	5


General Surgery

DUROGRIP™ TC Needle Holders

GENERAL
SURGERY


DUROGRIP TC MATHIEU Needle Holder

Catalog Number	Jaw	Length (in.)
BM149R	SERR	5 1/2
BM151R	SERR	6 1/2
BM154R	SERR	7 7/8



DUROGRIP TC MATHIEU Needle Holder, Heavy

Catalog Number	Jaw	Length (in.)
BM156R	SERR	7 7/8
BM158R	SERR	9 5/8



INCHES

General Surgery

Needle Cases



BL907R

Needle Case, Rectangle

Catalog Number	Dimensions (mm)
BL907R	70 x 53 x 9



BL930R-BL933R



BL931R-BL933R

Needle Case, Round

Catalog Number	Item	Dimensions (mm)
BL930R	COMPLETE NEEDLE CASE W/ INTERIOR BOX	65 x 16
BL931R	INTERIOR BOX ONLY	62 x 12
BL933R	INTERIOR BOX ONLY	90 x 12




BL941R

Needle Case, With 7 Sections


Catalog Number	Dimensions (mm)
BL941R	150 x 90 x 10



BM822R-BM824R



BM822R




BM824R




BM749R, BM781R, BM807R, BM809R,
BM811R, BM813R




BM781R



BM749R




BM809R



BM807R



BM811R



BM813R

DESCHAMPS Needle



Catalog Number	Hand Orientation	Tip	Length (in.)
BM822R	RIGHT	BLUNT, 3.5 MM	9 1/2
BM824R	RIGHT	BLUNT, 5 MM	11
BM781R	LEFT	SHARP, 2 MM	8 1/4
BM749R	LEFT	SHARP, 3.5 MM	8 1/2
BM809R	LEFT	BLUNT, 2.5 MM	8 1/4
BM807R	LEFT	BLUNT, 2.5 MM	8 1/2
BM811R	LEFT	BLUNT, 3.5 MM	8 1/2
BM813R	LEFT	STEEP, BLUNT, 2 MM	8 7/8

General Surgery

Ligature Needles




BM770R



BM770R



BM771R




BM771R

HURD Needle


Catalog Number	Hand Orientation	Tip	Length (in.)
BM770R	RIGHT	SHARP	8 1/4
BM771R	LEFT	SHARP	8 1/4





BM801R



BM801R



BM802R





BM802R

KRONECKER Needle

Catalog Number	Hand Orientation	Tip	Length (in.)
BM801R	RIGHT	BLUNT	6 1/8
BM802R	LEFT	BLUNT	6 1/8



BM855R




BM855R

SYME Suture Needle


Catalog Number	Tip	Length (in.)
BM855R	BLUNT	6 3/4

Monopolar Diathermi Forceps,
Insulated


Catalog Number	Length (in.)
GK033R	13 3/4




BD043R-BD053R




BD043R




BD045R




BD047R




BD049R




BD050R



BD051R



BD052R




BD053R

Dressing Forceps, Standard Patterns

Catalog Number	Tip	Length (in.)
BD043R	SERR	4 1/2
BD045R	SERR	5 1/8
BD047R	SERR	5 3/4
BD049R	SERR	6 1/4
BD050R	SERR	7
BD051R	SERR	7 7/8
BD052R	SERR	9 3/4
BD053R	SERR	11 3/4



BD077R



BD077R

Dressing Forceps, Extra Wide


Catalog Number	Tip	Length (in.)
BD077R	SERR	5 3/4

General Surgery


Dressing and Tissue Forceps



BD023R-BD033R



BD023R




BD025R




BD027R



BD029R




BD033R



BD035R-BD037R



BD035R




BD037R

Dressing Forceps, Delicate

Catalog Number	Tip	Length (in.)
BD023R	STR, SERR	4 1/2
BD025R	STR, SERR	5 1/8
BD027R	STR, SERR	5 3/4
BD029R	STR, SERR	6 1/4
BD033R	STR, SERR	12
BD035R	CVD, SERR	5 1/8
BD037R	CVD, SERR	5 1/2



BD210R





BD210R

Dressing Forceps, Extra Delicate


Catalog Number	Tip	Length (in.)
BD210R	SERR	4 1/2




BD533R-BD543R



MB803R




BD575R-BD581R



BD587R-BD591R



BD527R




Tissue Forceps, Delicate

Catalog Number	Tip	Teeth	Length (in.)
BD533R	STR	1 x 2	4 $\frac{1}{2}$
BD535R	STR	1 x 2	5 $\frac{1}{8}$
BD537R	STR	1 x 2	5 $\frac{3}{4}$
BD539R	STR	1 x 2	6 $\frac{1}{4}$
BD540R	STR	1 x 2	7
BD541R	STR	1 x 2	7 $\frac{7}{8}$
BD542R	STR	1 x 2	9 $\frac{3}{4}$
BD543R	STR	1 x 2	11 $\frac{3}{4}$
MB803R	STR	2 x 3	4 $\frac{1}{2}$
BD575R	STR	2 x 3	5 $\frac{1}{8}$
BD577R	STR	2 x 3	5 $\frac{3}{4}$
BD579R	STR	2 x 3	6 $\frac{1}{4}$
BD580R	STR	2 x 3	7
BD581R	STR	2 x 3	7 $\frac{7}{8}$
BD587R	STR	3 x 4	5 $\frac{3}{4}$
BD589R	STR	3 x 4	6 $\frac{1}{4}$
BD591R	STR	3 x 4	7 $\frac{7}{8}$
BD527R	CVD	1 x 2	5 $\frac{3}{4}$



BD547R-BD549R



Tissue Forceps, Medium


Catalog Number	Teeth	Length (in.)
BD547R	1 x 2	5 $\frac{3}{4}$
BD549R	1 x 2	6 $\frac{1}{4}$

General Surgery


Dressing and Tissue Forceps




BD553R-BD563R



BD553R




BD555R




BD557R




BD559R




BD560R



BD561R



BD562R




BD563R

Tissue Forceps, Standard Patterns


Catalog Number	Teeth	Length (in.)
BD553R	1 x 2	4 1/2
BD555R	1 x 2	5 1/8
BD557R	1 x 2	5 3/4
BD559R	1 x 2	6 1/4
BD560R	1 x 2	7
BD561R	1 x 2	7 7/8
BD562R	1 x 2	9 3/4
BD563R	1 x 2	11 3/4



BD690R-BD692R



BD690R



BD692R

STRASSBURGER Tissue Forceps

Catalog Number	Tip	Teeth	Length (in.)
BD690R	X-SERR	1 x 2	6
BD692R	X-SERR	1 x 2	7 7/8



MG077R



MG077R

CURTIS Tissue Forceps

Catalog Number	Tip	Length (in.)
MG077R	SERR	8 3/4



BD693R



BD693R



MG357R




MG357R

BONNEY Tissue Forceps

Catalog Number	Tip	Teeth	Length (in.)
BD693R	SERR	1 x 2	6 7/8
MG357R	SERR	2 x 3	6 3/4



BD733R




BD733R

MARTIN Tissue Forceps

Catalog Number	Teeth	Length (in.)
BD733R	7 x 8	6 1/8


General Surgery

Dressing and Tissue Forceps




HEANEY Tissue Forceps

Catalog Number	Tip	Teeth	Length (in.)
MG464R	X-SERR	1 x 2	7 1/4
MB224	X-SERR	2 x 3	5 1/2
MG465R	X-SERR	3 x 4	7 1/2



STILLE-BARRAYA CV Forceps

Catalog Number	Tip	Teeth	Length (in.)
BD725R	X-SERR	3 x 4	7
BD726R	X-SERR	3 x 4	7 7/8
BD727R	X-SERR	3 x 4	9 3/4



RUSSIAN Tissue Forceps

Catalog Number	Length (in.)
BD755R	6
BD760R	7 7/8
XG006R	9 3/4
BD761R	10 1/4
MB921R	12
MB922R	14





MB208R-MB923R




MAYO-RUSSIAN Tissue Forceps

Catalog Number	Tip	Length (in.)
MB208R	FENESTRATED	9
MB923R	FENESTRATED	14



BD750R



BD750R




MB795R



MB795R



MG078R





MG078R

SELMAN Tissue Forceps

Catalog Number	Length (in.)
BD750R	6 1/8
MB795R	7 7/8
MG078R	9 3/4


General Surgery

Dressing and Tissue Forceps





LUCAE Bayonet Forceps

Catalog Number	Tip	Teeth	Length (in.)
BD878R	SERR	-	5 1/2
BD879R	SMOOTH	1 x 2	5 1/2




GRUENWALD Bayonet Forceps

Catalog Number	Tip	Length (in.)
BD883R	SERR	7 7/8



CUSHING Bayonet Forceps

Catalog Number	Tip	Teeth	Length (in.)
US147	SERR	-	7 3/4
MB231	SMOOTH	1 x 2	7 1/4



BROWN-CUSHING Bayonet Forceps

Catalog Number	Teeth	Length (in.)
MB227	7 x 7	7 1/4

General Surgery

Dressing and Tissue Forceps

BROWN-BAHNSON Bayonet Forceps

Catalog Number	Tip	Length (in.)
BD876R	SERR	6 1/8
MB804R	SERR	8

BD876R, MB804R



BD876R, MB804R

JANSEN Bayonet Forceps

Catalog Number	Tip	Length (in.)
BD880R	SERR	6 1/4
MB236	SERR	6 1/2

BD880R, MB236




BD880R

MB236

Bayonet Forceps

Catalog Number	Tip	Teeth	Length (in.)
BD885R	SMOOTH	-	7 1/4
BD886R	SMOOTH	1 x 2	7 7/8

BD885R



BD885R

BD886R



BD886R

GERALD Bayonet Forceps

Catalog Number	Tip	Teeth	Length (in.)
US105	SERR	-	7 1/2
US104	SMOOTH	1 x 2	7 1/2

US104-US105




US105

US104

General Surgery

Dressing and Tissue Forceps




Bayonet Forceps

Catalog Number	Tip	Length (in.)
BD887R		6 1/4
FD110R	FINE	7 1/4
FD111R	BROAD	7 1/4




BD900R, BD902R, BD910R-BD911R



TROELTSCH Bayonet Forceps

Catalog Number	Tip	Length (in.)
BD900R	ANG/SERR	4
BD902R	ANG/SERR	4 1/2
BD910R	ANG/SERR	5 1/2
BD911R	ANG/SERR	6 3/4




WILDE Bayonet Forceps

Catalog Number	Tip	Length (in.)
BD904R	ANG/SERR	4 1/2

General Surgery

Dressing and Tissue Forceps



BD225R, MB291R



BD225R



MB291R

HUDSON Forceps


Catalog Number	Tip	Teeth	Length (in.)
BD225R	SERR	-	5
MB291R	SMOOTH	1 x 2	5



BD610R-BD611R



BD610R




BD611R

FERRIS-SMITH Tissue Forceps


Catalog Number	Tip	Teeth	Length (in.)
BD610R	X-SERR	1 x 2	6 7/8
BD611R	X-SERR	2 x 3	6 7/8



BD215R-BD216R




BD215R




BD216R




BD665R-BD667R



BD665R



BD667R




BD666R

SEMKEN Forceps, Delicate

Catalog Number	Tip	Teeth	Length (in.)
BD215R	SERR	-	5
BD216R	SERR	-	6
BD665R	SMOOTH	1 x 2	5
BD667R	SMOOTH	1 x 2	6
BD666R	SMOOTH	2 x 3	5




BD220R



BD220R



BD223R




BD223R

Micro ADSON Forceps, Delicate

Catalog Number	Tip	Length (in.)
BD220R	SERR, W/ PLATFORM	4 3/4
BD223R	SERR	6

General Surgery

Dressing and Tissue Forceps




ADSON Forceps, Delicate

Catalog Number	Tip	Teeth	Handle	Length (in.)
BD222R	SERR	-	-	4 3/4
MA156R	SERR, X-LONG	-	-	4 3/4
BD511R	SMOOTH	1 x 2	-	4 3/4
MG275R	SERR	1 x 2	-	4 3/4
BD512R	X-SERR W/ PLATFORM	1 x 2	-	4 3/4
MA536R	X-LONG	1 x 2	-	4 3/4
MB454R	SMOOTH	2 x 3	-	4 3/4
BD506R	SMOOTH	1 x 2	FENE-STRATED	4 3/4
BD507R	SMOOTH	1 x 2	FENE-STRATED	4 3/4
BD508R	SERR	1 x 2	FENE-STRATED	4 3/4
BD513R	X-SERR W/ PLATFORM	1 x 2	FENE-STRATED	4 3/4
BD521R	SMOOTH	1 x 2	FENE-STRATED	6


General Surgery

Dressing and Tissue Forceps




MUSTARDE Suture Removal Forceps

Catalog Number	Tip	Length (in.)
MG358R	LONGITUDINAL, SERR	5




CALLISON-ADSON Tissue Forceps, W/ Platform

Catalog Number	Teeth	Length (in.)
MA537R	1 x 2	5




ADSON-BROWN Tissue Forceps

Catalog Number	Teeth	Length (in.)
BD700R	7 x 8	4 3/4



BROWN Tissue Forceps

Catalog Number	Teeth	Length (in.)
BD701R	7 x 8	6
BD702R	7 x 8	7 7/8



MICRO ADSON Tissue Forceps, Delicate


Catalog Number	Teeth	Length (in.)
BD510R	1 x 2	4 3/4
BD520R	1 x 2	6

General Surgery

Dressing and Tissue Forceps



BD662R-BD664R




GERALD Tissue Forceps, Delicate

Catalog Number	Tip	Teeth	Length (in.)
BD662R	STR	1 x 2	6 7/8
BD664R	STR	1 x 2	9
BD663R	CVD	1 x 2	6 7/8



MG422R




GERALD-DEBAKEY Forceps, Narrow

Catalog Number	Tip	Length (in.)
MG422R	SERR	9



BD522R



ADSON Tissue Forceps, Delicate

Catalog Number	Teeth	Length (in.)
BD522R	1 x 2	7



BD523R




JEFFERSON Tissue Forceps, Delicate


Catalog Number	Teeth	Length (in.)
BD523R	1 x 2	7




BD228R-BD230R



BD228R



BD230R



BD229R

GERALD Dressing Forceps, Delicate

Catalog Number	Tip	Length (in.)
BD228R	STR, SERR	6 7/8
BD230R	STR, SERR	9
BD229R	CVD, SERR	6 7/8



BD236R



BD236R

MCINDOE Thumb Forceps, Delicate

Catalog Number	Tip	Length (in.)
BD236R	SERR	6



BD661R




BD661R

MCINDOE Tissue Forceps, Delicate

Catalog Number	Tip	Teeth	Length (in.)
BD661R	SERR	1 x 2	6



BD257R



BD257R

Delicate Forceps, Angled

Catalog Number	Tip	Length (in.)
BD257R	SERR	7 7/8

General Surgery


Dressing and Tissue Forceps




BD238R-BD243R




BD238R



BD239R



BD240R




BD243R

CUSHING Forceps, Delicate

Catalog Number	Tip	Length (in.)
BD238R	STR, SERR	7
BD239R	STR, SERR	7 7/8
BD240R	STR, SERR	9 3/4
BD243R	CVD, SERR	7 7/8



MD463




MD463

CUSHING Tissue Forceps

Catalog Number	Teeth	Length (in.)
MD463	1 x 2	7



BD246R




BD246R

BROPHY Forceps, Delicate


Catalog Number	Tip	Length (in.)
BD246R	SERR	7 7/8



BD231R-BD232R



BD231R




BD232R

Forceps, Delicate

Catalog Number	Tip	Length (in.)
BD231R	SMOOTH	7
BD232R	SMOOTH	7 7/8




BD301R-BD305R




BD301R



BD302R



BD303R



BD304R



BD305R




BD311R-BD315R




BD311R




BD312R



BD313R



BD314R




BD315R

Micro Dressing Forceps

Catalog Number	Tip	Length (in.)
BD301R	STR, SERR	3 1/2
BD302R	STR, SERR	4 1/8
BD303R	STR, SERR	4 1/2
BD304R	STR, SERR	5 1/8
BD305R	STR, SERR	5 3/4
BD311R	CVD, SERR	3 1/2
BD312R	CVD, SERR	4 1/8
BD313R	CVD, SERR	4 1/2
BD314R	CVD, SERR	5 1/8
BD315R	CVD, SERR	5 3/4


General Surgery

Dressing and Tissue Forceps




Micro Tissue Forceps

Catalog Number	Teeth	Length (in.)
BD500R	1 x 2	4
BD501R	1 x 2	4 1/4
BD502R	2 x 3	4 1/4




Micro Forceps #3 Jewelers

Catalog Number	Tip Size (mm)	Length (in.)
BD335R	0.3	5 1/4




Micro Forceps #4 Jewelers

Catalog Number	Tip Size (mm)	Length (in.)
BD330R	0.2	4 1/4




Micro Forceps #5 Jewelers

Catalog Number	Tip Size (mm)	Length (in.)
BD331R	0.15	4 1/4




Micro Forceps #7 Jewelers

Catalog Number	Tip Size (mm)	Length (in.)
BD333R	0.2	4 1/2




GILLIES Tissue Forceps, Delicate

Catalog Number	Tip	Teeth	Length (in.)
BD660R	X-SERR	1 x 2	6 1/8




SEMKEN Tissue Forceps, Delicate

Catalog Number	Teeth	Length (in.)
BD669R	1 x 2	5 3/4



WAUGH Tissue Forceps, Delicate

Catalog Number	Tip	Teeth	Length (in.)
BD670R	X-SERR	1 x 2	7
BD671R	X-SERR	1 x 2	7 7/8
BD672R	X-SERR	1 x 2	9 3/4
BD679R	X-SERR	1 x 2	12




DOCKTOR Forceps

Catalog Number	Tip	Teeth	Length (in.)
MB211R	X-SERR	1 x 2	4 3/4

General Surgery

Dressing and Tissue Forceps




DEBAKEY Atraumatic Ultra Light Forceps

Catalog Number	Tip Size (mm)	Length (in.)
FC411R	1.2	6
FC412R	1.5	7 7/8
FC413R	1.5	9 1/2
FC418R	1.5	11 3/4



DEBAKEY Atraumatic Vascular Tissue Forceps

Catalog Number	Tip Size (mm)	Length (in.)
FB394R	1	6
FB395R	1	7 3/4
FB396R	1	9 1/2



DEBAKEY Atraumatic Forceps

Catalog Number	Tip Size (mm)	Length (in.)
FB411R	1.5	6
FB412R	1.5	7 7/8
FB418R	1.5	9 1/2
FB419R	1.5	11 3/4




FB400R-FB409R




FB400R



FB402R




FB404R




FB405R




FB414R



FB415R



FB416R



FB417R




FB406R



FB407R



FB408R



FB409R

DEBAKEY Atraumatic Forceps


Catalog Number	Tip Size (mm)	Length (in.)
FB400R	2	6
FB402R	2	7 7/8
FB404R	2	9 1/2
FB405R	2	11 3/4
FB414R	2.8	6
FB415R	2.8	7 7/8
FB416R	2.8	9 1/2
FB417R	2.8	11 3/4
FB406R	3.5	7 7/8
FB407R	3.5	9 1/2
FB408R	3.5	11 3/4
FB409R	3.5	13 3/4

General Surgery


Dressing and Tissue Forceps




FB401R-FB410R



FB401R



FB403R



FB410R

DEBAKEY Atraumatic Forceps, Angled

Catalog Number	Tip Size (mm)	Length (in.)
FB401R	2	6
FB403R	2	7 7/8
FB410R	2	9 1/2



FB784R



FB784R



FB785R



FB785R




FB786R



FB786R

COOLEY Atraumatic Vascular Tissue Forceps

Catalog Number	Tip Size (mm)	Length (in.)
FB784R	1.5	6
FB785R	2.0	7 3/4
FB786R	2.3	9 1/2



BD250R




BD250R

COOLEY Vascular Tissue Forceps

Catalog Number	Length (in.)
BD250R	7 7/8



BD770R




MILLS Mammmary/Coronary Tissue
Forceps

Catalog Number	Length (in.)
BD770R	7



FB130R




Dilator Forceps

Catalog Number	Tip Size (mm)	Length (in.)
FB130R	0.2	4 1/2



BD980R




HUNTER Splinter Forceps

Catalog Number	Length (in.)
BD980R	4 1/8


General Surgery

DUROGRIP™ TC Tissue Forceps



DUROGRIP Tungsten Carbide Tissue Forceps

- Cross-Serrated Platform tips for firm and non-traumatic grip
- Long working life
- Gold-plated end




BD151R $\frac{1}{2}$

BD151R $\frac{1}{1}$

US109-US110 $\frac{1}{2}$

US109 $\frac{1}{1}$

US110 $\frac{1}{1}$




MB484R-MB485R $\frac{1}{2}$

MB484R $\frac{1}{1}$

MB485R $\frac{1}{1}$

Catalog Number	Tip	Teeth	Length (in.)
BD151R	X-SERR	-	4 3/4
US109	X-SERR	1 x 2	4 3/4
US110	DEL, X-SERR	1 x 2	4 3/4
MB484R	X-SERR	1 x 2	6
MB485R	DEL, X-SERR	1 x 2	6



US112 $\frac{1}{2}$

US112 $\frac{1}{1}$



Catalog Number	Tip	Teeth	Length (in.)
US112	DEL, SMOOTH	1 x 2	4 3/4

General Surgery

DUROGRIP™ TC Tissue Forceps


DUROGRIP TC DEBAKEY Needle Forceps

Catalog Number	Tip	Tip Size (mm)	Length (in.)
MB782R	-	1.5	6
MB783R	-	1.5	7 3/4
US119	NARROW	1.5	7 3/4
MB779R	-	2	6
MB780R	-	2	8
MB781R	-	2	10 1/2
US117	-	2	6
US118	-	2	8
US120	NARROW	2	9



GERALD Needle Forceps

Catalog Number	Tip	Teeth	Length (in.)
MG356R	X-SERR	1 x 2	6 7/8




DUROGRIP TC (POTTS-SMITH) Thumb Forceps

Catalog Number	Tip	Length (in.)
BD154R	X-SERR	5 3/4
BD156R	X-SERR	7
BD157R	X-SERR	7 7/8
BD158R	X-SERR	9
BD159R	X-SERR	9 3/4




General Surgery

DUROGRIP™ TC Tissue Forceps




DUROGRIP TC CUSHING Forceps,
W/ Scraper End

Catalog Number	Tip	Length (in.)
BD160R	X-SERR	7
BD198R	X-SERR	9 3/4



DUROGRIP TC CUSHING Forceps

Catalog Number	Tip	Length (in.)
MB308R	ANG, X-SERR	6 1/4




DUROGRIP TC SEMKEN Forceps

Catalog Number	Tip	Tip Size (mm)	Teeth	Length (in.)
US111	X-SERR	1	1 x 2	6 1/4

General Surgery
DUROGRIP™ TC Tissue Forceps



BD164R-BD166R




DUROGRIP TC STANDARD Forceps

Catalog Number	Tip	Length (in.)
BD164R	X-SERR	5 3/4
BD165R	X-SERR	7 7/8
BD166R	X-SERR	9 3/4



BD168R




DUROGRIP TC GRUENWALD Forceps

Catalog Number	Tip	Length (in.)
BD168R	X-SERR	7 1/4



BD171R-BD173R




DUROGRIP TC WANGENSTEEN Forceps



Catalog Number	Tip	Length (in.)
BD171R	X-SERR	6
BD173R	X-SERR	9

General Surgery



Abdominal and Intestinal Grasping Forceps




BF309R



BF310R, MG360R



BF311R



LAHEY Forceps


Catalog Number	Teeth	Length (in.)
BF309R	2 x 2	6 1/8
BF310R	3 x 3	6 1/8
MG360R	3 x 3	8
BF311R	4 x 4	6 1/8

General Surgery

Abdominal and Intestinal Grasping Forceps

BOYS-ALLIS Forceps

Catalog Number	Teeth	Length (in.)
EA010R	5 x 6	6 1/8




EA010R

EA010R

Baby ALLIS Forceps

Catalog Number	Teeth	Length (in.)
EA012R	4 x 5	5 1/8




EA012R

EA012R


ALLIS Forceps

Catalog Number	Teeth	Length (in.)
EA015R	4 x 5	6 1/8
EA016R	5 x 6	6 1/8
EA017R	5 x 6	7 1/2
FB938R	5 x 6	8 3/4
EA018R	5 x 6	10
MB910R	5 x 6	16



EA015R

EA015R




EA016R-MB910R

EA016R-MB910R

General Surgery


Abdominal and Intestinal Grasping Forceps




MG235R



MG235R




MG236R-MG237R



MG236R-MG237R

ALLIS Forceps, Heavy Version

Catalog Number	Teeth	Length (in.)
MG235R	4 x 5	6 1/8
MG236R	5 x 6	6 1/8
MG237R	5 x 6	7 7/8



EA020R



EA020R

ALLIS-THOMS Forceps




Catalog Number	Teeth	Length (in.)
EA020R	6 x 7	7 7/8

General Surgery

Abdominal and Intestinal Grasping Forceps


JUDD-ALLIS Forceps

Catalog Number	Teeth	Length (in.)
EA025R	3 x 4	6 1/8
EA026R	3 x 4	7 7/8




ALLIS-ADAIR Clamp

Catalog Number	Teeth	Length (in.)
MD568	9 x 10	6 1/4




LOCKWOOD-ALLIS Intestinal Clamp

Catalog Number	Teeth	Length (in.)
MG455R	4 x 5	6 1/2




General Surgery


Abdominal and Intestinal Grasping Forceps



EA092R-EA099R




EA092R




EA095R




EA096R




EA097R



EA099R



EA100R-EA101R



EA100R



EA101R

ALLIS Atraumatic Forceps


Catalog Number	Length (in.)
EA092R	6 1/4
EA095R	7 7/8
EA096R	8 3/4
EA097R	10
EA100R	11 5/8
EA099R	11 3/4
EA101R	14

General Surgery

Abdominal and Intestinal Grasping Forceps


GLASSMAN-ALLIS Common Duct Forceps Non-Crushing

Catalog Number	Length (in.)
MG435R	8 1/2



BABCOCK Tissue Holding Forceps

Catalog Number	Length (in.)
MB136R	12
MB137R	14



General Surgery

Abdominal and Intestinal Grasping Forceps



EA030R-MG359R




MG359R




EA030R




EA031R




EA032R



FB934R, MB911R, MG300R



FB934R, MB911R, MG300R



MB911R

BABCOCK Forceps








Catalog Number	Tip	Length (in.)
MG359R	-	5 1/2
EA030R	-	6 1/8
EA031R	-	6 7/8
EA032R	-	8 1/2
FB934R	-	8 3/4
MB300R	-	14
MB911R	DEL, PETITE	16

General Surgery

Abdominal and Intestinal Grasping Forceps



BABCOCK Atraumatic Forceps

Catalog Number	Length (in.)
EA050R	6 1/4
EA052R	8 1/2
FB936R	8 3/4
MD470	9 1/2





General Surgery

Abdominal and Intestinal Grasping Forceps



DUROGRIP™ TC BABCOCK Forceps

Catalog Number	Length (in.)
MG364R	6 1/4
MB338R	8 1/4



ALLIS Intestinal Grasping Forceps





Catalog Number	Length (in.)
EA022R	7 1/2
EA023R	11 1/2

General Surgery

Abdominal and Intestinal Grasping Forceps




COLLIN Intestinal Forceps

Catalog Number	Length (in.)
EA039R	5 1/2
EA040R	7




COLLIN Atraumatic Forceps

Catalog Number	Length (in.)
EA061R	7 3/4
EA062R	9




General Surgery


Abdominal and Intestinal Grasping Forceps



EA041R-EA042R




EA041R



EA042R

DUVAL-COLLIN Intestinal Forceps


Catalog Number	Length (in.)
EA041R	7 3/4
EA042R	7 7/8




BH902R



BH902R



EA037R




EA037R

PENNINGTON Forceps

Catalog Number	Length (in.)
BH902R	6 1/8
EA037R	7 7/8

BAINBRIDGE Intestinal Forceps, Longitudinal Serrations

Catalog Number	Jaw	Length (in.)
EA106R	STR	6 1/8
US114	STR	7 1/4
EA107R	CVD	6 1/8
US115	CVD	7 1/4




EA106R-EA107R, US114-US115



EA106R



US114




EA107R

US115

ALLEN Intestinal Forceps, Longitudinal Serrated


Catalog Number	Teeth	Length (in.)
EA114R	1 x 2	6 1/8
MB609R	1 x 2	7 7/8



EA114R



EA114R






MB609R



MB609R

General Surgery

Abdominal and Intestinal Forceps




DOYEN Intestinal Forceps, Elastic
Blades

Catalog Number	Jaw	Length (in.)
EA120R	STR	6 1/2
EA126R	STR	9 1/4
EA121R	CVD	6 1/2
EA127R	CVD	9 1/4

**DOYEN Intestinal Forceps, Elastic
Blades, Longitudinal Serations**

Catalog Number	Jaw	Length (in.)
EA130R	STR	9 1/2
EA131R	CVD	9 1/2




EA130R-EA131R

EA130R

EA131R

**HARTMANN Intestinal Clamp, Elastic
Blades**

Catalog Number	Jaw	Length (in.)
EA124R	STR	7 7/8
EA125R	CVD	7 7/8





EA124R-EA125R

EA124R

EA125R

General Surgery

Abdominal and Intestinal Forceps



KOCHER Intestinal Forceps,
Elastic Blades

Catalog Number	Jaw	Length (in.)
EA116R	STR	6 1/8
EA140R	STR	8 1/4
EA145R	STR	9
EA117R	CVD	5 1/4
EA141R	CVD	8 1/4
EA146R	CVD	9

KOCHER Intestinal Forceps, Elastic
Blades

Catalog Number	Jaw	Length (in.)
EA150R	STR	8 3/4
EA151R	STR	9 3/4
EA152R	STR	11
EA155R	CVD	8 3/4
EA156R	CVD	9 3/4



EA150R-EA156R



EA150R



EA151R



EA152R



EA155R






EA156R

General Surgery

Abdominal and Intestinal Forceps

MAYO-ROBSON Intestinal Forceps,
Hard Elastic Blades

Catalog Number	Jaw	Length (in.)
EA170R	STR	9 3/4
EA171R	CVD	9 3/4



General Surgery


Abdominal and Intestinal Forceps




EA305R-EA450R

PAYR Intestinal Clamp


Catalog Number	Length (in.)
EA305R	7
EA306R	7 7/8
EA309R	10 3/4
EA450R	12 1/2




EA305R



EA309R



EA306R



EA450R



General Surgery

Abdominal and Intestinal Forceps



EA331R-EA332R

DEMARTEL-WOLFSON Anastomosis Clamp


Catalog Number	Description	Length (in.)
EA331R	SET OF 3	4
EA332R	SET OF 3	4 3/4



EA334R

Clamp Holder for Anastomosis Clamp


Catalog Number	Description	Length (in.)
EA334R	3" FOR USE WITH EA331R-EA332R	3 1/8



EA334R


KLEINSCHMIDT Appendectomy Forceps

Catalog Number	Jaw	Length (in.)
EA397R	CVD	5 1/8




BENSON Pylorus Spreader

Catalog Number	Length (in.)
EA466R	6 1/8



TAUFIC Cholangiography Clamp

Catalog Number	Length (in.)
MG454R	8 1/2



General Surgery

Abdominal and Intestinal Forceps



MD582

1/2

1/1

MD582

PRATT-SMITH Hemostatic "T" Forceps

Catalog Number	Length (in.)
MD582	6 1/8

EA323R

1/2

1/1


EA323R

Purse String Clamp

Catalog Number	Jaw Length (mm)	Length (in.)
EA323R	50	11


LOWER Gall Duct Forceps

Catalog Number	Jaw	Length (in.)
EB001R	HEAVY-CVD	7




LAHEY (SWEET) Gall Duct Forceps

Catalog Number	Jaw	Length (in.)
EB003R	CVD	7 7/8



DESJARDINS Gall Duct Forceps

Catalog Number	Jaw	Length (in.)
EB008R	HEAVY-CVD	8 3/4



General Surgery
Abdominal and Intestinal Forceps



EB004R

1/2

EB004R

1/1

LAHEY Gall Duct Forceps

Catalog Number	Jaw	Length (in.)
EB004R	HEAVY-CVD	7 1/2

EB011R-EB012R

1/2

EB011R

1/1

EB012R


1/1

GRAY Cystic Gall Duct Forceps

Catalog Number	Jaw	Length (in.)
EB011R	SHORT, CVD	8 3/4
EB012R	LONG, CVD	8 3/4

MIXTER Gall Duct Forceps

Catalog Number	Jaw	Length (in.)
MB120R	HEAVY-CVD	6 1/4




MB120R



MB120R

SHALLCROSS Cystic Duct Forceps


Catalog Number	Jaw	Length (in.)
MB126R	CVD	6 7/8
US152	FULL-CVD	8



MB126R, US152




MB126R



US152

BLAKE Gall Stone Forceps


Catalog Number	Jaw	Length (in.)
EB220R	STR	8
EB221R	CVD	8



EB220R-EB221R




EB220R




EB221R

General Surgery

Abdominal and Intestinal Forceps



EB231R-EB234R



EB231R



EB232R



EB233R



EB234R


DESJARDINS Gall Stone Forceps

Catalog Number	Tip	Length (in.)
EB231R	OVAL	8 7/8
EB232R	OVAL	8 7/8
EB233R	TAPERED	8 7/8
EB234R	TAPERED	8 7/8




PARKER-KERR Forceps

Catalog Number	Jaw	Length (in.)
EA365R	SLIGHT-CVD	9 3/4




HEIMING Kidney Stone Forceps

Catalog Number	Jaw	Length (in.)
EF047R	SLIGHT-CVD	6 3/4
EF048R	HEAVY-CVD	6 1/4



General Surgery


Abdominal and Intestinal Forceps




EF051R-EF054R

RANDALL Kidney Stone Forceps

Catalog Number	Jaw	Length (in.)
EF051R	SLIGHT-CVD	8 7/8
EF052R	CVD	8 7/8
EF053R	HEAVY-CVD	8 7/8
EF054R	U-CVD	8 7/8




EF051R



EF052R




EF053R




EF054R

**RANDALL Kidney Stone Forceps
(continued)**

Catalog Number	Jaw	Length (in.)
EF055R	CVD	7 1/2
EF056R	HEAVY-CVD	7 1/2
EB244R	90° ANG	9 5/8
EB246R	HEAVY-CVD, DEL	9




EF055R-EF056R



EF055R




EF056R



EB244R




EB244R



EB246R

General Surgery

Atraumatic Intestinal Forceps and Clamps



EA202R-EA207R

KOCHER Atraumatic Forceps,
Elastic Blades

Catalog Number	Jaw	Length (in.)
EA202R	STR	8 3/4
EA204R	STR	9 3/4
EA206R	STR	11
EA203R	CVD	8 3/4
EA205R	CVD	9 3/4
EA207R	CVD	11



EA202R



EA204R



EA206R



EA203R



EA205R





EA207R

General Surgery

Atraumatic Intestinal Forceps and Clamps


HARTMANN Atraumatic Forceps, Elastic Blades

Catalog Number	Jaw	Length (in.)
EA210R	STR	7 7/8
EA211R	CVD	7 7/8




MAYO-ROBSON Atraumatic Forceps

Catalog Number	Jaw	Length (in.)
EA216R	STR	9 3/4
EA217R	CVD	9 3/4



General Surgery

Atraumatic Intestinal Forceps and Clamps



EA212R-EA223R

DOYEN Atraumatic Forceps,
Elastic Blades

Catalog Number	Jaw	Length (in.)
EA222R	STR	6 3/4
EA212R	STR	9
EA223R	CVD	6 3/4
EA213R	CVD	9



EA222R



EA212R




EA223R



EA213R

SCUDDER Atraumatic Forceps

Catalog Number	Jaw	Length (in.)
EA240R	STR	13



EA240R




EA240R

General Surgery


Atraumatic Intestinal Forceps and Clamps

GLASSMAN Anterior Resection Clamp


Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
MB931R	90°, SERR	25	9 3/4
MB130R	90°, SERR	50	9 3/4
MB932R	90°, SERR	75	10




MB130R-MB932R



MB931R



MB130R






MB932R

General Surgery

Atraumatic Intestinal Forceps and Clamps

CARTER-GLASSMAN Resection Clamp

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
MG460R	STR	110	11
MG439R	45° ANG	103	10 1/4




General Surgery

Atraumatic Intestinal Forceps and Clamps

GLASSMAN Intestinal and Liver Holding Clamp, Non-Crushing

Catalog Number	Jaw	Length (in.)
MG436R	STR	9 3/4
MG437R	CVD	9 5/8
MG438R	FULL CVD	9



MG436R-MG438R



MG436R



MG437R



MG438R


General Surgery

Atraumatic Intestinal Forceps and Clamps



GLASSMAN Gastroenterostomy Clamp


Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
MG433R	STR	116	11 1/2
MG434R	ANG	115	11 1/2



MG433R-MG434R



MG433R




MG434R

General Surgery

Atraumatic Intestinal Forceps and Dressing and Sponge Forceps

GLASSMAN Gastrointestinal Clamp, Non-Crushing


Catalog Number	Jaw	Length (in.)
MB138R	STR	9 1/4
MB139R	ANG	9



MB138R-MB139R




MB138R




MB139R

LISTER Sponge Forceps


Catalog Number	Jaw	Length (in.)
BF004R	STR	5 1/4
BF006R	STR	6 1/4
BF008R	STR	7 1/4




BF004R-BF008R



BF004R




BF006R



BF008R

GROSS Forceps, W/ Ratchet


Catalog Number	Jaw	Length (in.)
BF046R	STR	7 7/8
BF047R	CVD	7 7/8



BF046R-BF047R




BF046R




BF047R

General Surgery

Dressing and Sponge Forceps






SIMS-MAIER Sponge Forceps		
Catalog Number	Jaw	Length (in.)
BF078R	STR	13 3/4



ULRICH-AESCULAP Sponge Forceps		
Catalog Number	Jaw	Length (in.)
BF062R	STR	10 1/2
BF063R	CVD	10 1/2
BF065R	CVD	13 3/4

RAMPLEY Sponge Forceps

Catalog Number	Tip	Jaw	Length (in.)
BF116R	SERR	STR	7
BF118R	SERR	STR	9 3/4




General Surgery


Dressing and Sponge Forceps

FOERSTER Sponge Forceps

Catalog Number	Tip	Jaw	Length (in.)
MB098R	SERR	STR	12
MB099R	SERR	STR	14
MB904R	SERR	STR	16



MB098R-MB904R



MB098R









MB099R



MB904R


FOERSTER-BALLENGER Sponge Forceps

Catalog Number	Tip	Jaw	Length (in.)
BF112R	SERR	STR	7
BF120R	SMOOTH	STR	9 5/8
BF122R	SERR	STR	9 5/8
BF113R	SERR	CVD	7
BF121R	SMOOTH	CVD	9 5/8
BF123R	SERR	CVD	9 5/8





General Surgery

Dressing and Sponge Forceps




BALLENGER Sponge Forceps

Catalog Number	Tip	Jaw	Length (in.)
MB232	SMOOTH	STR	7



KELLY Sponge Forceps

Catalog Number	Tip	Jaw	Handle	Length (in.)
BF135R	SERR	CVD	W/ OUT RATCHET	12 1/2
BF136R	SERR	CVD	W/ RATCHET	12 1/2




General Surgery

Tongue Forceps

MAGILL Forceps


Catalog Number	Size	Length (in.)
AN379R	INFANT	6
AN380R	PEDIATRIC	7 7/8
AN381R	ADULT	9 3/4



AN379R-AN381R



AN379R




AN380R




AN381R

COLLIN Forceps W/ Ratchet


Catalog Number	Jaw	Length (in.)
BH906R	OVAL	6 1/4
BH908R	TEARDROP	6 1/4



BH906R-BH908R




BH906R




BH908R

General Surgery


Tongue Forceps




AN710R



AN710R



AN711R




AN711R

HEYWOOD-SMITH Forceps

Catalog Number	Tip	Length (in.)
AN710R	SERR	7 7/8
AN711R	SERR	8 1/4

YEOMAN Rectal Biopsy Forceps

Catalog Number	Tip Size (mm)	Shaft Length (in.)
EA903R	4	16 1/2




General Surgery

Biopsy Forceps

AESCULAP Rectal Swab Holder


Catalog Number	Shaft Length (in.)
EA943R	15 3/4



PRICE Muscle Biopsy Clamp

Catalog Number	Jaw	Jaw Size (mm)	Length (in.)
MB200R	DOUBLE-PRONG	8	4
MG098R	DOUBLE-PRONG	15	4 3/8





BH332R-BH339R



BH332R



BH334R



BH336R



BH338R



BH333R



BH335R




BH337R



BH339R

SPENCER-WELLS Forceps

Catalog Number	Jaw	Length (in.)
BH332R	STR	5 1/8
BH334R	STR	6
BH336R	STR	6 7/8
BH338R	STR	7 7/8
BH333R	CVD	5 1/8
BH335R	CVD	6
BH337R	CVD	6 7/8
BH339R	CVD	7 7/8



BH412R-BH415R



BH412R



BH414R



BH413R





BH415R

PEAN Forceps, Delicate

Catalog Number	Jaw	Length (in.)
BH412R	STR	5 1/8
BH414R	STR	5 1/2
BH413R	CVD	5 1/8
BH415R	CVD	5 1/2



General Surgery

Hemostatic Forceps and Clamps




PEAN Forceps

Catalog Number	Jaw	Length (in.)
BH410R	STR	5
BH422R	STR	5 1/8
BH424R	STR	5 1/2
BH425R	CVD	5 1/2



PEAN Artery Forceps

Catalog Number	Jaw	Length (in.)
MB180R	CVD	12
MB181R	CVD	14



BH442R-BH451R

ROCHESTER-PEAN Forceps

Catalog Number	Jaw	Length (in.)
BH442R	STR	6 1/4
BH444R	STR	7 1/4
BH446R	STR	7 7/8
BH448R	STR	8 7/8
BH450R	STR	9 1/2
BH443R	CVD	6 1/4
BH445R	CVD	7 1/4
BH447R	CVD	7 7/8
BH449R	CVD	8 7/8
BH451R	CVD	9 1/2



BH442R



BH444R



BH446R



BH448R



BH450R



BH443R




BH445R



BH447R




BH449R



BH451R

General Surgery

Hemostatic Forceps and Clamps



BH470R-BH475R, MB128R

PEAN Forceps, Long

Catalog Number	Jaw	Length (in.)
BH470R	STR	9 1/2
BH472R	STR	10 1/4
BH474R	STR	11
MB128R	STR	12
BH471R	CVD	9 1/2
BH473R	CVD	10 1/4
BH475R	CVD	11



BH470R



BH472R



BH474R




MB128R



BH471R




BH473R



BH475R

KOCHER Forceps, Delicate

Catalog Number	Jaw	Teeth	Length (in.)
BH614R	STR	1 x 2	5 1/2
BH615R	CVD	1 x 2	5 1/2




KOCHER Forceps

Catalog Number	Jaw	Teeth	Length (in.)
BH630R	STR	1 x 2	6
BH631R	CVD	1 x 2	6

General Surgery

Hemostatic Forceps and Clamps



BH618R-BH619R




BH618R



BH619R

ROCHESTER-OCHSNER Forceps

Catalog Number	Jaw	Teeth	Length (in.)
BH618R	STR	1 x 2	5 1/2
BH619R	CVD	1 x 2	5 1/2




MB172R



MB172R

ROCHESTER-OCHSNER Artery Forceps

Catalog Number	Jaw	Teeth	Length (in.)
MB172R	STR	1 x 2	14



BH642R-BH655R

KOCHER-OCHSNER Forceps

Catalog Number	Jaw	Teeth	Length (in.)
BH642R	STR	1 x 2	6 1/4
BH644R	STR	1 x 2	7 1/4
BH646R	STR	1 x 2	7 7/8
BH648R	STR	1 x 2	8 7/8
BH650R	STR	1 x 2	9 1/2
BH652R	STR	1 x 2	10 1/4
BH654R	STR	1 x 2	11
BH643R	CVD	1 x 2	6 1/4
BH645R	CVD	1 x 2	7 1/4
BH647R	CVD	1 x 2	7 7/8
BH649R	CVD	1 x 2	8 7/8
BH651R	CVD	1 x 2	9 1/2
BH653R	CVD	1 x 2	10 1/4
BH655R	CVD	1 x 2	11



BH642R



BH644R



BH646R



BH648R



BH650R



BH652R



BH654R



BH643R



BH645R



BH647R



BH649R



BH651R



BH653R




BH655R

General Surgery

Hemostatic Forceps and Clamps

ROCHESTER-OCHSNER Forceps


Catalog Number	Jaw	Teeth	Length (in.)
MB807R	STR	1 x 2	14
MB909R	STR	1 x 2	16
MB908R	CVD	1 x 2	13 3/4



MB807R, MB909R




MB807R



MB908R



MB908R



BH800R, MG432R

ROCHESTER-CARMALT Forceps

Catalog Number	Jaw	Length (in.)
BH800R	STR	6 1/4
BH804R	STR	7 7/8
MG196R	STR	9
MG431R	STR	10 1/4
BH801R	CVD	6 1/4
BH805R	CVD	7 7/8
MG197R	CVD	9
MG432R	CVD	10 1/4



BH800R



BH804R



MG196R



MG431R



BH801R




BH805R



MG197R



MG432R



MG353R-MG354R

STOREY Forceps

Catalog Number	Tip	Length (in.)
MG353R	STRONG-CVD	8 3/4
MG354R	RIGHT-ANG	8 3/4



MG353R



MG354R



General Surgery

Hemostatic Forceps and Clamps





JONES Hemostatic Forceps

Catalog Number	Jaw	Length (in.)
MG426R	STR	5



MOYNIHAN Forceps

Catalog Number	Jaw	Length (in.)
BH831R	CVD	9
BH833R	CVD	9 1/2




NEGUS Tonsil Forceps

Catalog Number	Jaw	Length (in.)
BH951R	HEAVY-CVD	7 1/2
BH952R	U-CVD	7 1/2


CARMALT Artery Forceps, Longitudinal
Serrations

Catalog Number	Jaw	Length (in.)
MB182R	SLIGHT-CVD	12




SCHNIDT Tonsil Clamp

Catalog Number	Ring Handle	Jaw	Length (in.)
MA172R	OPEN RING	STR	7 1/2
BH961R	OPEN-RING	SLIGHT-CVD	7 1/4
BH963R	OPEN-RING	CVD	7 1/4
BH957R	CLOSED-RING	SLIGHT-CVD	7 1/4
BH959R	CLOSED-RING	CVD	7 1/4



General Surgery

Hemostatic Forceps and Clamps



MD450


1/2

1/1

MD450

OCHSNER-DIXON Artery Forceps

Catalog Number	Jaw	Teeth	Length (in.)
MD450	STR	1 x 2	6 1/4



BH104R-BH115R

1/2

1/1

BH104R

1/1

BH114R

1/1


BH105R

1/1

BH115R

Baby HARTMANN Mosquito Artery Forceps

Catalog Number	Jaw	Teeth	Length (in.)
BH104R	STR	-	4
BH114R	STR	1 x 2	4
BH105R	CVD	-	4
BH115R	CVD	1 x 2	4



BH108R-BH119R

1/2

1/1

BH108R

1/1

BH118R

1/1

BH109R

1/1







BH119R

Micro HALSTED Mosquito Forceps, Delicate

Catalog Number	Jaw	Teeth	Length (in.)
BH108R	STR	-	5
BH118R	STR	1 x 2	5
BH109R	CVD	-	5
BH119R	CVD	1 x 2	5



HALSTED Mosquito Forceps, Delicate

Catalog Number	Jaw	Teeth	Length (in.)
BH110R	STR	-	5
BH120R	STR	1 x 2	5
BH210R	STR	-	7 7/8
BH214R	STR	1 x 2	7 7/8
BH111R	CVD	-	5
BH121R	CVD	1 x 2	5
BH211R	CVD	-	7 7/8
BH215R	CVD	1 x 2	7 7/8



JACOBSON Mosquito Forceps


Catalog Number	Jaw	Length (in.)
BH106R	STR	5 1/8
MB229R	CVD	5



General Surgery

Hemostatic Forceps and Clamps





MG085R


1/2



MG085R 1/1

Mosquito Forceps

Catalog Number	Jaw	Length (in.)
MG085R	DEL, CVD.	6



BJ002R


1/2



BJ002R 1/1

Suture Forceps

Catalog Number	Tip	Jaw	Length (in.)
BJ002R	SMOOTH	STR	3 3/4




BH187R

1/2



BH187R 1/1



BH189R

1/2



BH189R 1/1

JACOBSON Hemostatic Forceps


Catalog Number	Jaw	Length (in.)
BH187R	SLIGHT-CVD	7 1/2
BH189R	HEAVY-CVD	7 1/2

General Surgery

Hemostatic Forceps and Clamps


ADSON Hemostatic Forceps

Catalog Number	Jaw	Length (in.)
BH205R	CVD	8 1/2




HALSTED Forceps Delicate

Catalog Number	Jaw	Teeth	Length (in.)
BH202R	STR	-	7 1/4
BH212R	STR	1 x 2	7 1/4
BH203R	CVD	-	7 1/4
BH213R	CVD	1 x 2	7 1/4




PROVIDENCE Hospital Forceps, Delicate

Catalog Number	Jaw	Length (in.)
BH130R	STR	5 1/2
BH131R	CVD	5 1/2




General Surgery

Hemostatic Forceps and Clamps




KELLY Forceps

Catalog Number	Jaw	Length (in.)
BH134R	STR	5 1/2
BH135R	CVD	5 1/2







Baby CRILE Forceps

Catalog Number	Jaw	Teeth	Length (in.)
BH140R	STR	-	5 1/2
BH150R	STR	1 x 2	5 1/2
BH141R	CVD	-	5 1/2
BH151R	CVD	1 x 2	5 1/2




CRILE Forceps

Catalog Number	Jaw	Teeth	Length (in.)
BH144R	STR	-	5 1/2
BH154R	STR	1 x 2	5 1/2
BH166R	STR	-	6 1/4
BH176R	STR	1 x 2	6 1/4
BH145R	CVD	-	5 1/2
BH155R	CVD	1 x 2	5 1/2
BH167R	CVD	-	6 1/4
BH177R	CVD	1 x 2	6 1/4



General Surgery

Hemostatic Forceps and Clamps



MD464-MD469



MD464



MD466R



MD468



MD465R




MD467R



MD469

COLLER Artery Forceps

Catalog Number	Jaw	Length (in.)
MD464	STR	5 3/4
MD466R	STR	6 1/2
MD468	STR	7
MD465R	CVD	5 3/4
MD467R	CVD	6 1/2
MD469	CVD	7



BH160R-BH171R



BH160R



BH170R



BH161R




BH171R

LERICHE Forceps, Delicate

Catalog Number	Jaw	Teeth	Length (in.)
BH160R	STR	-	6
BH170R	STR	1 x 2	6
BH161R	CVD	-	6
BH171R	CVD	1 x 2	6


RANKIN Forceps, Delicate

Catalog Number	Jaw	Length (in.)
BH164R	STR	6 1/4
BH165R	CVD	6 1/4




BIRKETT Hemostatic Forceps

Catalog Number	Jaw	Length (in.)
BH197R	CVD	7 1/4




NISSEN Tonsil Forceps, Delicate

Catalog Number	Jaw	Length (in.)
BH199R	CVD	7 1/4




General Surgery

Hemostatic Forceps and Clamps



BH200R-BH201R




BH200R


BH201R

ADSON Tonsil Forceps, Delicate

Catalog Number	Jaw	Length (in.)
BH200R	STR	7 1/4
BH201R	CVD	7 1/4



BH206R-BH217R



BH206R

BH216R



BH207R

BH217R



BH209R

HEISS Forceps, Delicate

Catalog Number	Jaw	Teeth	Length (in.)
BH206R	STR	-	7 7/8
BH216R	STR	1 x 2	7 7/8
BH207R	CVD	-	7 7/8
BH217R	CVD	1 x 2	7 7/8
BH209R	HVY-CVD	-	7 3/4



BENGOLEA Forceps, Delicate

Catalog Number	Jaw	Length (in.)
BH228R	STR	9 5/8
BH229R	CVD	9 5/8




BRIDGE Forceps, Delicate

Catalog Number	Jaw	Length (in.)
BH230R	STR	10 3/4
BH231R	CVD	10 3/4



General Surgery

Dissecting and Right Angle Forceps



BJ100R-BJ106R

GEMINI Right Angle Forceps

Catalog Number	Length (in.)
BJ100R	5 1/2
BJ102R	7
BJ103R	7 7/8
BJ104R	9
BJ106R	11



BJ100R



BJ104R




BJ102R



BJ106R



BJ103R



BJ037R-BJ039R

HOLLE Forceps, Delicate

Catalog Number	Length (in.)
BJ037R	5 3/4
BJ038R	7 1/4
BJ039R	9 1/4



BJ037R



BJ038R



BJ039R








General Surgery

Dissecting and Right Angle Forceps

GENERAL
SURGERY

KANTROWITZ Right Angle Forceps


Catalog Number	Length (in.)
MB058R	6 3/8
BJ057R	7 7/8
BJ058R	9 5/8
BJ059R	10 3/4
MB061R	12





General Surgery


Dissecting and Right Angle Forceps



BJ013R-BJ014R




BJ013R $\frac{1}{1}$



BJ014R $\frac{1}{1}$

Baby ADSON Forceps

Catalog Number	Length (in.)
BJ013R	5 1/2
BJ014R	7




BJ008R $\frac{1}{2}$



BJ008R $\frac{1}{1}$

Micro ADSON Suture Forceps

Catalog Number	Length (in.)
BJ008R	5 1/2



MD569-MD570R



MD569 $\frac{1}{1}$



MD570R $\frac{1}{1}$

MIXTER (MEEKER) Right Angle Forceps

Catalog Number	Length (in.)
MD569	5 1/2
MD570R	7

General Surgery

Dissecting and Right Angle Forceps




BJ020R-BJ035R

OVERHOLT-GEISSENDOERFER Forceps


Catalog Number	Jaw Curve (#)	Length (in.)
BJ020R	0	7 3/4
BJ021R	1	8
BJ022R	2	8 3/4
BJ023R	3	8 7/8
BJ024R	4	8 3/4
BJ025R	5	9
BJ026R	6	8 7/8
BJ031R	1	10 1/4
BJ032R	2	10 5/8
BJ033R	3	10 3/4
BJ034R	4	10 5/8
BJ035R	5	11



BJ020R




BJ021R,
BJ031R




BJ022R,
BJ032R



BJ023R,
BJ033R



BJ024R,
BJ034R




BJ025R,
BJ035R



BJ026R

General Surgery


Dissecting and Right Angle Forceps




BJ080R-BJ096R




BJ080R



BJ081R



BJ082R



BJ083R




BJ086R



BJ090R



BJ091R





BJ096R

OVERHOLT Forceps

Catalog Number	Jaw Curve (#)	Length (in.)
BJ080R	0	8 1/2
BJ081R	1	8 1/4
BJ082R	2	8 3/4
BJ083R	3	8 7/8
BJ086R	6	8 7/8
BJ090R	0	11 5/8
BJ091R	1	11 5/8
BJ096R	6	11 3/4







MEEKER Forceps

Catalog Number	Length (in.)
MB459R	11



General Surgery

Dissecting and Right Angle Forceps



ZENKER Ligature Forceps

Catalog Number	Jaw	Length (in.)
BJ047R	SLIGHT-CVD	11 3/4
BJ048R	HEAVY-CVD	11 1/2
BJ130R	SLIGHT-CVD	13 3/4
BJ131R	CVD	13 3/4
BJ132R	HEAVY-CVD	13 3/4

General Surgery

Dissecting and Right Angle Forceps

GENERAL
SURGERY


MIXTER-O'SHAUGNESSY Forceps

Catalog Number	Length (in.)
BJ050R	7 1/2



BJ050R

1/2




BJ050R

1/1

FINOCHIETTO Forceps

Catalog Number	Length (in.)
BJ052R	9 5/8



BJ052R

1/2




BJ052R

1/1

INCHES


General Surgery

Dissecting and Right Angle Forceps



LAHEY Forceps

Catalog Number	Length (in.)
BJ060R	8 7/8



MIXTER Thoracic Forceps


Catalog Number	Length (in.)
MB059R	14

General Surgery

Dissecting and Right Angle Forceps

MIXTER Forceps, Longitudinal Serrations

Catalog Number	Length (in.)
BJ061R	8 3/4
BJ062R	9 3/4
BJ073R	11 1/2




BJ061R-BJ062R



BJ061R



BJ062R






BJ073R



BJ073R


General Surgery

Dissecting and Right Angle Forceps



MIXTER Forceps

Catalog Number	Length (in.)
BJ055R	9
MB808R	14
MB924R	15 3/4



Baby MIXTER Forceps

Catalog Number	Length (in.)
BJ011R	5 1/2
BJ012R	7


General Surgery

Dissecting and Right Angle Forceps


GENERAL
SURGERY

HARRINGTON-MIXTER Forceps,
Longitudinal Serrations

Catalog Number	Length (in.)
BJ072R	12 1/2




BJ072R



BJ072R

WIKSTROEM Forceps

Catalog Number	Length (in.)
BJ051R	8



BJ051R



BJ051R


INCHES

General Surgery

Dissecting and Right Angle Forceps


LAWRENCE Forceps

Catalog Number	Length (in.)
BJ070R	10 5/8



REINHOF SWAN NECK CLAMP

Catalog Number	Tip	Jaw	Handle	Length (in.)
US116	SERR	CVD	CVD	10 1/4




General Surgery

Dissecting and Right Angle Forceps


KIEBACK Forceps

Catalog Number	Jaw	Length (in.)
BJ075R	CVD, LONGITUDINAL SERR	9 1/4




COLLIER-ANDERSON Forceps

Catalog Number	Jaw	Length (in.)
MB897R	CVD, LONGITUDINAL SERR	10
MB898R	CVD, LONGITUDINAL SERR	11



DEBAKEY-COLLIER Forceps

Catalog Number	Jaw	Length (in.)
MB899R	CVD, SERR	8




General Surgery

Towel and Tubing Clamps




BF407R




BF407R

JONES Towel Clamp

Catalog Number	Length (in.)
BF407R	3 1/4



BF426R



BF426R

BERNHARD Towel Clamp


Catalog Number	Length (in.)
BF426R	6 1/4




BF431R-BF433R



BF431R



BF432R




BF433R

BACKHAUS Towel Clamp


Catalog Number	Length (in.)
BF431R	3 1/2
BF432R	4 1/4
BF433R	5 1/4

General Surgery


Towel and Tubing Clamps



BF436R-BF437R



BF436R




BF437R

BACKHAUS Towel Clamp, W/ Screw

Catalog Number	Length (in.)
BF436R	3 1/8
BF437R	4 1/8



BF444R




BF444R

BACKHAUS-KOCHER Towel Clamp, W/ Screw

Catalog Number	Length (in.)
BF444R	5 1/4



BF448R



BF448R

ROEDER Towel Clamp, W/ Ball On Jaws

Catalog Number	Length (in.)
BF448R	5 1/4

General Surgery

Towel and Tubing Clamps




BF464R

1/2




BF464R

1/1



BF460R

1/2




BF460R


1/1

Towel Clamp, Paper Drape

Catalog Number	Length (in.)
BF464R	5 1/2
BF460R	6




BF458R-BF459R



BF458R

1/1



BF459R

1/1

Non-Perforating Towel Forceps



Catalog Number	Length (in.)
BF458R	4
BF459R	5 1/2

General Surgery

Towel and Tubing Clamps





LORNA Non-Perforating Towel Clamp

Catalog Number	Length (in.)
BF461R	3 3/4
BF463R	5 1/8





Tubing Clamp, W/ Guard

Catalog Number	Jaw	Length (in.)
MD652	SERR	6
MD653	SERR	7
MD654	SERR	7 7/8



General Surgery


Skin Hooks




Dimensions for blade size of retractors
1st number = depth
2nd number = width
Example: 25 x 38 mm




OL620R, MA185R, OL622R-OL624R, MD581




OL620R




MA185R



OL622R



OL623R



OL624R



MD581

JOSEPH Skin Hook

Catalog Number	Tip	Prongs	Tip Size (mm)	Length (in.)
OL620R	SHARP	1	-	6 1/2
MA185R	SHARP	2	2	6 1/4
OL622R	SHARP	2	3	6 1/2
OL623R	SHARP	2	6	6 1/2
OL624R	SHARP	2	8	6 1/2
MD581	SHARP	2	10	6 1/2



OL621R




OL621R

FREER Double Hook


Catalog Number	Tip	Tip Size (mm)	Length (in.)
OL621R	SHARP	2.0 x 3.2	6



OA332R-OA333R



OA332R



OA333R

GRAEFE Hook

Catalog Number	Tip	Length (in.)
OA332R	SHARP	5
OA333R	BLUNT	5



OA334R-OA335R



OA334R




OA335R

GUTHRIE Double Hook


Catalog Number	Tip	Length (in.)
OA334R	SHARP	5
OA335R	SHARP	5 1/8




OA336R-OA338R



OA336R



OA337R




OA338R

ROLLET Hook Retractor

Catalog Number	Tip	Blade End	Tip Size (mm)	Length (in.)
OA336R	SHARP	2 PRONGS	2.2 x 4.5	5 1/8
OA337R	SHARP	4 PRONGS	2.5 x 7.0	5 1/4
OA338R	BLUNT	4 PRONGS	2.2 x 7.0	5 1/4



OL615R




OL615R

MCINDOE Hook

Catalog Number	Tip	Blade End	Tip Size (mm)	Length (in.)
OL615R	SHARP	2 PRONGS	9	6



OA340R



OA340R

PAUL Tendon Hook

Catalog Number	Tip	Blade End	Tip Size (mm)	Length (in.)
OA340R	BLUNT	3 PRONGS	11 x 8	4 3/4



OA351R



OA351R

JAEGER Hook Retractor, Fine

Catalog Number	Tip	Length (in.)
OA351R	SHARP	5 1/4



OL610R



OL610R

KILNER Skin Hook

Catalog Number	Tip	Length (in.)
OL610R	SHARP	6 1/4



OL611R-OL612R



OL611R




OL612R

GILLIES Skin Hook

Catalog Number	Tip	Hook End	Length (in.)
OL611R	SHARP	SMALL	7
OL612R	SHARP	LARGE	7




General Surgery

Skin Hooks



COTTLE Hook

Catalog Number	Tip	Hook End	Length (in.)
OL617R	SHARP	ANG	6 1/4
OL619R	SHARP	HVY-CVD	6 1/4



KLEINERT-KUTZ Skin Hook

Catalog Number	Size	Tip	Tip Size (mm)	Length (in.)
OL601R	1	SHARP	2.5	6 1/4
OL602R	2	SHARP	4.2	6 1/4
OL603R	3	SHARP	4.8	6 1/4



Prong Micro Hooks

Catalog Number	Tip	Tip Size (mm)	Length (in.)
MB885R	SHARP	1	5 5/8
MB886R	SHARP	2	5 5/8
MB887R	BLUNT	2	5 5/8
MB888R	SHARP	2	5 5/8
MB889R	SHARP	5	5 5/8
MB890R	SHARP	7	5 5/8
MB891R	SHARP	10	5 5/8



BT075R



BT075R

BROM Vein Hook

Catalog Number	Length (in.)
BT075R	6 3/4



BT100R



BT100R

BERGMANN Hook

Catalog Number	Tip	Tip Size (mm)	Length (in.)
BT100R	SHARP	7	5 1/2



BT161R



BT161R

EMMET Fistula Hook

Catalog Number	Tip	Length (in.)
BT161R	SHARP	8 3/4




BT166R



BT166R



BT167R



BT167R

SIMON Fistula Hook


Catalog Number	Tip	Blade End	Length (in.)
BT166R	SHARP	1 PRONG	8 3/4
BT167R	SHARP	2 PRONGS	8 3/4

General Surgery

Hand Held Retractors



BT006R-BT007R



BT006R


BT007R

**Baby SENN-MILLER Retractor,
Flat Handle**

Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
BT006R	SHARP	3 PRONGS	8 x 7 / 22 x 7	6 1/8
BT007R	BLUNT	3 PRONGS	8 x 7 / 22 x 7	6 1/8




BT008R




BT008R

SENN-MILLER Retractor


Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
BT008R	SHARP	3 PRONGS	8 x 7 / 18 x 5.5	6 1/2



MG276R-MG277R



MG276R



MG277R

SENN Retractor, Round Tapered Handle

Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
MG276R	SHARP	3 PRONGS	9 x 10 / 22 x 7	6 1/4
MG277R	BLUNT	3 PRONGS	9 x 10 / 22 x 7	6 1/4

General Surgery

Hand Held Retractors


GENERAL
SURGERY



BT004, MG274R



BT004R




MG274R

RAGNELL Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT004R	3 x 5 / 6.3 x 14	6
MG274R	3 x 11 / 5.4 x 17	6



MG475R



MG475R

DAVIS Double Ended Retractor

Catalog Number	Blade Size (mm)	Length (in.)
MG475R	19 x 9 / 23 x 6.4	6



BT187R

Double Ended Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT187R	20 x 15 / 20 x 19	8



MD497

GOELET Retractor

Catalog Number	Blade Size (mm)	Length (in.)
MD497	30 x 32 / 35 x 40	7 1/2

INCHES

General Surgery

Hand Held Retractors



BT010R




BT010R

LUER Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT010R	10 x 8 / 12 x 11	4 1/4




BT020R



BT020R



BT021R



BT021R

FARABEUF Retractor, 2/Set

Catalog Number	Blade Size (mm)	Length (in.)
BT020R	25 x 10 / 32 x 12 22 x 10 / 29 x 12	5
BT021R	23 x 16 / 28 x 16 19 x 16 / 24 x 16	6




MB984R

MAYO-COLLINS Retractor, 2/Set

Catalog Number	Blade Size (mm)	Length (in.)
MB984R	31 x 20 / 30 x 18 28 x 19 / 26 x 17	5 7/8

General Surgery

Hand Held Retractors




COLLIN Retractor, 2/Set

Catalog Number	Blade Size (mm)	Length (in.)
BT022R	21 x 14 / 32 x 16 17 x 14 / 28 x 16	6



CRILE Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT012R	11 x 4 / 19 x 6	4 1/4



MB727R-MB733R



CRILE Retractor


Catalog Number	Size	Blade Size (mm)	Length (in.)
MB727R	1	3 x 11 / 6 x 21	4 1/4
MB728R	2	6 x 13 / 10 x 21	6 1/8
MB729R	3	16 x 25 / 17 x 35	6 1/4
MB730R	4	22 x 29 / 29 x 44	7
MB731R	5	29 x 41 / 35 x 51	7
MB732R	6	32 x 64 / 35 x 70	7 7/8
MB733R	7	41 x 76 / 51 x 79	8 7/8

General Surgery

Hand Held Retractors



BT011R



BT011R

ROSE Tracheal Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT011R	10 x 12 / 11 x 15	5 1/8



MG474R




MG474R

Children's Hospital Retractor

Catalog Number	Blade Size (mm)	Length (in.)
MG474R	11 x 22 / 14 x 35	4 1/2



BT025R



BT025R




BT026R

PARKER Retractor, 2/Set


Catalog Number	Blade Size (mm)	Length (in.)
BT025R	22 x 15 / 22 x 15 18 x 15 / 18 x 15	5 1/4
BT026R	31 x 25 / 27 x 25 31 x 25 / 27 x 25	7




BT034R-BT036R



BT034R



BT035R



BT036R

ROUX Retractor


Catalog Number	Blade Size (mm)	Length (in.)
BT034R	20 x 22 / 25 x 29	6
BT035R	25 x 26 / 30 x 38	6 1/4
BT036R	28 x 29 / 33 x 44	6 3/4

General Surgery


Hand Held Retractors




BT030R-BT033R



BT031R



BT032R




BT033R

ROUX Retractor


Catalog Number	Blade Size (mm)	Length (in.)
BT030R	COMPLETE SET BT031R-BT033R	-
BT031R	21 x 22 / 27 x 28	5 3/4
BT032R	26 x 26 / 32 x 26	6 1/8
BT033R	30 x 30 / 36 x 44	6 1/2



MB584R-MB585R



MB584R




MB585R

ROUX Retractor

Catalog Number	Blade Size (mm)	Length (in.)
MB584R	25 x 26 / 30 x 37	6 1/2
MB585R	28 x 28 / 29 x 40	6 1/2



BT029R



BT029R

ROUX Retractor, Long


Catalog Number	Blade Size (mm)	Length (in.)
BT029R	28 x 28 / 46 x 28	7 1/2

MATHIEU Retractor, 2/Set

Catalog Number	Blade Size (mm)	Length (in.)
BT040R	42 x 13 / 41 x 13	7 7/8
	45 x 25 / 40 x 25	



BT040R



BT040R

US Army Retractor, 2/Set

Catalog Number	Blade Size (mm)	Length (in.)
BT041R	26 x 15 / 22 x 15	8 3/4
	43 x 15 / 39 x 15	




BT041R




BT041R

General Surgery
Hand Held Retractors



Hook, Small Curve

Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
BT111R	SHARP	1 PRONG	-	6 1/2
BT112R	SHARP	2 PRONGS	5.0 x 4.2	6 1/2
BT113R	SHARP	3 PRONGS	5.0 x 7.0	6 1/2
BT114R	SHARP	4 PRONGS	5.0 x 9.5	6 1/2
BT116R	BLUNT	1 PRONG	-	6 1/2
BT117R	BLUNT	2 PRONGS	4.5 x 4.2	6 1/2
BT118R	BLUNT	3 PRONGS	4.5 x 7.0	6 1/2
BT119R	BLUNT	4 PRONGS	4.5 x 9.0	6 1/2



Hook, Large Curve

Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
BT121R	SHARP	1 PRONG	-	6 1/2
BT122R	SHARP	2 PRONGS	7.5 x 7.5	6 1/2
BT123R	SHARP	3 PRONGS	7.5 x 10.0	6 1/2
BT124R	SHARP	4 PRONGS	7.5 x 14.0	6 1/2
BT126R	BLUNT	1 PRONG	-	6 1/2
BT127R	BLUNT	2 PRONGS	6.8 x 7.5	6 1/2
BT128R	BLUNT	3 PRONGS	6.8 x 10.0	6 1/2
BT129R	BLUNT	4 PRONGS	6.8 x 14.0	6 1/2

General Surgery

Hand Held Retractors


CRILE Hook and Dissector

Catalog Number	Length (in.)
MB060R	11 3/4



KLAPP Tendon Hook

Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
BT172R	BLUNT	4 PRONGS	15 x 11	6 7/8




Vein Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT179R	10 x 13	6 1/4



KOENIG Vein Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT182R	13 x 13	7
BT183R	12 x 16	7




General Surgery


Hand Held Retractors




BT193R-BT198R



BT193R



BT194R



BT195R



BT196R



BT197R



BT198R

MEYERDING Thumb Retractor

Catalog Number	Blade Size	Size (#)	Blade Size (mm)	Length (in.)
BT193R	A	1	18 x 7	7
BT194R	B	2	16 x 5	7
BT195R	C	3	7 x 4	7
BT196R	D	4	5 x 10	7
BT197R	E	5	11 x 4	7
BT198R	F	6	6 x 4	7



BT190R



BT190R

CUSHING-KOCHER Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT190R	10 x 10	6 7/8



BT184R-BT185R



BT184R



BT185R

CUSHING Vein Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT184R	10 x 13	8
BT185R	14 x 18	8

CUSHING Vein Retractor
(continued)

Catalog Number	Blade Size (mm)	Length (in.)
BT200R	14 x 11	9 3/4
BT202R	22 x 15	9 3/4



BT200R-BT202R



CUSHING Vein Retractor

Catalog Number	Blade Size (mm)	Length (in.)
MB912R	14	12 1/2
MB913R	14	14
MB914R	14	16



MB912R-MB914R



GIL VERNET Retractor, Malleable

Catalog Number	Blade Size (mm)	Length (in.)
MB766R	9	11 1/4
MB767R	11	11 1/4
MB768R	13	11 1/4
MB769R	18	11 1/4




MB766R-MB769R



General Surgery
Hand Held Retractors




BT221R




BT221R



BT222R-BT230R



BT222R



BT227R



BT223R



BT228R



BT224R



BT229R



BT225R



BT230R

VOLKMANN Retractor


Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
BT221R	SHARP	1 PRONG	8.5 x 8	8 3/4
BT222R	SHARP	2 PRONGS	9 x 8	8 7/8
BT223R	SHARP	3 PRONGS	9 x 13	8 7/8
BT224R	SHARP	4 PRONGS	9 x 19	8 7/8
BT225R	SHARP	6 PRONGS	9 x 29	8 7/8
BT227R	BLUNT	2 PRONGS	8.5 x 12	8 7/8
BT228R	BLUNT	3 PRONGS	8.5 x 13	8 7/8
BT229R	BLUNT	4 PRONGS	8.5 x 19	8 7/8
BT230R	BLUNT	6 PRONGS	8.5 x 29	8 7/8

General Surgery


Hand Held Retractors



BT240R-BT241R




BT240R




BT241R




BT242R-BT265R



BT242R



BT257R



BT252R



BT243R



BT259R



BT253R



BT244R



BT260R



BT254R



BT246R



BT261R



BT256R



BT248R



BT265R








BT258R

VOLKMANN Retractor, Small

Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
BT240R	SHARP	1 PRONG	-	8 3/4
BT241R	BLUNT	1 PRONG	-	8 3/4
BT242R	SHARP	2 PRONGS	9 x 8	8 3/4
BT243R	SHARP	3 PRONGS	9 x 13	8 3/4
BT244R	SHARP	4 PRONGS	9 x 19	8 3/4
BT246R	SHARP	6 PRONGS	9 x 29	8 3/4
BT248R	SHARP	8 PRONGS	9 x 40	8 3/4
BT257R	SEMI-S	2 PRONGS	8.5 x 8	8 3/4
BT259R	SEMI-S	3 PRONGS	8.5 x 13	8 3/4
BT260R	SEMI-S	4 PRONGS	8.5 x 19	8 3/4
BT261R	SEMI-S	6 PRONGS	8.5 x 29	8 3/4
BT265R	SEMI-S	8 PRONGS	8.5 x 40	8 3/4
BT252R	BLUNT	2 PRONGS	8 x 8	8 3/4
BT253R	BLUNT	3 PRONGS	8 x 13	8 3/4
BT254R	BLUNT	4 PRONGS	8 x 19	8 3/4
BT256R	BLUNT	6 PRONGS	8 x 29	8 3/4
BT258R	BLUNT	8 PRONGS	8 x 40	8 3/4

General Surgery

Hand Held Retractors



VOLKMANN Retractor, US Version


Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
MG260R	SHARP	2 PRONGS	20 x 11	8 3/4
MG261R	SHARP	3 PRONGS	20 x 18	8 3/4
MG262R	SHARP	4 PRONGS	20 x 26	8 3/4
MG264R	SHARP	6 PRONGS	20 x 40	8 3/4
MG265R	BLUNT	2 PRONGS	20 x 11	8 3/4
MG266R	BLUNT	3 PRONGS	20 x 18	8 3/4
MG267R	BLUNT	4 PRONGS	20 x 26	8 3/4
MG269R	BLUNT	6 PRONGS	20 x 40	8 3/4

General Surgery

Hand Held Retractors



BT262R-BT276R



BT262R




BT263R



BT264R



BT266R



BT272R



BT273R



BT274R




BT276R

KOCHER Retractor, Large

Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
BT262R	SEMI-S	2 PRONGS	16 x 9	8 3/4
BT263R	SEMI-S	3 PRONGS	16 x 14	8 3/4
BT264R	SEMI-S	4 PRONGS	16 x 20	8 3/4
BT266R	SEMI-S	6 PRONGS	16 x 30	8 3/4
BT272R	BLUNT	2 PRONGS	15.5 x 9	8 3/4
BT273R	BLUNT	3 PRONGS	15.5 x 14	8 3/4
BT274R	BLUNT	4 PRONGS	15.5 x 20	8 3/4
BT276R	BLUNT	6 PRONGS	15.5 x 30	8 3/4



BT278R-BT279R



BT279R



BT278R

KOERTE Retractor

Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
BT279R	SEMI-S	8 PRONGS	26 x 40	9 1/2
BT278R	BLUNT	8 PRONGS	26 x 40	9 1/2





General Surgery
Hand Held Retractors

ISRAEL Retractor

Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
BT284R	BLUNT	4 PRONGS	39 x 40	10
BT285R	BLUNT	5 PRONGS	47 x 51	10
BT286R	BLUNT	6 PRONGS	53 x 62	10
BT287R	BLUNT	6 PRONGS	70 x 70	10



BT284R-BT287R




General Surgery

Hand Held Retractors

ISRAEL Retractor, T-Handle

Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
MG284R	BLUNT	4 PRONGS	46 x 40	7 1/2



MG284R



MG284R

MURPHY Retractor


Catalog Number	Tip	Blade End	Length (in.)
BT293R	SHARP	3 PRONGS	7 3/4
BT294R	SHARP	4 PRONGS	7 3/4
BT296R	BLUNT	4 PRONGS	7 3/4



BT293R-BT296R



BT293R




BT296R



BT294R

KOERTE-WAGNER Retractor

Catalog Number	Tip	Blade End	Blade Size (mm)	Length (in.)
BT290R	SEMI-S	8 PRONGS	28 x 40	11



BT290R



BT290R

General Surgery



Hand Held Retractors



LANGENBECK Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT311R	30 x 10	9
BT315R	40 x 10	9
BT316R	43 x 13	9


BT311R-BT316R



SAUERBRUCH Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT351R	46 x 16	8 7/8
BT352R	55 x 22	8 7/8
BT353R	72 x 19	8 7/8

BT351R-BT353R



General Surgery


Hand Held Retractors



BT354R-BT365R

KOCHER-LANGENBECK Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT354R	25 x 6	8 1/2
BT356R	35 x 8	8 1/2
BT357R	35 x 11	8 1/2
BT358R	35 x 15	8 1/2
BT359R	41 x 11	8 1/2
BT361R	55 x 11	8 1/2
BT363R	70 x 14	8 1/2
BT364R	80 x 12	8 1/2
BT365R	80 x 16	8 1/2



BT354R



BT356R



BT357R




BT358R



BT359R



BT361R



BT363R



BT364R



BT365R



MB785R

LAHEY Goiter Retractor

Catalog Number	Blade Size (mm)	Length (in.)
MB785R	24 x 8	7 7/8



MB785R

General Surgery

Hand Held Retractors



MB748R

MCBURNEY Thyroid Retractor

Catalog Number	Blade Size (mm)	Length (in.)
MB748R	20 x 10	7 1/2



MB748R




BT404R-BT406R

MIDDELDORPF Retractor


Catalog Number	Blade Size (mm)	Length (in.)
BT404R	15 x 15	8 3/4
BT405R	20 x 22	8 1/2
BT406R	28 x 28	9 1/4



BT404R



BT406R




BT405R



BT408R

GREENE Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT408R	28 x 22	8 7/8



BT408R

DOCKHORN Cheek Retractor

Catalog Number	Blade Size (mm)	Length (in.)
D0438R	40 x 12	8 1/4
D0439R	55 x 14	8 1/4



D0438R-D0439R



D0438R




D0439R

General Surgery

Abdominal Retractors

HOESEL Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT442R	100 x 30	9 3/4
BT446R	220 x 40	9 3/4




KOCHER Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT461R	80 x 42	9 3/4
BT462R	80 x 50	9 3/4
BT463R	80 x 60	9 3/4



BT461R-BT463R




General Surgery

Abdominal Retractors



BT464R



BT464R

KRAEMER Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT464R	66 x 10	8 1/4



BT465R



BT465R

CASPAR Retractor


Catalog Number	Blade Size (mm)	Length (in.)
BT465R	67 x 17	8 1/4



BT450R-BT469R

MIASPAS® Mini ALIF LANGENBECK Retractor


Catalog Number	Blade Size (mm)	Length (in.)
BT450R	125 x 20	8 1/4
BT451R	125 x 35	8 1/4
BT466R	80 x 20	8 1/4
BT467R	110 x 20	8 1/4
BT468R	80 x 35	8 1/4
BT469R	110 x 35	8 1/4



BT470R-BT472R

RICHARDSON-EASTMAN Retractor

Catalog Number	Description	Tip	Length (in.)
BT470R	2/SET (BT471R, BT472R)	-	-
BT471R	-	SM	9 3/4
BT472R	-	LRG	10 5/8



General Surgery


Abdominal Retractors




BT475R-BT477R

RICHARDSON Retractor


Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT475R	23 x 20	1 x 3/4	9 1/2
BT476R	30 x 29	1 1/4 x 1 1/8	9 1/2
BT477R	37 x 37	1 1/2 x 1 1/2	9 1/2
BT478R	47 x 21	2 x 3/4	9 1/2



BT475R




BT476R



BT477R




BT478R




BT478R

RICHARDSON Retractor

Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT479R	65 x 41	2 1/2 x 1 5/8	9 3/4



BT479R




BT479R

General Surgery

Abdominal Retractors

KELLY Retractor

Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT481R	50 x 40	2 x 1 1/2	10 1/4
BT482R	64 x 54	2 1/2 x 2	10 1/4
BT483R	70 x 65	2 3/4 x 2 1/2	10 1/4
BT484R	90 x 76	3 1/2 x 3	10 1/4




BT481R-BT484R




BT481R



BT482R



BT483R




BT484R

General Surgery


Abdominal Retractors

SAUERBRUCH Retractor, Open

Catalog Number	Blade Size (mm)	Length (in.)
BT514R	130 x 38	10 1/4





BT514R



BT514R



General Surgery
Abdominal Retractors



HARRINGTON Lung Spatula

Catalog Number	Blade Size (mm)	Blade End	Length (in.)
BT551R	123 x 45	FLEX	11 5/8
BT552R	123 x 64	FLEX	11 5/8
MB930R	190 x 115	FLEX	12 1/2

BT551R, BT552R, MB930R

1/1
BT551R

1/1
BT552R


1/1
MB930R

HARRINGTON Retractor, Large Blade

Catalog Number	Blade Size (mm)	Blade End	Length (in.)
MG472R	195 x 115	FLEX	14 3/8



MG472R




MG472R


General Surgery
Abdominal Retractors

HARRINGTON-PEMBERTON Lung
Spatula

Catalog Number	Blade Size (mm)	Blade End	Length (in.)
BT553R	175 x 64	FLEX	11 5/8



BT553R







ZENKER Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT594R	36 x 8	9 1/2
BT595R	60 x 12	10
BT596R	75 x 12	10 5/8
BT597R	120 x 12	11 1/2





BT594R-BT597R




General Surgery

Abdominal Retractors



DEEVER Retractor

Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT600R	19	3/4	7



DEEVER Retractor, Double Ended

Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT606R	25	1	9 1/4

General Surgery


Abdominal Retractors

DEAVER Retractor

Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT601R	22	7/8	8 1/2
BT602R	25	1	11 3/4
BT603R	25	1	12 3/8



BT601R




BT601R



BT602R-BT603R



BT602R




BT603R

General Surgery

Abdominal Retractors


DEAVER Retractor
(continued)

Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT604R	38	1 1/2	11 3/4



DEAVER Retractor
(continued)


Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT605R	50	2	12 1/4
BT607R	25	1	14 3/8



BT605R-BT607R



BT605R




BT607R

General Surgery


Abdominal Retractors

DEAVER Retractor
(continued)

Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT608R	75	3	12 1/2



BT608R




BT608R

DEAVER Retractor, W/ Grip


Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT611R	25	1	12 1/4
BT613R	38	1 1/2	12 1/4
BT614R	50	2	12
BT615R	75	3	12 1/4




BT611R-BT615R




BT611R



BT614R




BT613R



BT615R

General Surgery


Abdominal Retractors



MD071R-MD076

HASSON "S" Retractor, 2/Set

Catalog Number	Blade	Blade Size (mm)	Length (in.)
MD076	MALLEABLE	10	6 1/4
MD077	MALLEABLE	13	7
MD071R	RIGID	10	7 1/8
MD072R	RIGID	13	7 1/8



MD071R-MD076



BT750R

TUFFIER Spatula, Malleable

Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT750R	17 / 25	3/4 / 1	7 7/8



BT751R-BT763R, MD498, MG480R-MG484R

Ribbon Spatula, Malleable

Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
MG480R	5	1/8	8
BT751R	12	1/2	7 7/8
BT752R	17	3/4	7 7/8
MG481R	25	1	8
BT755R	25	1	9 3/4
MD498	75	3	13
MG482R	20	3/4	13
BT756R	25	1	13
BT761R	30	1 1/4	13
BT762R	40	1 1/2	13
BT763R	50	2	13
MG483R	65	2 1/2	13
MG484R	100	4	13



BT758R

HABERER Spatula, Malleable

Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT758R	40 / 50	1 1/2 / 2	12



MG374R

Spatula and Packer, Double Ended

Catalog Number	Length (in.)
MG374R	5 3/4



MG374R

General Surgery
Abdominal Retractors

DOYEN Abdominal Retractor


Catalog Number	Blade Size (mm)	Length (in.)
BT667R	53 x 80	9 3/4
BT669R	72 x 82	10 1/4




BT667R



BT669R



BT667R



BT669R


WEINBERG Vagotomy Retractor, Joe's
Hoe

Catalog Number	Blade Size (mm)	Length (in.)
MG491R	47 x 57 / 93 x 157	9



MG491R


1/2




General Surgery
Abdominal Retractors

ST. MARKS Pelvis Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT670R	124 x 60	11 1/2
BT672R	174 x 60	13




BT670R



BT672R

ST. MARKS Pelvis Retractor
(continued)

Catalog Number	Blade Size (mm)	Length (in.)
BT674R	178 x 60	13




BT674R

General Surgery

Abdominal Retractors

MUELLER Retractor


Catalog Number	Blade Size (mm)	Length (in.)
BT626R	250 x 30	10 5/8



BT626R


MAYO-COLLINS Appendectomy Retractor

Catalog Number	Length (in.)
MB129R	6 3/4



MB129R

1/2




BT708R-BT709R

**ROCHARD Retractor Mounting Frame,
F/ BT711R-BT715R**

Catalog Number	Mounting Arms
BT708R	DOUBLE
BT709R	SINGLE


Mounting frame for fastening to the lateral bars of the operating tables, can be dismantled for decontamination and sterilization.



BT710R

**ROCHARD Fixation Device,
F/ BT711R-BT715R**

Catalog Number	Length (mm)	Length (in.)
BT710R	190	7 1/2



BT711R-BT715R

ROCHARD Blades

Catalog Number	Blade Size (mm)	Length (in.)
BT711R	48 x 90	12 1/4
BT712R	48 x 105	12 1/4
BT713R	48 x 120	12 1/4
BT714R	48 x 135	12 1/4
BT715R	48 x 155	12 1/4

General Surgery


Abdominal Retractors

Abdominal Retractor, Malleable,
W/ Handle

Catalog Number	Blade Size (mm)	Blade Size (in.)	Length (in.)
BT773R	30	1 1/4	16 1/8
BT774R	40	1 1/2	16 1/8
BT775R	50	2	16 1/8
BT776R	60	2 1/2	16 1/8
BT777R	70	2 3/4	16 1/8



BT773R-BT777R




BT773R-BT777R



MB226, MG582R-MG586R

Laparotomy Ring

Catalog Number	Diameter (in.)
MG582R	1 1/2
MG585R	1 3/4
MB226	2
MG586R	3




BV512R

HERTZLER Abdominal Retractor

Catalog Number	Blade Size (mm)
BV512R	20 x 25

General Surgery


Abdominal Retractors



BV520R-BV530R



BV551R



BV551R

COLLIN Abdominal Retractor W/ Ball Snap Closure


Catalog Number	Description
BV530R	SET COMPLETE
BV520R	SPREADER ONLY
BV551R	LATERAL BLADES, Pair, 37 x 62 MM



BV399R

Blade Ejector


Catalog Number	Description	Length (in.)
BV399R	FOR CHANGING BLADES WITH BALL SNAP CLOSURE	4 1/2



BV500R

GOSSET Abdominal Retractor, Small

Catalog Number	Blade Depth (mm)
BV500R	38




BV502R

GOSSET Abdominal Retractor, Medium

Catalog Number	Blade Depth (mm)
BV502R	55

General Surgery

Abdominal Retractors




BV505R-BV506R

GOSSET Abdominal Retractor, Large

Catalog Number	Description
BV505R	SET COMPLETE
BV506R	CENTER BLADE ONLY 59 x 57 MM

BV505R Includes:
1 Center Blade BV506R




BV601R-BV602R

BALFOUR Baby Abdominal Retractor

Catalog Number	Description
BV601R	SET COMPLETE
BV602R	CENTER BLADE ONLY, 21 x 24 MM

BV601R Includes:
1 Center Blade BV602R




BALFOUR Abdominal Retractor

Catalog Number	Description	Lateral Blades Depth (mm)
BV605R	BALFOUR RETRACTOR, 200 x 155 MM	60
BV608R	BALFOUR RETRACTOR, 200 x 235 MM	60
BV610R	BALFOUR RETRACTOR, 200 x 265 MM	90
BV606R	CENTER BLADE, 47 x 80 MM	-
BV609R	CENTER BLADE, 62 x 76 MM	-
BV611R	CENTER BLADE, 80 x 80 MM	-

BV605R includes BV606R


BV608R includes BV609R

BV610R includes BV611R




General Surgery

Abdominal Retractors



BV615R-BV621R



BV632R-BV638R



BV622R-BV628R

BALFOUR Abdominal Retractor (continued)

Catalog Number	Description
BV615R	BALFOUR RETRACTOR, 200 x 255 MM COMPLETE W/ BV616R
BV632R	BALFOUR RETRACTOR, COMPLETE W/ BV606R, BV621R, BV622R
BV634R	BALFOUR RETRACTOR, COMPLETE W / BV609R, BV621R, BV624R
BV636R	BALFOUR RETRACTOR, COMPLETE W/ BV611R, BV621R, BV626R
BV638R	BALFOUR RETRACTOR, COMPLETE W/ BV611R, BV621R, BV628R
BV616R	CENTER BLADE, 107 x 59 MM
BV621R	SPREADER ONLY, 200 x 250 MM
BV622R	LATERAL BLADE, PAIR 50 MM
BV624R	LATERAL BLADE, PAIR 60 MM
BV626R	LATERAL BLADE, PAIR 75 MM
BV628R	LATERAL BLADE, PAIR 90 MM

**BALFOUR Side Blade Solid
F/ Retractor MB221R, MB222R**

Catalog Number	Description
MB637R	PAIR, 2 1/2" (64 MM) DEEP
MB638R	PAIR, 3 1/2" (89 MM) DEEP
MB639R	PAIR, 4" (102 MM) DEEP



MB637R



MB638R



MB639R

**BALFOUR Side Blade Wire
F/ Retractor MB221R, MB222R**


Catalog Number	Description
MB577R	PAIR, 2 1/2" (64 MM) DEEP
MB578R	PAIR, 3 1/2" (89 MM) DEEP
MB579R	PAIR, 4" (102 MM) DEEP



MB577R




MB578R



MB579R

General Surgery


Abdominal Retractors



MB221R-MB222R

BALFOUR Retractor W/ Exchangeable Lateral Blades

Catalog Number	Total Width (in.)	Opening Width (in.)	Description
MB221R	16	10 1/2	2 PAIRS OF DETACHABLE BLADES 2 1/2" AND 3 1/2", 1 CENTER BLADE 2 3/4" x 3 1/4"
MB222R	16	8 1/2	2 PAIRS OF DETACHABLE BLADES 2 1/2" AND 3 1/2", 1 CENTER BLADE 2 1/2" x 2 1/2"



MB225

BALFOUR Retractor

Catalog Number	Total Width (in.)	Opening Width (in.)	Description
MB225	13	11	2 LATERAL BLADES 2" x 2" AND 1 1/4" x 2", 1 CENTER BLADE 3" x 10 1/2"

General Surgery

Abdominal Retractors



MB573R-MB574R

DIXON Center Blade

Catalog Number	Blade	Blade Size (mm)
MB573R	MALLEABLE	171 x 51
MB574R	MALLEABLE	267 x 57



BV617R

BALFOUR Center Blade For All BALFOUR Retractors

Catalog Number	Blade Size (mm)
BV617R	54 x 80




BV618R


Center Blade F/ BV615R

Catalog Number	Blade Size (mm)
BV618R	75 x 80

General Surgery

Abdominal Retractors





BV646R

DENIS-BROWNE Abdominal Retractor

Catalog Number	Description
BV650R	COMPLETE SET
BV648R	2 BLADES ONLY, 40 x 30 MM
BV649R	2 BLADES ONLY, 40 x 40 MM
BV646R	BLADE ONLY, 20 x 30 MM

BV650R Includes:

1 Frame

2 Blades of each BV648R and BV649R


General Surgery

Abdominal Retractors

SEMM Abdominal Retractor

Catalog Number	Description	Blade Size (mm)
BV654R	COMPLETE SET	-
BV653R	FRAME	-
BV655R	2 BLADES	46 x 64
BV656R	BLADE	53 x 54
BV657R	BLADE	82 x 97
BV658R	BLADE	98 x 51
BV659R	BLADE	107 x 36


BV654R Includes:
1 Frame BV653R
1 Blade of each BV656R, BV657R,
BV658R and BV659R
2 Blades BV655R



General Surgery

Universal Ring Retractor System

■ For more information about Universal Ring Retractor System see our brochure number DOC380.




MB499R

Universal Ring Retractor System

Catalog Number	Description
MB499R	COMPLETE SET

Complete System Includes:	
MB497	Segmented Medium Oval Ring
MB516	Medium Oval Ring
MB506	Standard Table Clamp 18"
MB488R	Horizontal Flex Bar w/ 360° Rotation 19 1/4"
MB487R	Horizontal Fixed Bar 15 1/2"
MB504	Post Clamp for 1" Bar
MB503	Non-Tilting Blade Clamp (2 each)
MB510	Tilting Blade Clamp (6 each)
MB524	KELLY Retractor Blade 2" x 3" (2 each)
MB525	KELLY Retractor Blade 2" x 4"
MB526	KELLY Retractor Blade 2" x 5"
MB527	KELLY Retractor Blade 2" x 6"
MB529	Malleable Retractor Blade 1 1/2" x 6"
MB530	Malleable Retractor Blade 2" x 6"
MB532	Malleable Retractor Blade 3" x 6"
MB534	BALFOUR Bladder Blade 3" x 2 1/2"
MB544	HARRINGTON Sweetheart Blade 2 1/2" x 5"
MB546	GELPI Point Retractor Blade




MB515-MB517

Oval Ring

Catalog Number	Size	Width x Length (in.)
MB515	SMALL	8 1/2 x 12
MB516	MEDIUM	10 1/2 x 15
MB517	LARGE	12 1/2 x 18 1/2

General Surgery

Universal Ring Retractor System



MB496-MB462R

Segmented Oval Ring

Catalog Number	Size	Description	Width x Length (in.)
MB496	SMALL	INCLUDES REMOVABLE 4" SEGMENTS	8 1/2 x 14
MB497	MEDIUM	INCLUDES REMOVABLE 8" SEGMENTS	10 1/2 x 19 1/2
MB498	LARGE	INCLUDES REMOVABLE 6" SEGMENTS	12 1/2 x 19 1/2
MB462R	X-LARGE	INCLUDES REMOVABLE 6" SEGMENTS	16 x 22



MB471R-MB473R

Segments for Segmented Oval Ring, Pair

Catalog Number	Length (in.)
MB471R	4
MB472R	6
MB473R	8



MB564

URRS Segmented Ring Screw

Catalog Number
MB564


General Surgery

Universal Ring Retractor System

Round Ring

Catalog Number	Size	Diameter (in.)
MB519	SET	-
MB520	SMALL	8 1/2
MB521	MEDIUM	10 1/2
MB522	LARGE	12 1/2


MB519 Includes:
1 Ring of each MB520-MB522



MB519-MB522

Segmented Round Ring


Catalog Number	Size	Diameter (in.)
MB825R	SMALL	8 1/2
MB826R	MEDIUM	10 1/2
MB827R	LARGE	12 1/2
MB828R	X-LARGE	16



MB825R-MB828R

General Surgery

Universal Ring Retractor System



MB489

Vaginal Ring


Catalog Number	Description	Width x Length (in.)
MB489	20° ANGLE	8 x 10 1/2



MG511R

Tilt and Swivel Clamp

Catalog Number
MG511R



MB973R

Multi Directional Blade Clamp

Catalog Number
MB973R



MB503

Non-Tilting Blade Clamp


Catalog Number
MB503



MB510

Tilting Blade Clamp

Catalog Number
MB510



MB504

Post Clamp


Catalog Number
MB504

General Surgery

Universal Ring Retractor System

Table Clamp

Catalog Number	Working Length (in.)	Length (in.)
MB506	15	18
MB460R	21	24



MB460R-MB506

Flat Table Clamp

Catalog Number	Working Length (in.)	Total Length (in.)
MB670R	11 1/2	18
MB821R	17	24



MB670R-MB821R

Horizontal Fixed Bar

Catalog Number	Length (in.)
MB487R	15 1/2
MB823R	22



MB487R-MB823R

General Surgery

Universal Ring Retractor System

Horizontal Flex Bar W/ 360° Rotation

Catalog Number	Length (in.)
MB488R	19
MB824R	22 1/2



MB488R-MB824R

Vertical Extension Bar

Catalog Number	Length (in.)
MB508	12
MB822R	14



MB508-MB822R

Adapter For Attaching Handheld Retractors

Catalog Number
MB509




MB509



General Surgery

Universal Ring Retractor System



MB469R-MB555

KELLY Retractor Blade

Catalog Number	Description	Width x Depth (in.)
MB555	-	1 x 1 1/2
MB556	-	1 x 2 1/2
MB936R	-	2 x 1 1/2
MB937R	-	2 x 2 1/2
MB524	-	2 x 3
MB525	-	2 x 4
MB526	-	2 x 5
MB527	-	2 x 6
MB469R	BARIATRIC	2 x 8



MB938A-MB943A

Aluminum KELLY Retractor Blade

Catalog Number	Description	Width x Depth (in.)
MB938A	-	2 x 2 1/2
MB939A	-	2 x 3
MB940A	-	2 x 4
MB941A	-	2 x 5
MB942A	-	2 x 6
MB943A	BARIATRIC	2 x 8

General Surgery

Universal Ring Retractor System



MB467R-MB860R

Malleable Retractor Blade

Catalog Number	Description	Width x Depth (in.)
MB860R	-	3/4 x 2
MB861R	-	3/4 x 3
MB862R	-	1 x 1
MB863R	-	1 x 4
MB864R	-	1 x 5
MB865R	-	1 1/2 x 3
MB866R	-	1 1/2 x 4
MB529	-	1 1/2 x 6
MB867R	-	1 1/2 x 8
MB530	-	2 x 6 1/2
MB531	-	2 x 8
MB463R	BARIATRIC	2 x 10
MB532	-	3 x 6
MB464R	BARIATRIC	3 x 10
MB466R	BARIATRIC	3 x 12
MB869R	BARIATRIC	4 x 8
MB465R	BARIATRIC	4 x 10
MB467R	BARIATRIC	4 x 12



MB954A-MB959A

Aluminum Malleable Retractor Blade, Flat

Catalog Number	Description	Width x Depth (in.)
MB954A	-	1 1/2 x 6
MB955A	-	2 x 6
MB956A	-	2 x 8
MB957A	-	3 x 6
MB958A	BARIATRIC	2 x 10
MB959A	BARIATRIC	3 x 10

General Surgery

Universal Ring Retractor System



MB557-MB559

RICHARDSON Retractor


Catalog Number	Width x Depth (in.)
MB557	3/4 x 1
MB558	3/4 x 2
MB559	1 x 1 1/4



MB470R-MB534

BALFOUR Bladder Blade

Catalog Number	Description	Width x Depth (in.)
MB534	-	3 x 2 1/2
MB535	-	3 x 3 1/4
MB470R	BARIATRIC	3 x 3 1/2



MB829A-MB830A


Aluminum BALFOUR Bladder Blade

Catalog Number	Width x Depth (in.)
MB829A	2 x 2 3/8
MB830A	3 x 2 1/2

General Surgery

Universal Ring Retractor System


GENERAL
SURGERY



MB968R-MB969R

MAYO Retractor Blade

Catalog Number	Width x Depth (in.)
MB968R	1 1/4 x 1 1/4
MB969R	2 x 2



MB970A

Aluminum MAYO Retractor Blade

Catalog Number	Width x Depth (in.)
MB970A	2 3/4 x 2



MB537-MB539

DEAVER Retractor Blade

Catalog Number	Width x Depth (in.)
MB537	1 x 4 1/2
MB538	1 1/2 x 4 1/2
MB539	2 x 4 1/2



MB492-MB493


DEAVER Vaginal Retractor Blade

Catalog Number	Description	Width x Depth (in.)
MB492	SMALL	1 x 3 1/2
MB493	LARGE	1 x 5


INCHES

General Surgery

Universal Ring Retractor System



Vaginal Retractor Posterior Blade	
Catalog Number	Width x Depth (in.)
MB490	1 x 4 1/2



Vaginal Retractor Lateral Blade	
Catalog Number	Width x Depth (in.)
MB491	1 x 3




ST. MARK'S Deep Pelvic Retractor Blade	
Catalog Number	Width x Depth (in.)
MB871R	2 3/8 x 7 1/2



MB494-MB495

Rectal Blade

Catalog Number	Description	Width x Depth (in.)
MB494	SMALL	7/8 x 2 7/8
MB495	LARGE	7/8 x 3 7/8




MB541

Skin Incision Blade

Catalog Number	Width x Depth (in.)
MB541	1 1/8 x 3 1/4

General Surgery


Universal Ring Retractor System



MB543-MB832R

HARRINGTON Sweetheart Blade

Catalog Number	Description	Width x Depth (in.)
MB543	-	1 1/2 x 5
MB544	-	2 1/2 x 5 5/8
MB468R	BARIATRIC	3 1/2 x 10
MB832R	BARIATRIC	4 1/2 x 10



MB833A-MB835A

Aluminum HARRINGTON Sweetheart Blade

Catalog Number	Description	Width x Depth (in.)
MB833A	-	2 1/2 x 4
MB834A	-	3 1/2 x 8
MB835A	BARIATRIC	3 1/2 x 10




MB546

GELPI Point Retractor Blade

Catalog Number
MB546

General Surgery

Universal Ring Retractor System



MB548-MB549

SCOVILLE Hook Retractor Blade

Catalog Number	Description	Depth (in.)
MB549	SMALL	2
MB548	LARGE	3 1/4



MB566-MB568

Rake Retractor

Catalog Number	Description	Width x Depth (in.)
MB566	3 PRONGS	2/3 x 3/4
MB567	4 PRONGS	3/4 x 1/3
MB568	6 PRONGS	1 x 1/3




MB551-MB553

MEYERDING Retractor Blade

Catalog Number	Width x Depth (in.)
MB551	1 x 3 1/3
MB552	2 x 3 3/4
MB553	2/3 x 2 1/3

General Surgery

Universal Ring Retractor System



MB962A-MB965A

Aluminum Renal Retractor Blade

Catalog Number	Width x Depth (in.)
MB962A	1 x 4
MB963A	1 x 6
MB964A	1 x 8
MB965A	1 x 10




MB966A-MB967A

Aluminum Splanchnic Retractor Blade

Catalog Number	Width x Depth (in.)
MB966A	2 x 10
MB967A	3 x 8

Perforated Super Malleable Retractor
Blade


Catalog Number	Width x Depth (in.)
MB960R	2 x 9
MB961R	3 x 9



MB960R-MB961R

YU-HOLTGREWE Malleable Retractor
Blade

Catalog Number	Width x Depth (in.)
MB870R	2 1/2 x 6



MB870R



General Surgery

Self Retaining Retractors



BV001R




BV001R

FINSEN Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV001R	SHARP	2 x 3 PRONGS	2



BV002R



BV002R

FINSEN Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV002R	SHARP	3 x 4 PRONGS	2 3/4



BV003R



BV003R

FINSEN Modified Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV003R	BLUNT	2 x 3 PRONGS	3



BV005R



BV005R

FINSEN Retractor


Catalog Number	Tip	Blade End	Length (in.)
BV005R	SHARP	3 x 4 PRONGS	2 3/4



BV014R

PARSONNET Epicardial Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV014R	SHARP	3 x 3 PRONGS	2






BV006R-BV007R



BV006R




BV007R



MD471

Spring Wire Retractor

Catalog Number	Tip	Prongs	Length (in.)
BV006R	BLUNT	3 x 3	1 3/4
BV007R	BLUNT	3 x 3	2 1/4
MD471	SHARP	1 x 1	3



MB196R-MB777R

Spring Wire Retractor

Catalog Number	Size	Length (in.)
MB777R	SMALL	5
MB778R	MEDIUM	5 1/8
MB196R	LARGE	9

General Surgery

Self Retaining Retractors



BV010R-BV011R



BV010R



BV011R

ALM Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV010R	SHARP	4 x 4 PRONGS	2 3/4
BV011R	SHARP	4 x 4 PRONGS	4



BV008R-BV015R




BV008R



BV015R

HEISS Self-Retaining Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV015R	SHARP	4 x 4 PRONGS	4
BV008R	BLUNT	4 x 4 PRONGS	4









General Surgery

Self Retaining Retractors


WEITLANER Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV073R	BLUNT	2 x 3 PRONGS	4 1/4
MA438R	BLUNT	3 x 3 PRONGS	4 1/4
BV074R	BLUNT	3 x 4 PRONGS	5 1/8
BV075R	BLUNT	3 x 4 PRONGS	6 1/2
BV104R	BLUNT/DEEP	3 x 4 PRONGS	5 1/8



General Surgery


Self Retaining Retractors



BV070R




BV070R




BV071R




BV071R




BV072R



BV072R



BV101R




BV101R




BV066R




BV066R




BV067R



BV067R



BV068R



BV068R

WEITLANER Retractor




Catalog Number	Tip	Blade End	Length (in.)
BV070R	SHARP	2 x 3 PRONGS	4 1/4
BV071R	SHARP	3 x 4 PRONGS	5 1/8
BV072R	SHARP	3 x 4 PRONGS	6 1/2
BV101R	SHARP/DEEP	3 x 4 PRONGS	5 1/8
BV066R	SEMI-S	2 x 3 PRONGS	4 1/4
BV067R	SEMI-S	3 x 4 PRONGS	5 1/8
BV068R	SEMI-S	3 x 4 PRONGS	6 1/2

General Surgery

Self Retaining Retractors


WULLSTEIN Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV076R	SHARP	3 x 3 PRONGS	5
BV077R	BLUNT	3 x 3 PRONGS	5




Baby ADSON Retractor, Hinged

Catalog Number	Tip	Blade End	Length (in.)
BV085R	SEMI-S	3 x 4 PRONGS	5 $\frac{1}{2}$
BV094R	SEMI-S	3 x 4 PRONGS	6 $\frac{1}{2}$



Baby ADSON Retractor, Angled


Catalog Number	Tip	Blade End	Length (in.)
MG007R	BLUNT	4 x 4 PRONGS	7 $\frac{1}{2}$



MG007R














General Surgery

Self Retaining Retractors



WEITLANER Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV200R	SHARP	3 x 4 PRONGS	7 7/8
BV206R	SHARP	3 x 4 PRONGS	9 1/2
BV208R	SHARP	5 x 6 PRONGS	10
BV201R	BLUNT	3 x 4 PRONGS	7 7/8
BV207R	BLUNT	3 x 4 PRONGS	9 1/2
BV209R	BLUNT	5 x 6 PRONGS	10




MISKIMON Cerebellar Retractor

Catalog Number	Tip	Blade End	Length (in.)
MB476R	SHARP	7 x 7 PRONGS	7 7/8




Thyroid Retractor

Catalog Number
BV086R




JOLL (FRIEDMANN) Thyroid/Vaginal Retractor

Catalog Number	Length (in.)
BV087R	6 1/8




MOLLISON Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV088R	SHARP	4 x 4 PRONGS	6 1/8



General Surgery


Self Retaining Retractors




BV112R

MAYO-ADAMS Retractor

Catalog Number	Length (in.)
BV112R	6 1/2




BV112R 1/1




MB929R

WILSTE-GELPI Retractor

Catalog Number	Tip	Tip Size (mm)	Length (in.)
MB929R	SHARP	70	11



MB929R



MB199R

GELPI Perineal Retractor

Catalog Number	Tip	Length (in.)
MB199R	SHARP WITH BALL STOPS	7




MB199R

General Surgery


Self Retaining Retractors

GELPI Vaginal Retractor

Catalog Number	Tip	Length (in.)
MG176R	SHARP	3 1/2
BV996R	SHARP	5 1/4
BV997R	SHARP	6 7/8
MD651	SHARP	7 1/4
MG177R	BLUNT	4 1/2



MG176R-MG177R




MG176R




MG177R




BV996R




BV996R



BV997R



BV997R



MD651




MD651

General Surgery

Self Retaining Retractors




BV998R



BV998R

DEBAKEY Femur Retractor for Children

Catalog Number	Tip	Length (in.)
BV998R	SHARP	4




BV999R




BV999R

DEBAKEY Femur Retractor for Adults

Catalog Number	Tip	Length (in.)
BV999R	SHARP	5 1/2



BV993R




BV993R

JACKSON-BURROWS Retractor

Catalog Number	Tip	Length (in.)
BV993R	BLUNT	7

General Surgery

Self Retaining Retractors



HENLY Retractor

Catalog Number	Description	Blade Size (mm)	Blade Size (in.)	Length (in.)
BV120R	COMPLETE SET	-	-	7
MB842R	CENTER BLADE, SMALL	19 x 16	3/4 x 5/8	-
MB843R	CENTER BLADE, MEDIUM	25 x 16	1 x 5/8	-
MB844R	CENTER BLADE, LARGE	32 x 16	1 1/4 x 5/8	-
MB848R	CENTER BLADE, X-LARGE	76 x 16	3 x 5/8	-
MB845R	SIDE BLADE, PAIR (4 PRONGS)	25 x 25	1 x 1	-
MB846R	SIDE BLADE, PAIR (4 PRONGS)	51 x 25	2 x 1	-
MB847R	SIDE BLADE, PAIR (4 PRONGS)	76 x 22	3 x 1	-
MB849R	FRAME ONLY	-	-	7


BV120R Complete Set Includes:

1 Frame,


1 each of 4 Center Blades,

1 Pair of each of 3 Side Blades


Materials in set sold separately: MB842R-MB849R



MB842R



MB843R



MB844R



MB848R




MB845R



MB846R




MB847R




MB849R

General Surgery


Rectal Specula Instruments




BV585R




BV581R




BV583R



BV586R



BV584R




BV589R

PARKS Rectal Spreader

Catalog Number	Item	Length (mm)
BV585R	COMPLETE SET	-
BV581R	SPREADER ONLY	-
BV582R	CENTER BLADE	-
BV583R	SIDE BLADE, PAIR	70
BV586R	SIDE BLADE, PAIR	90
BV584R	SIDE BLADE, PAIR	93
BV589R	SIDE BLADE, PAIR	150


Set Includes:
1 Spreader BV581R
1 Center Blade BV582R
2 Blades of each BV583R and BV584R



BV906R


SMITH-BUIE Rectal Retractor

Catalog Number	Blade Size (mm)
BV906R	58 x 23




ST. MARKS MODIFIED GIRONA Perineal Retractor

Catalog Number
BV907R



CZERNY Rectal Speculum

Catalog Number	Length (in.)
EA800R	8 7/8




SAWYER Retractor

Catalog Number	Size	Blade Size (mm)	Length (in.)
MB773R	SMALL	70 x 19	10 3/4
MB772R	MEDIUM	82 x 29	10 3/4
MG178R	LARGE	97 x 38	11 5/8

General Surgery

Rectal Specula Instruments



MB774R-MB776R

HILL-FERGUSSON Rectal Retractor

Catalog Number	Blade Size (mm)	Length (in.)
MB775R	64 x 24	8 1/4
MB774R	76 x 35	8 1/4
MB776R	76 x 32	8 1/4




EA820R



EA821R

SIMS Rectal Speculum

Catalog Number	Blade	Length (in.)
EA820R	SOLID	6 1/4
EA821R	FENESTRATED	6 1/4



EA831R

PRATT Rectal Speculum


Catalog Number	Length (in.)
EA831R	8

General Surgery

Rectal Specula Instruments

MATHIEU Rectal Speculum

Catalog Number	Length (in.)
EA843R	8 1/2



EA843R

IVES-FANSLER Anuscope

Catalog Number	Length (in.)
EA880R	5 3/4



EA880R

Rectal Speculum

Catalog Number	Length (in.)
MB274R	2 3/8



MB274R

General Surgery


Rectal Specula Instruments



MG160R

CHELSEA-EATON Anal Speculum

Catalog Number	Size	Diameter (mm)	Finish	Length (in.)
MG160R	SMALL	30	MIRROR FINISH	2 1/2



MB179R, MB203R, MB288R

HIRSCHMANN Rectal Speculum

Catalog Number	Size	Blade Size (mm)
MB203R	SMALL	60 x 16
MB288R	MEDIUM	72 x 19
MB179R	LARGE	89 x 23

General Surgery

Rectal Specula Instruments



MB284R, MG162R-MG163R

BRINKERHOFF Rectal Speculum

Catalog Number	Size	Diameter (in.)	Length (in.)
MG162R	SMALL	1/2 x 1	3 3/4
MG163R	MEDIUM	5/8 x 1	3 3/4
MB284R	LARGE	3/4 x 1 1/4	4 1/2




EA882R-EA887R

KELLY Sphincteroscope

Catalog Number	Diameter (mm)	Working Length (mm)	Length (in.)
EA882R	21	54	6 1/4
EA887R	22	140	5 1/2
EA883R	25	50	6 1/4

General Surgery

Rectal Specula Instruments



MCGIVNEY Hemorrhoidal Forceps


Catalog Number	Length (in.)
EA990R	7 1/4



RUDD Hemorrhoidal Ligator

Catalog Number	Description	Length (in.)
EA991R	COMPLETE SET	12 3/8
EA992R	LIGATOR, INCLUDES CONE	12 1/2
EA992208	LIGATOR CONE ONLY	-
EA993	LATEX FREE RINGS F/ EA991R, 100/PKG	-

EA991R Complete Set Includes:
EA992R and EA993



General Surgery
Rectal Specula Instruments

GENERAL
SURGERY



STEWART Crypt Hook

Catalog Number	Length (in.)
MB892R	9



INCHES

General Surgery


Probes, Dilators and Scoops

BAKES Bile-Duct Dilator

Catalog Number	Tip Size (mm)	Length (in.)
EB061R	1	12 3/8
EB062R	2	12 3/8
EB063R	3	12 3/8
EB064R	4	12 3/8
EB065R	5	12 3/8
EB066R	6	12 1/2
EB067R	7	12 1/2
EB068R	8	12 1/2
EB069R	9	12 1/2
EB070R	10	12 1/2
EB071R	11	12 1/2
EB072R	12	12 1/2
EB073R	13	12 1/2
EB074R	SET, 3-11 mm	-



EB061R-EB074R




EB061R



EB065R



EB069R



EB062R



EB066R



EB070R



EB063R



EB067R



EB071R



EB064R



EB068R



EB072R




EB073R

LUER-KOERTE Scoops, Malleable Shaft


Catalog Number	Tip Size (mm)	Length (in.)
EB168R	2.8	12 3/4
EB169R	3.4	12 3/4
EB170R	4.3	12 3/4
EB171R	5.5	12 3/4
EB172R	6.7	12 3/4
EB173R	8.0	12 3/4
EB174R	11.0	12 3/4
EB175R	15.0	12 3/4




EB168R-EB175R




EB168R




EB169R



EB170R



EB171R



EB172R



EB173R




EB174R



EB175R



BN003R-BN006R




BN003R-BN006R

Director W/ Tongue

Catalog Number	Blade Shape	Length (in.)
BN003R	STR	5 1/8
BN004R	STR	5 3/4
BN006R	STR	6 1/4



BN010R



BN010R

DOYEN Director W/ Tongue

Catalog Number	Blade Shape	Length (in.)
BN010R	CVD	5 3/4

General Surgery


Probes, Dilators and Scoops



MB336R

Grooved Director W/ Tongue


Catalog Number	Length (in.)
MB336R	4 1/2




BN012R

NELATON Director

Catalog Number	Blade Shape	Length (in.)
BN012R	CVD	6 1/4




BN012R




BN015R

Director, W/ Probe Tip & Tongue

Catalog Number	Diam. of Tip (mm)	Blade Shape	Length (in.)
BN015R	3.3	STR	5 1/2




BN015R




BN020R

BRODIE Director, W/ Probe Tip & Tongue

Catalog Number	Diam. of Tip (mm)	Blade Shape	Length (in.)
BN020R	3.5	STR	7 7/8



BN020R







BN170R

Probe

Catalog Number	Diam. of Tip (mm)	Length (in.)
BN170R	1.4	6 1/4

LOCKHART-MUMMERY Fistula Probe

Catalog Number	Blade Shape	Length (in.)
BN021R	#1-STR	6 1/2
BN022R	#2-SLT-CVD	6 1/2
BN023R	#3-CVD	6 1/2
BN024R	#4-HVY-CVD	6 1/2



INCHES

General Surgery

Probes, Dilators and Scoops

MG168R-MG169R



BUIE Fistula Probe

Catalog Number	Description	Length (in.)
MG168R	LIGHT	6 1/2
MG169R	HEAVY	6 1/2


BN034R-BN040R, MB603R



Probe W/ Eye

Catalog Number	Tip Size (mm)	Length (in.)
MB603R	1.5	5
BN034R	2.5	5 3/4
BN036R	2.5	6 1/4
BN040R	3	7 7/8

MG571



Probe W/ Eye, Silver

Catalog Number	Tip Size (mm)	Length (in.)
MG571	1.75	4 1/2

BN054R



BN063R-BN070R



Myrtle Leaf Probe

Catalog Number	Tip Size (mm)	Length (in.)
BN054R	1.5	5 3/4
BN063R	2	5 1/8
BN064R	2	5 3/4
BN066R	2	6 1/4
BN068R	2	7
BN070R	2	7 7/8



BN104R-BN150R

Double Ended Probe

Catalog Number	Tip Size (mm)	Length (in.)
BN104R	1.0	5 3/4
BN113R	1.5	5 1/8
BN114R	1.5	5 3/4
BN116R	1.5	6 1/4
BN118R	1.5	7
BN120R	1.5	7 7/8
BN125R	1.5	9 3/4
BN134R	2.0	5 3/4
BN136R	2.0	6 1/4
BN138R	2.0	7
BN140R	2.0	7 7/8
BN150R	2.0	11 3/4



BN090R

Double Ended Probe W/ Eye, Malleable


Catalog Number	Tip Size (mm)	Length (in.)
BN090R	3.0	7 1/2
BN095	3.5	9 3/4
BN096	3.5	11 3/4



BN095-BN096

PRATT Rectal Probe, Malleable

Catalog Number	Tip Size (mm)	Length (in.)
MB486	2.0	10



MB486

Probe

Catalog Number	Tip Size (mm)	Length (in.)
BN101R	2.0	11 3/4



BN101R



General Surgery

Probes, Dilators and Scoops

A long, thin surgical instrument with a double-ended design. It features a long, straight handle with a textured grip in the middle. At each end, there is a shallow, spoon-like scoop. The scoops are slightly curved and have a smooth, polished surface.

FERGUSSON Gall Stone Scoop, Double Ended

Catalog Number	Size	Length (in.)
MB697R	SMALL	9 1/2
MB698R	MEDIUM	9 1/2
MB699R	LARGE	9 1/2

MB697R-MB699R

1/2

1/1

1/1

MB697R-MB699R

A long, thin surgical instrument with a double-ended design. It features a long, straight handle with a textured grip in the middle. At each end, there is a thin, curved probe. The probes are slightly curved and have a smooth, polished surface.

BARR Fistula Probe, Double Ended

Catalog Number	Length (in.)
MG170R	6 5/8

MG170R

1/2


A long, thin surgical instrument with a single-ended design. It features a long, straight handle with a textured grip in the middle. At one end, there is a thin, curved probe. The probe is slightly curved and has a smooth, polished surface.

PRATT Fistula Probe, Crypt Hook

Catalog Number	Length (in.)
EA890R	9 1/2


EA890R

1/2




ROSSER Crypt Hook

Catalog Number	Length (in.)
MG165R	6 5/8



DESJARDINS Gall Duct Probe, Flexible

Catalog Number	Tip Size (mm)	Length (in.)
EB101R	3.8	11 1/4
EB102R	4.3	11 1/4
EB103R	6.8	11 1/4



DESJARDINS Gall Stone Scoop

Catalog Number	Tip Size (mm)	Length (in.)
EB111R	6	11 3/4
EB112R	7	11 3/4
EB113R	8	11 3/4
EB114R	9	11 3/4



MAYO Gall Duct Probe

Catalog Number	Tip Size (mm)	Tip	Length (in.)
E0005R	5	MALLEABLE	13

General Surgery

Scalpel Handles


Micro Scalpel Handle

Catalog Number	Size	Length (in.)	Handle
MB104R	3H	4	Hexagonal
MB141R	3HL	6	Hexagonal
MB280R	3KM	3	Round Knurled
MB997R	3K	4	Round Knurled
MB281R	3KD	5	Round Knurled
BB046R	3K	5 1/4	Round Knurled
MB271R	3KL	6	Round Knurled
BB048R	3KL	6 1/8	Round Knurled
MG221R	3	7	Round Knurled
BB054R	3	7 7/8	Round Knurled
MG222R	3	8	Round Knurled
BB060R	3	17 3/4	Round Knurled
MB475R	COLLET		

MB475R
Replacement collet for MB104R, MB141R, MB271R, MB280R, MB281R and MB997R



MB104R



MB141R



MB280R



MB997R



MB281R



BB046R



MB271R



BB048R




MG221R



BB054R



MG222R



BB060R



MB475R

General Surgery

Scalpel Handles

CORWIN Knife Handle

Catalog Number	Length (in.)
MA203R	5 1/2



MA203R

CASPAR Scalpel Handle

Catalog Number	Length (in.)
BB052R	7



BB052R

LANDOLT Scalpel Handle, F/ Microsurgery Blades


Catalog Number	Length (in.)
BB053R	8 1/4
BB057R	9 3/4



BB053R



BB053R



BB057R



BB331R



BB332R

YASARGIL Arachnoid Blades F/ LANDOLT Scalpel Handles


Catalog Number	Tip Size (mm)	Packaging
BB331R	1.5	10/BOX
BB332R	2.5	10/BOX

General Surgery

Scalpel Handles


Transssphenoidal Scalpel Handle,
Bayonet, F/ Vertical and Lateral Blades

Catalog Number	Length (in.)
BB071R	8 1/4



Sterile Microsurgery Blades

Catalog Number	Tip	Packaging
BB361R	SMALL CHISEL #61	10/BOX
BB362R	WIDE CHISEL #62	10/BOX
BB363R	SPEAR #63	1 EACH
BB364R	ROUND #64	10/BOX
BB365R	SMALL, SHARP #65	10/BOX
BB367R	CURVED, POINTED #67	10/BOX
BB368R	SMALL SICKLE #68	10/BOX



BARRON Scalpel Handle #3 F/ Blades
nos. 10-15, 40, 42

Catalog Number	Length (in.)
BB068R	5 1/8



Scalpel Handle #3, F/ Blades
nos. 10-15, 40, 42

Catalog Number	Length (in.)
BB070R	6 1/8



General Surgery

Scalpel Handles

GENERAL
SURGERY

Scalpel Handle #3, F/ Blades
nos. 10-15, 40, 42

Catalog Number	Size	Length (in.)
BB074R	#3	5
BB075R	#3L	8 1/4
BB175R	#3XL	9 3/4
MB917R	#3	16



Scalpel Handle #3L, F/ Blades
nos. 10-15, 40, 42

Catalog Number	Blade Shape	Length (in.)
MB127R	ANG	8



Scalpel X-Long Handle, #3, F/ Blades
nos. 10-15, 40, 42

Catalog Number	Blade Shape	Length (in.)
BB078R	STR	8 7/8

INCHES

General Surgery

Scalpel Handles




Scalpel X-Long Handle, #3, F/ Blades nos. 10-15, 40, 42

Catalog Number	Blade Shape	Length (in.)
BB079R	ANG	8 7/8



Scalpel Handle #4, F/ Blade nos. 18-36

Catalog Number	Length (in.)
BB080R	8 7/8



Scalpel Handle #4, F/ Blades nos. 18-36

Catalog Number	Length (in.)
BB084R	5 1/4



Scalpel Handle #4L, F/ Blades nos. 18-36

Catalog Number	Length (in.)
BB085R	8 1/2

Scalpel Handle #4XL, F/ Blades
nos. 18-36

Catalog Number	Length (in.)
BB087R	9 3/4



BB087R

2/3

Scalpel Handle #7, F/ Blades
nos. 10-15, 40, 42

Catalog Number	Length (in.)
BB077R	6 1/4



BB077R

1/1

Scalpel Handle #9XL, F/ Blades
nos. 18-36

Catalog Number	Length (in.)
MA202R	4 5/8
BB089R	5 3/4



MA202R

1/1



BB089R

1/1

General Surgery

Scalpel Blades



BB510, 16600495



BB515, 16600525




BB511, 16600509



BB540 double cutting



BB512, 16600517



BB542 double cutting

Sterile Scalpel Blade, For #3 Handle

Catalog Number		Blade Number	Packaging
Carbon	Stainless Steel		
BB510	16600495	#10	100/BOX
BB511	16600509	#11	100/BOX
BB512	16600517	#12	100/BOX
BB515	16600525	#15	100/BOX
BB540	-	#40	100/BOX
BB542	-	#42	100/BOX



BB520, 16600550



BB521, 16600568



BB522, 16600576



BB523, 16600584



BB524, 16600592



BB525, 16616634



BB536, 16615522

Sterile Scalpel Blade, For #4 Handle

Catalog Number		Blade Number	Packaging
Carbon	Stainless Steel		
BB520	16600550	#20	100/BOX
BB521	16600568	#21	100/BOX
BB522	16600576	#22	100/BOX
BB523	16600584	#23	100/BOX
BB524	16600592	#24	100/BOX
BB525	16616634	#25	100/BOX
BB536	16615522	#36	100/BOX




BB500, 16600487

Sterile Stitch Cutter Blade #3 Handle

Catalog Number		Packaging
Carbon	Stainless Steel	
BB500	16600487	100/BOX

General Surgery

Sterile Scalpels With Plastic Handle



Sterile Scalpels W/ Plastic Handle and Blade Protector

Catalog Number		Size	Packaging
Carbon	Stainless Steel		
BA210	5518059	10	10/BOX
BA211	5518040	11	10/BOX
BA212	5518105	12	10/BOX
BA215	5518032	15	10/BOX
BA218	-	18	10/BOX
BA219	-	19	10/BOX
BA220	5518067	20	10/BOX
BA221	5518075	21	10/BOX
BA222	5518083	22	10/BOX
BA223	5518016	23	10/BOX
BA224	5518091	24	10/BOX

General Surgery

Amputation and Autopsy Instruments

LISTON Amputating Knife

Catalog Number	Blade Length (mm)	Blade Size (in.)
FH100R	130	5 1/8
FH101R	160	6 1/4
FH102R	190	7 1/2
FH103R	220	8 3/4



FH100R-FH103R

CATLIN Amputating Knife, Double Cutting

Catalog Number	Blade Length (mm)	Blade Size (in.)
FH116R	160	6 1/4




FH116R

1/2

CHARRIÈRE Amputation Saw, 2 Blades

Catalog Number	Item	Blade Size (in.)	Blade Width (mm)	Total Length (in.)
FH300R	SAW COMPLETE W/2 BLADES	8 1/2	-	14
FH301C	CARBON STEEL BLADE	8 1/2	12	-
FH302C	CARBON STEEL BLADE	8 1/2	18	-



FH300R



FH301C-FH302C

CHARRIÈRE Amputation Saw


Catalog Number	Item	Blade Size (in.)	Total Length (in.)
FH320R	SAW COMPLETE	7 1/4	10 1/2



FH320R

SATTERLEE Amputation Saw

Catalog Number	Item	Total Length (in.)
FH325R	SAW COMPLETE	11 1/2
FH326R	CARBON STEEL BLADE F/ FH325R	-



FH325R-FH326R

General Surgery


Amputation and Autopsy Instruments



FH335R

LANGENBECK Metacarpal Saw

Catalog Number	Length (in.)
FH335R	9 1/4



FH414-FH415, FH430R



FH414-FH415, FH430R

GIGLI Wire Saw

Catalog Number	Diameter (mm)	Wires	Sales Unit	Length (in.)
FH414	1.4	6	BOX OF 5	15 3/4
FH415	1.4	6	BOX OF 5	19 3/4
FH430R	1.8	4	BOX OF 5	19 3/4



MB986R

GIGLI Saw Handle

Catalog Number	Packaging
MB986R	2/SET



FH480R

GIGLI Hook Handle, F/ Wire Saws

Catalog Number	Packaging
FH480R	2/SET

POOLE Suction Tube

Catalog Number	Handle	Tip Size (FR)	Hole in Tip Inner/Outer	Tube Diameter Inner/Outer	Length (in.)
MG272R	STR	30	YES/YES	5 mm/10 mm	8 1/2
GF862R	STR	30	YES/NO	5 mm/10 mm	8 7/8
MB927R	STR	30	YES/YES	4 mm/10 mm	14
MB928R	STR	30	YES/YES	4 mm/10 mm	16
MG271R	CVD	24	YES/NO	4 mm/8 mm	8 3/4
GF860R	CVD	21	YES/YES	4 mm/8 mm	8 7/8



MG272R



GF862R




MB927R-MB928R



MB927R



MB928R



MG271R



GF860R



General Surgery

Suction Tubes



GF917R

ADSON Suction Tube

Catalog Number	Tip Size (mm)	Tip Size (FR)	Length (in.)
GF917R	4	11	7 1/2



GF919R

Nasal Suction Tube W/ Luer Hub, Tapered

Catalog Number	Tip Size (mm)	Length (in.)
GF919R	5	7



GF821R

COOLEY Suction Tube

Catalog Number	Tip Size (mm)	Length (in.)
GF821R	8.0	13 1/4


General Surgery

Suction Tubes

GENERAL
SURGERY

COOLEY Suction Tube W/ Detachable Tip

Catalog Number	Tip Size (mm)	Shaft	Length (in.)
GF822R	8.0	SLT-CVD	12 1/4
GF823R	9.0	CVD	11 5/8
GF824R	8.0	ANG	12 1/2
GF825R	9.0	90°	10 1/2



INCHES

General Surgery

Suction Tubes

Grafting Suction Tube

Catalog Number	Tip Size (mm)	Length (in.)
MB286R	3	9 1/2
MB287R	4	9 1/2



MB286R-MB287R

DEBAKEY Suction Tube

Catalog Number	Tip Size (mm)	Length (in.)
GF946R	3	10 3/4
GF947R	5	6 1/4
GF948R	5	10 5/8
GF946205	REPLACEMENT TIP FOR GF946R	-



GF946R



Diam. 9 mm

GF946R




GF947R



Diam. 5 mm

GF947R



GF948R



Diam. 5 mm


GF948R



MF333R-MF347R

Atraumatic Suction Tube

Catalog Number	Tip Size (FR)	Working Length (in.)	Length (in.)
MF333R	3	3	6 1/4
MF334R	5	3	6 1/4
MF335R	7	3	6 1/4
MF336R	10	3	6 1/4
MF337R	12	3	6 1/4
MF338R	3	4	7 1/4
MF339R	5	4	7 1/4
MF340R	7	4	7 1/4
MF341R	10	4	7 1/4
MF342R	12	4	7 1/4
MF343R	3	5	8 1/4
MF344R	5	5	8 1/4
MF345R	7	5	8 1/4
MF346R	10	5	8 1/4
MF347R	12	5	8 1/4



MF106R-MF111R



MF106R-MF108R



MF109R-MF111R

Atraumatic Suction Tube


Catalog Number	Tip	Tip Size (FR)	Length (in.)
MF106R	90° DOWN	5	8 1/2
MF107R	90° DOWN	7	8 1/2
MF108R	90° DOWN	10	8 1/2
MF109R	90° UP	5	8 1/2
MF110R	90° UP	7	8 1/2
MF111R	90° UP	10	8 1/2


General Surgery

Suction Tubes

YANKAUER Suction Tube

Catalog Number	Description	Tip Opening (mm)	Length (in.)
GF944R	-	2.0	11 1/4
OM671R	-	4.8	11 1/4
MB919R	-	2.0	12
MG092R	W/ VALVE	2.0	13 1/4
MB051R	-	2.0	14
MB920R	W/ VALVE	2.0	16
GF944200	REPLACEMENT TIP F/ GF944R	-	-
OM671200	REPLACEMENT TIP F/ OM671R	-	-
MB051200	REPLACEMENT TIP F/ MB051R, MB053R, MG919R, MG920R	-	-
MG092200	REPLACEMENT TIP F/ MG092R	-	-






MB052R

YANKAUER Suction W/ Thumb Control

Catalog Number	Description	Length (in.)
MB052R	HANDLE W/ THUMB CONTROL	6 7/8
MB053R	SUCTION TUBE, CHROME PLATED	6 5/8



MB053R




GF933R

FERGUSON Suction Tube

Catalog Number	Tip	Tip Size (mm)	Length (in.)
GF933R	ROUND	4	7

General Surgery


Suction Tubes



GF350R-GF385R

FERGUSSON Suction Tube, Luer Hub

Catalog Number	Diameter (mm)	Working Length (in.)	Length (in.)
GF350R	1.5	4 1/4	7 3/4
GF351R	2.0	4 1/4	7 3/4
GF352R	2.5	4 1/4	7 3/4
GF353R	3.0	4 1/4	7 3/4
GF354R	4.0	4 1/4	7 3/4
GF355R	5.0	4 1/4	7 3/4
GF360R	1.5	5 1/8	8 1/2
GF361R	2.0	5 1/8	8 1/2
GF362R	2.5	5 1/8	8 1/2
GF363R	3.0	5 1/8	8 1/2
GF364R	4.0	5 1/8	8 1/2
GF365R	5.0	5 1/8	8 1/2
GF370R	1.5	6 1/4	9 5/8
GF371R	2.0	6 1/4	9 5/8
GF372R	2.5	6 1/4	9 5/8
GF373R	3.0	6 1/4	9 5/8
GF374R	4.0	6 1/4	9 5/8
GF375R	5.0	6 1/4	9 5/8
GF380R	1.5	7 7/8	11
GF381R	2.0	7 7/8	11
GF382R	2.5	7 7/8	11
GF383R	3.0	7 7/8	11
GF384R	4.0	7 7/8	11
GF385R	5.0	7 7/8	11




MD124R-MD125300

MA501R, MA501300, MD127R, MD127300

FRAZIER Suction Tube

Catalog Number	Tip Size (FR)	Description	Working Length (in.)
MD124R	7		4 1/4
MD124300	7	SUCTION STYLET W/ OLIVE TIP FOR MD124R	
MD125R	7		4 1/4
MD125300	8	SUCTION STYLET W/ OLIVE TIP FOR MD125R	
MA501R	9		4 1/4
MA501300		SUCTION STYLET W/ OLIVE TIP FOR MA501R	
MD127R	10		4 1/4
MD127300	10	SUCTION STYLET W/ OLIVE TIP FOR MD127R	
MA499R	11		4 1/4
MA499300	11	SUCTION STYLET W/ OLIVE TIP FOR MA499R	
MD129R	12		4 1/4
MD129300	12	SUCTION STYLET W/ OLIVE TIP FOR MD129R	




MA501300

MA499R, MA499300, MD129R, MD129300



MF479R-MF489R, MF499R-MF508R




MF489R-MF498R, MF509R-MF518R

FRAZIER Slotted Suction Tube

Catalog Number	Tip Size (FR)	Shaft	Length (in.)
MF479R	3	STR	5 1/2
MF480R	4	STR	5 1/2
MF481R	5	STR	5 1/2
MF482R	6	STR	5 1/2
MF483R	7	STR	5 1/2
MF484R	8	STR	5 1/2
MF485R	9	STR	5 1/2
MF486R	10	STR	5 1/2
MF487R	11	STR	5 1/2
MF488R	12	STR	5 1/2
MF499R	3	STR	7
MF500R	4	STR	7
MF501R	5	STR	7
MF502R	6	STR	7
MF503R	7	STR	7
MF504R	8	STR	7
MF505R	9	STR	7
MF506R	10	STR	7
MF507R	11	STR	7
MF508R	12	STR	7
MF489R	3	ANG	5 1/2
MF490R	4	ANG	5 1/2
MF491R	5	ANG	5 1/2
MF492R	6	ANG	5 1/2
MF493R	7	ANG	5 1/2
MF494R	8	ANG	5 1/2
MF495R	9	ANG	5 1/2
MF496R	10	ANG	5 1/2
MF497R	11	ANG	5 1/2
MF498R	12	ANG	5 1/2
MF509R	3	ANG	7
MF510R	4	ANG	7
MF511R	5	ANG	7
MF512R	6	ANG	7
MF513R	7	ANG	7
MF514R	8	ANG	7
MF515R	9	ANG	7
MF516R	10	ANG	7
MF517R	11	ANG	7
MF518R	12	ANG	7

General Surgery

Suction Tubes



MD091-MD099

FRAZIER Long-Suction Tube

Catalog Number	Tip Size (FR)	Length (in.)
MD091	3	6 1/2
MD092	5	6 1/2
MD093	6	6 1/2
MD094	7	6 1/2
MD095	8	6 1/2
MD096	9	6 1/2
MD097	10	6 1/2
MD098	11	6 1/2
MD099	12	6 1/2




GF923R-GF925R

FRAZIER Suction Tube W/ Luer Hub

Catalog Number	Tip Size (FR)	Length (in.)
GF923R	9	7
GF925R	15	7

- Teardrop shaped thumb control for a very precise suction regulation
- Malleable material for individual forming of the suction hose
- Conical design of the tube prevents plugging
- Delivered with stylet



GF391R-GF429R

Suction Tube, Fukushima Style,
W/ Tapered Teardrop Control

Catalog Number	Size	Tip Size (mm)	Tip Size (FR)	Working Length (in.)	Length (in.)
GF401R	S	1.0	3	4	6 1/2
GF402R	S	1.4	4	4	6 1/2
GF403R	S	1.7	5	4	6 1/2
GF404R	S	2.0	6	4	6 1/2
GF405R	S	2.3	7	4	6 1/2
GF406R	S	2.7	8	4	6 1/2
GF407R	S	3.0	9	4	6 1/2
GF408R	S	3.3	10	4	6 1/2
GF409R	S	4.0	12	4	6 1/2
GF391R	M	1.0	3	4 1/2	7
GF392R	M	1.4	4	4 1/2	7
GF393R	M	1.7	5	4 1/2	7
GF394R	M	2.0	6	4 1/2	7
GF395R	M	2.3	7	4 1/2	7
GF396R	M	2.7	8	4 1/2	7
GF397R	M	3.0	9	4 1/2	7
GF398R	M	3.3	10	4 1/2	7
GF399R	M	4.0	12	4 1/2	7
GF411R	L	1.0	3	5 1/2	8
GF412R	L	1.4	4	5 1/2	8
GF413R	L	1.7	5	5 1/2	8
GF414R	L	2.0	6	5 1/2	8
GF415R	L	2.3	7	5 1/2	8
GF416R	L	2.7	8	5 1/2	8
GF417R	L	3.0	9	5 1/2	8
GF418R	L	3.3	10	5 1/2	8
GF419R	L	4.0	12	5 1/2	8
GF421R	XL	1.0	3	6 1/2	9
GF422R	XL	1.4	4	6 1/2	9
GF423R	XL	1.7	5	6 1/2	9
GF424R	XL	2.0	6	6 1/2	9
GF425R	XL	2.3	7	6 1/2	9
GF426R	XL	2.7	8	6 1/2	9
GF427R	XL	3.0	9	6 1/2	9
GF428R	XL	3.3	12	6 1/2	9
GF429R	XL	4.0	12	6 1/2	9

General Surgery

Suction Tubes




GF980R

olive diam. 6 - 9 mm

1/2

LANDOLT Dissecting Suction Tube

Catalog Number	Tip	Tip Size (mm)	Length (in.)
GF980R	ANG, SHARP	3.2	7 7/8



MB635R-MB636R

1/2

SACHS Suction Tube

Catalog Number	Tip Size	Length (in.)
MB635R	SMALL	8 1/2
MB636R	LARGE	8 1/2



US106-US107

1/2

PRIBRAM Suction Tube with Thumb Control

Catalog Number	Tip	Length (in.)
US106	16	5 1/2
US107	18	5 1/2



MG100R-MG110R

1/2

OCHSNER Trocar


Catalog Number	Tip Size (FR)	Length (in.)
MG100R	10	6
MG101R	12	6
MG102R	14	6
MG103R	16	6
MG104R	18	6
MG105R	20	6
MG106R	22	6
MG107R	24	6
MG108R	26	6
MG109R	27	6
MG110R	28	6



AA804R-AA806R

Steel Ruler, Graduated In Millimeters and Inches



Catalog Number	Length (in.)
AA804R	6
AA805R	11 3/4
AA806R	19 3/4



AA808R

Ruler, F/ X-Ray Work

Catalog Number	Length (in.)
AA808R	19 3/4



AD110R

Finger-Goniometer

Catalog Number
AD110R

General Surgery

Accessories



MG001R-MG004R

Surgical Spoon

Catalog Number	Size	Blade Size (mm)	Length (in.)
MG003R	SMALL	25	4 1/4
MG002R	MEDIUM	30	5 1/4
MG004R	MEDIUM-LARGE	35	6 1/4
MG001R	LARGE	40	7 1/2




783101-783102

1/2

Staple Remover

Catalog Number	Description
783102	REUSABLE
783101	SINGLE USE (BOX OF 6)



MG377R-MG378R

TAYLOR Percussion Hammer

Catalog Number	Size	Length (in.)
MG377R	SMALL	7
MG378R	LARGE	7 3/4

Metal Sponge Bowl

Catalog Number	Capacity (qts.)
US061	1
US060	1 1/2
US062	1 7/8



US060-US062

1/2

Metal Cup

Catalog Number	Item	Capacity (oz.)
US066	MEDICINE CUP	2
US063	IODINE CUP	6
US067	SOLUTION BASIN	12
US065	IODINE CUP	14
US064	GRADUATED MEASURE	32



US063-US067

1/2

Kidney Tray

Catalog Number	Length (in.)
JG504R	6 3/4
JG506R	9 3/4



JG504R-JG506R

General Surgery

